



# MADE BY SCHICK

THE BEST FOR YOUR LABORATORY



**SCHICK**  
DENTAL

# create the future

## List of contents

HANDPIECES		
<b>Q3</b>	▶ User-friendly & practice-oriented	4
	▶ Technical data	6
	▶ Q3 Powerstick	7
	▶ Handpiece for left-handed	8
	▶ ISO-motor	8
	▶ Q Basic slim	9
<b>Q Mobile</b>	▶ Highest precision without grid connection	10
<b>Accessories</b>	▶ Riveting hammer for micromotors	12
	▶ Selection - practical processing sets	14
WORKPLACE LAYOUT		
<b>Absorberbox</b>	▶ Compact, quiet and top illumination	16
<b>L Protect</b>	▶ The ultimate combination of protection and illumination	18
	▶ The one of a kind LED swan neck lamp – L Protect SH-D	19
<b>Loupe</b>	▶ Outstanding sharpness and clarity	20
<b>Air puffer system</b>	▶ Clear choice	21
MILLING TECHNOLOGY		
<b>S1 Basic</b>	▶ The compact milling unit	22
<b>S2 Profi</b>	▶ The all-round milling unit	24
<b>S3 Premium</b>	▶ The masterclass	26
<b>S Technology</b>	▶ Technical data	28
<b>S Ceramic</b>	▶ The ceramic milling set	30
	▶ Technical data	32
<b>S W1</b>	▶ Watercooling for milling spindle	34
MODEL PROCESSING		
<b>G2 Concept</b>	▶ Time-saving model saw	36
	▶ Technical data	38

All for my Lab

# Craftsmanship

at its best

Schick, a family-owned company, has developed in around 60 years to become one of the world's leading premium manufacturers of micromotors, milling machines and devices for model processing. Driven by creativity, quality consciousness and German inventive spirit, the development engineers continue to set new standards in the industry in collaboration with leading dental technicians.

## Your benefits at a glance

- ▶ Longlasting premium products guaranteed developed and assembled in Germany
- ▶ Personal point of contact
- ▶ Reliable service
- ▶ Own service center
- ▶ Longtime availability of spare parts
- ▶ 1 year warranty and 5 years warranty on the electric motor (stator and rotor magnet) of the current handpieces and milling spindles





User-friendly & practice-oriented

## Outstanding performance

The new Q3 series from Schick offers unbeatable performance parameters with up to 60,000 rpm and 11 Ncm. Thanks to new technology, the speed can be controlled extremely sensitively, making it ideal for processing even the most demanding dental materials.

## Multifunctional and reliable

All important parameters are clearly displayed on the bright OLED display and are easy to read from all angles. Additional functions such as clockwise/anticlockwise rotation and cruise control leave nothing to be desired and can be conveniently activated via direct selection keys.

## Three power classes

With three different performance classes, the Q3 covers the requirements of a wide range of applications. The three unit types are supplied with the following handpieces as standard:

- ▶ The **Q3 Basic** with Q Basic handpiece, which is perfectly suited for universal use in everyday laboratory work thanks to its robust design.
- ▶ The **Q3 Profi** with Q Profi motor handpiece, which thanks to the soft-grip surface lies comfortably in the hand and is therefore very popular and widespread.
- ▶ The **Q3 Premium** with Q Premium handpiece - the ultra-light and ergonomically perfectly balanced motor handpiece, which with its robust interior made of corrosion-resistant stainless steel is also suitable for the heaviest grinding, cutting and milling work.

**With Q3 you get absolute top technology  
Made by Schick.**



Basic





Profi

Premium



### Q3 Basic

Knee control complete Art. No. 10765  
Foot control complete Art. No. 10768  
**Table-top control complete Art. No. 10773**

### Q3 Profi

Knee control complete Art. No. 10775  
**Foot control complete Art. No. 10778**  
Table-top control complete Art. No. 10783

### Q3 Premium

**Knee control complete Art. No. 10785**  
Foot control complete Art. No. 10788  
Table-top control complete Art. No. 10793

As option available:



Footswitch (on/off)  
Dynamic foot control



Art. No. 6370/2  
Art. No. 9440





## Technical data

	Q3 Basic	Q3 Profi	Q3 Premium
Speed range	200 – 40,000 rpm	200 – 50,000 rpm	200 – 60,000 rpm
Torque	7 Ncm	9 Ncm	11 Ncm
Operating voltage	100 – 240 Volt	100 – 240 Volt	100 – 240 Volt
Operating frequency	50/60 Hz	50/60 Hz	50/60 Hz
Power consumption	160 Watt / max. 180 Watt	180 Watt / max. 240 Watt	200 Watt / max. 340 Watt
Overload protection	No waiting time when restarting	No waiting time when restarting	No waiting time when restarting

**Scope of delivery:** Control unit, handpiece with cable, handpiece rack, mains cable

Handpiece	Q Basic (9400/3)	Q Profi (9300)	Q Premium (10720)
Diameter	Max. 30 mm	Max. 28 mm	Max. 29 mm
Length	158 mm	158 mm	144 mm
Weight without cable	247 g	256 g	195 g
Weight with cable	378 g	384 g	330 g
Drive system	Brushless DC-Motor	Brushless DC-Motor	Brushless DC-Motor
Concentricity	< 0.02 mm	< 0.02 mm	< 0.02 mm
Chucks	Ø 2.35 mm (standard) Ø 3.00 mm (optional)	Ø 2.35 mm (standard) Ø 3.00 mm (optional)	Ø 2.35 mm (standard) Ø 3.00 mm (optional)
Tool change	Quick release	Quick release	Quick release

Control units			
Dimensions	Knee device	Foot device	Table-top device
Width	128 mm	175 mm	106 mm
Height	208 mm	95 mm	200 mm
Depth	207 mm	260 mm	207 mm
Weight	1,360 g	2,300 g	1,880 g





Unique solution for more performance

### More Power

With the newly developed Powerstick from Schick, the two power classes Q3 Basic and Q3 Profi can be retrofitted with additional power if required.

### Your benefits

- ▶ Increase in speed by 10,000 rpm and 2 Ncm torque
- ▶ Simple and fast use of the Powerstick- no conversions required
- ▶ Subsequently a powerful plus in performance without the need to purchase an additional handpiece
- ▶ Not device-specific, can be used flexibly at different laboratory workplaces



**Q3 Powerstick**  
Art. No. 10745

# Left-hand handpiece

Unrestricted working

## Optimal handpiece for left-handers

Schick offers convenient laboratory handpieces also for left-handers: The left-hand handpiece by Schick optimally supports dental technicians in their work. When used in conjunction with special anticlockwise milling cutters, the cutting or working direction corresponds to that of right-handed operators with clockwise handpieces and thus facilitates extraction. A special control unit is not required because the direction of rotation is switched by the special anticlockwise configuration on the handpiece.

## Your benefits

- Unlimited maximum speed
- adaptable use
- Non-fatiguing working
- Connects to all Q/QUBE/Q3 control units from Schick\*

Art. No. 9002/08



Technical data	
Diameter	Max. 28 mm
Length	158 mm
Weight without cable	256 g
Weight with cable	384 g
Drive system	Brushless DC-Motor
Concentricity	< 0.02 mm
Chucks	Ø 2.35 mm (standard) Ø 3.00 mm (optional)
Tool change	Quick release

# Q-Motor with ISO-Coupling



## For Schick control units

The ISO-motor enables the usage of dental angular handpieces with the powerful Schick control units. You can attach the regular dental contra angle handpieces or the fiber optics handpieces to this ISO-motor.

## Your benefits

- Connects to all Q/QUBE/Q3 control units from Schick\*

Art. No. 9002/04

Technical data	
Diameter	Max. 28 mm
Length	115 mm
Weight without cable	152 g
Weight with cable	282 g



# Basic slim

The slim alternative

## Slim and variable

The Q Basic slim motorhandpiece is an alternative to the standard handpieces from Schick. The lean design allows users to hold the handpiece as needed for the task in the laboratory. The unique design of the handpiece tip allows users to access hard-to-reach places.

## Your benefits

- ▶ Pen and fist grip enable flexible working
- ▶ Especially suited to small hands
- ▶ Connects to all Q/QUBE/Q3 control units from Schick\*

Art. No. 9400/13



### Technical data

Diameter	Max. 28 mm
Length	148 mm
Weight without cable	230 g
Weight with cable	360 g
Drive system	Brushless DC-Motor
Concentricity	< 0.02 mm
Chucks	Ø 2.35 mm (standard) Ø 3.00 mm (optional)
Tool change	Quick release

\*Older control units (C-Series) with adapter cable

# Mobile

Highest precision without grid connection

## Powerful and mobile

The rechargeable handpiece Q Mobile from Schick combines two worlds for the first time: unlimited mobility thanks to battery power and the quality of a Schick handpiece.

## Your benefits

- ▶ Highest precision and concentricity properties
- ▶ High performance of 5.5 Ncm and 30,000 rpm
- ▶ Quick-release clamping
- ▶ Speed and capacity display
- ▶ Can be switched between left- and right hand motion
- ▶ Optional: connecting cable for separation of battery and handpiece



### Q Mobile

Art. No. 9705

#### Scope of delivery:

Q Mobile motorhandpiece without cable,  
control unit with integrated battery,  
charging station, transport case



## Safe and Practical

Delivery and storage in transport case

## Optionally available

Thanks to the waist bag, Q Mobile together with the connection cable can be worn practically on the body. The bag offers space for the battery control unit and the handpiece.



### Waist bag (without content)

Art. No. 9654

### Connecting cable for Q Mobile handpiece

Art. No. 9604

### Connecting cable for Q Basic, Q Profi, QUBE LV

Art. No. 9604/1

Technical data	
Speed range	3,000 - 30,000 rpm
Torque	5.5 Ncm
Operating voltage	90 - 264 Volt
Operating frequency	50/60 Hz
Overload protection	Electrical

Handpiece	
Diameter	max. 29 mm
Length without cable	120 mm
Weight	194 g
Drive system	Brushless DC-Motor
Concentricity	< 0.02 mm
Chucks	Ø 2.35 mm (Standard) Ø 3.0 mm (Optional)
Tool change	Quick release
Battery	
Width	71 mm
Height	125 mm
Depth	43 mm
Weight	250 g
Capacity	3,500 mAh

# Riveting hammer


For handpieces

## Versatile

This micromotor extension has a multitude of applications and will even separate telescopic and conical crowns which are wedged together tightly. The vibrations created in the extension separate the inner and outer copings without harming their delicate marginal areas.

## Your benefits

- ▶ Quick and easy to handle
- ▶ Robust construction
- ▶ Tempered ram, will not distort
- ▶ Lifetime-lubricated ball bearing
- ▶ For 2.35 mm diameter chucks
- ▶ To be used with nearly all micromotors



**Riveting hammer**  
Art. No. 1850\*



**Riveting hammer with chisel**  
Art. No. 1850/1\*



**Riveting hammer QUBE**  
Art. No. fehlt: 1850/2  
For handpiece Art. No. 9001 (QUBE),  
9003 (QUBE Assist), 9500 (QUBE II),  
10720 (Q Premium), 9705 (Q Mobile)



**Twin-riveting hammer**  
with plastic cap to be used  
on sensitive surfaces  
Art. No. 1860\*



See user video on  
[www.schick-dental.de/rivetinghammer](http://www.schick-dental.de/rivetinghammer)

\*To be used with micromotors  
from other manufacturers

universal



# Selection

Exceptional quality


## Practical processing sets

The Selection processing sets are a collection of high-quality tools in whose development we concentrated on what is crucial. These highquality tools are tested and recommended by dental technicians.

## Less is more

The compact Selection processing sets allow performing all cutting and milling works effectively and with a minimum amount of tool changes – for less is more.

- 1 Conditioning Set Acrylics**  
Art. No. 3905
- 2 Conditioning Set Precious Metal**  
Art. No. 3910
- 3 Conditioning Set Framework**  
Art. No. 3920
- 4 Conditioning Set Plaster**  
Art. No. 3930
- 5 Conditioning Set Porcelain**  
Art. No. 3940



You will find the article numbers of single tools on [www.schick-dental.de](http://www.schick-dental.de)





# Absorberbox XL

## Compact, quiet and top illumination

The new Absorberbox XL from Schick is the cleanest solution for dustintensive work in the laboratory. The extra-large internal volume allows comfortable working with free moving space. Wide rubber cuffs allow easy insertion of handpieces and workpieces. An integrated suction with fine filter ensures always a clear view of the workpiece. For maximum extraction power, the Absorberbox XL can be connected to an external suction. For this purpose the filter unit can simply be exchanged by the connection plate for external extraction. The large-area LED lighting guarantees bright and shadow-free illumination of the working area.

## Your benefits

- ▶ Flexible use in the laboratory thanks to the integrated suction system
- ▶ Extra-large internal volume (21 litres)
- ▶ Simple conversion from internal to external extraction
- ▶ Large-area LED lighting





For external extraction the connection plate can be used on the rear side of the unit.



**Absorberbox XL**  
Art. No. 10835

**Scope of delivery:**  
Absorberbox XL with integrated suction, particle filter, connection plate for external extraction

Technical data	
Power	10 W
Airflow rate	184 m³/h
Sound pressure level	55 dB(A)
Operating voltage	12 V
Nominal voltage	100 - 230 V
Rated frequency	60 Hz
Dimensions (WxHxD)	475 x 275 x 445 mm
Weight (empty)	3.5 kg
Illuminance	2,500 Lux
Colour temperature	6,000 K



# L Protect

The ultimate combination of protection and illumination

## Illumination at a max

The LED illuminated safety screen offer a unique combination of energy-saving workplace illumination and a protective screen made of toughened safety glass (TSG). Equipped with highly efficient SMD-LEDs L Protect offers an outstanding illumination of the working area without irritating shadows. At the same time it ensures a maximum occupational safety.

With a luminosity of 2500 lux (brighter than two 60 watt incandescent bulbs) and minimum energy requirements of 8 watts L Protect is an economical alternative or supplement to conventional workplace lighting systems.

## Your benefits

- ▶ Bright and shadow-free illumination
- ▶ Long-lasting and energy-saving SMD-LED's
- ▶ Special protective glass pane (toughened safety glass, TSG)
- ▶ The LED illumination frame and protective glass can be separated easily
- ▶ Compatible with conventional suction devices
- ▶ It can also be used without a suction device with the optionally available bracket



**L Protect complete**  
Art. No. 4460

**Scope of delivery:**  
Protective glass pane,  
LED illumination frame,  
adapter unit



**Bracket**  
(without L Protect)  
Art. No. 4474

For mounting without  
suction hood

# L Protect SH-D

One of a kind LED swan neck lamp with a protective glass pane

## Individual positioning

L Protect SH-D is enabling a flexible positioning in any direction due to the cleverly constructed flexible arm. It can be fixed to any workspace with a table clamp or a table mount.



### L Protect SH-D

Table clamp  
Art. No. 4479

#### Scope of delivery:

Protective glass pane,  
LED illumination frame,  
Flexible arm,  
Table clamp



### L Protect SH-DTF

Table mount  
Art. No. 4479/1

#### Scope of delivery:

Protective glass pane,  
LED illumination frame,  
Flexible arm, Table mount



Technical data	
Illuminance	2,500 Lux
Colour temperature	6,000 K
Light colour	Daylight white
Nominal voltage	100 – 240 Volt
Operating frequency	50/60 Hz
Input capacity	8 Watt
Operating voltage	12 Volt
Weight SH-D	1405 g
Weight SH-DTF	1020 g

LED illumination frame	
Height	240 mm
Width	190 mm
Protective glass pane	
Height	205 mm
Width	160 mm
Thickness	4 mm
Single safety glass pane acc. to DIN-12150-1	
Flexible arm	
Length	500 mm

# Loupe

Outstanding sharpness  
and clarity

## Versatile use

The loupe from Schick is adaptable and easy to use. It can be attached to any smooth surface by means of suction cups. The cups are easy to remove while offering maximum retention at the same time. In combination with the illuminated L Protect protective glass from Schick, the high-quality loupe with its 3 x magnification is ideally suited to ultra-fine work.



### Loupe

Art. No. 4430

### Scope of delivery:

Lens, two suction cups, holder

Technical data	
Diameter	75 mm
Length	140 mm
Height	83 mm
Depth	17 mm
Weight	50 g
Magnification	3 x




# Air puffer system

A clear choice

## Quick and comfortable

This handy Schick device is used for cleaning and cooling of all materials and tools right at the workstation. It is activated via a foot switch – the inconvenience of having to reach for the air cleaning nozzle is a thing of the past.

The air puffer is attached to the suction hood or the workstation and features ball joints for flexible positioning. The grinding dust is therefore blown directly towards the extraction system. The air intensity can be easily adjusted via a knob placed on the connection block.



### Air puffer system complete

Art. No. 4400

#### Scope of delivery:

Connection block with articulated hose, foot switch, two connecting cables

#### Technical data

Inlet pressure	3 – 6 bar, 43 – 87 psi
Compressed air connection	NW 7.2
Weight	50 g
Articulated hose length	130 mm
Connecting cables	2 m



See user video on  
[www.schick-dental.de/airpuffer](http://www.schick-dental.de/airpuffer)

# S1 Basic

The compact milling unit

## Ergonomic work

The proven S1 Basic milling unit is not only suitable for milling, measuring, transferring or blocking out, but also ideally suited for the telescopic technology. The heightened base ensure excellent ergonomics. For even more working convenience, a foot switch can be connected to the control unit.

## Utmost precision

Another special feature is the optional LED light head for comfortable working. Other reliable features are the doublejointed arm for unrestricted flexibility and the ball bearing guides for highest precision during drilling.

## S1 Basic Parallelometer

The S1 Basic parallelometer for surveying, transferring and blocking out can be converted into a complete milling unit easily and at low cost. Individual components perfectly matched to your needs are also available for wet grinding applications.



S1 Basic Parallelometer

**High-performance  
milling spindle**  
with extremely quiet running

**Rotable grip**  
for comfortable work

**High-precision  
supporting surface**  
quick-release toggle  
clamp for simple fixing of  
the model table



## Utmost precision

during drilling and milling, thanks to the ball-bearing guides

## Double-jointed arm

for unrestricted flexibility

## Ergonomic working height

### Control unit:

- ▶ Latest technology with up to 40,000 rpm
- ▶ Compact dimensions
- ▶ Optional connection for a foot switch



The all-round milling device

## Safe and perfect millings

With the new S2 Profi, we succeeded in combining high quality and practical parts with a compact and state-of-the-art design. The extremely solid construction of the milling arm and column prevents any vibrations and thus enables the best milling results. For more work comfort the model table can be fixed magnetically to the milling table at the push of a button.

## Your benefits


- ▶ Powerful milling spindle running ultra-smooth and quiet
- ▶ LED light head for a perfect illumination
- ▶ Magnetic fixation table
- ▶ Solid milling arm and column for a vibration free milling process
- ▶ Optional footswitch for better working comfort



**Depth stop**  
for precise adjustment  
of the drilling depth

**Double-jointed arm**  
for unrestricted flexibility

**Magnetic  
fixation**  
for model table



# S3 Premium

The masterclass

## The number 1

The S3 Premium is the further development of the proven S3 Master milling machine. The new milling spindle with high-performance electronics permits speeds of up to 80,000 rpm and opens up additional machining further machining possibilities.

The S3 Premium uncompromisingly combines the precision and quiet running of an electric milling spindle with the high material removal rate of a turbine. Whether wax milling at 500 1/min, non-ferrous milling at 50,000 1/min, or zircon machining under water cooling at 80,000 rpm, the S3 Premium covers the entire spectrum of milling technology.

## High-performance milling spindle

- ▶ Speeds up to 80,000 1/min
- ▶ Low starting speed of 500 1/min for comfortable wax milling

## 3D milling arm

- ▶ Magnetic lock via foot via foot pedal switch or operating satellite
- ▶ Highest precision and mobility due to ball bearing guidance
- ▶ 25 mm drill travel

## Operating unit

- ▶ Flexible connection allows positioning to the optimum ergonomic, individual working position
- ▶ All adjustments are clearly visible and controllable during the working process
- ▶ Switching function for optional SW1 water cooling for milling spindle

## Height-adjustable milling table

- ▶ Ensures constant working at eye level even with different heights of models
- ▶ Magnetic fixation for model table

## Two footswitches

for better workflow

- ▶ Milling spindle On/Off
- ▶ Switchable between model table fixation and 3D milling arm lock

## Optional arm rests

- ▶ Individual three-dimensionally adjustable arm support for fatigue-free precision milling

**Depth stop**  
with micrometer spindle

**3D milling arm**  
with magnetic lock





\* Optional accessories,  
not included in scope of delivery

# S Technology

## Technical data

### S1 Basic

	milling unit	control unit
Width	220 mm	106 mm
Height	400 mm	200 mm
Depth	270 mm	207 mm
Weight	7.2 kg	1,880 g
Operating voltage		100 – 240V
Operating frequency		50/60 Hz
Power consumption		160 Watt / max. 180 Watt
Speed range		1,000 – 40,000 rpm
Concentricity	< 0.015 mm	
Collet chuck	Ø 2.35 mm series incl. depth stop for short burs Ø 3.00 mm on request	



**S1 Basic milling unit cpl.** Art. No. 2800

**Scope of delivery:**  
 S1 Basic unit  
 S1 Basic milling spindle  
 Q3 Basic control unit  
 Model table stainless steel

### S2 Profi

Width	240 mm
Height	507 mm
Depth	345 mm
Weight	13 kg
Operating voltage	100 – 240V
Operating frequency	50/60 Hz
Power consumption	80 Watt
Speed range	1,000 – 40,000 rpm
Concentricity	< 0.015 mm
Collet chuck	Ø 2.35 mm series incl. depth stop for short burs Ø 3.00 mm on request

### S3 Premium

Width	300 mm
Height	500 mm
Depth	420 mm
Weight	21 kg
Operating voltage	100 – 240V
Operating frequency	50/60 Hz
Power consumption	200 Watt / max. 340 Watt
Speed range	500 – 80,000 rpm
Concentricity	< 0.015 mm
Collet chuck	Ø 2.35 mm series incl. depth stop for short burs Ø 3.00 mm on request



**S1 Basic parallelometer** Art. No. 2810

**Scope of delivery:**  
 S1 Basic unit  
 Measuring spindle  
 Lead holder 2.35mm



**S2 Profi milling unit cpl.** Art. No. 2950

**Scope of delivery:**

S2 Profi milling unit,  
S2 Profi milling spindle with light head  
Model table is not included




**S3 Premium milling unit cpl.** Art. No. 12500

**Scope of delivery:**

S3 Premium milling unit,  
S3 Premium milling spindle with light head,  
footswitch magnetic control,  
footswitch motor on/off  
Model table is not included


## Accessories




**Measuring spindle**  
adjustable from 0.5 – 3.00 mm  
Ø 22 mm  
Art. No. 2052/1




**Measuring set**  
Ney 2.35 mm  
Art. No. 10450



**Model table**  
Art. No. 2407/9



**Adjustable angle**  
Art. No. 2506




**Lead holder**  
Ø 2.35 mm Art. No. 2268  
Ø 3.00 mm Art. No. 2268/1



**Holding tray**  
for magnetic fixation of the  
milling bowl on the milling unit  
Art. No. 2509



**Milling tray**  
Art. No. 2507/1



**Milling set 2.35 mm**  
10 pieces  
Art. No. 2530/1



**Transfer unit**  
Ø 2.35 mm Art. No. 2795  
Ø 3.00 mm Art. No. 2795/1



**Foot switch on/off**  
Art. No. 6370/2



**Light head**  
Art. No. 2510



**Power supply light head**  
Art. No. 2480/1



**Wooden armrest compl.**  
Art.-Nr. 12560

- Accessories for S1 Basic
- Accessories for S2 Profi
- Accessories for S3 Premium
- Accessories for S3 Premium and S3 Master
- For milling units of other manufacturers

# S Ceramic

The ceramic milling unit

## Safe and clean

All milling units from Schick are upgradable with the specially adapted accessories into a fully workable ceramic milling unit. This special ceramic milling accessory is also available for milling units of other manufacturer.

## Turbine


In place of a motor milling spindle a turbine with 300,000 rpm is used. Additionally the turbine is equipped with the proven LED lighting technology of Schick. To avoid micro cracks in the ceramic material, water cooling is used. A suction tub enables the technician to use this system in his standard milling device.

## Special tools

To complete the system one set of each diamond tools and special polishers for ceramic are available for the rework of ceramic primary parts and implant abutments.

## Suction

Arising dust and water spray is removed completely by the ergonomic suction tub in combination with the special water separator. User and milling unit are protected efficiently. By using the special separator, the system can be operated with all laboratory suctions or cental suctions. Alternatively you can connect the suction tub to a wet-dry vacuum cleaner, too.



300,000  
rpm

### Adapter

for turbines and milling units of other manufacturers available upon request

### Turbine


with LED light head

### Model table

made of stainless steel

### Suction tub

for connection to extraction system via the separator



See user video on  
[www.schick-dental.de/ceramic](http://www.schick-dental.de/ceramic)

# S Ceramic

## Technical Data

### Collection tub

Width x height x depth	260 x 245 x 110 mm
Weight	0.83 kg

### Suction tub

Width x height x depth	300 x 245 x 125 mm
Weight	1.3 kg
Intake socket (suction support)	31 mm – 38.1 mm (conical)

### Separator

Width x height x depth	405 x 405 x 540 mm
Weight	5.75 kg
Connecting adaptor	Inside Ø 31.4 mm Outside Ø 37.1 mm

Hose length	1 x 2 m, 1 x 3.5 m
Hose weight	2.0 kg incl. adaptor
Ø hose	Inside 32.7 mm Outside 40.0 mm

### T 100 turbine control unit

Width x height x depth	295 x 192 x 155 mm
Air pressure	5 - 8 bar
Tank capacity	1.4 l
Foot pedal	Pneumatik
Weight	3,200 g

### Turbine handpiece

Speed	300,000 rpm
Air volume	35 l/min average
Spray	Max. 100 ml/min
Collet chuck	Ø 1.60 mm (FG)
Weight	60 g

Subject to technical modification without prior notice



**Ceramic milling set for S1 and other manufacturers**  
Art. No. 2650/15

#### Scope of delivery:

Suction tub, separator, turbine T100, model table stainless steel, light head for turbine, power supply light head, diamond tool set for turbine 1.6 mm (8 pcs.), polishing set 2.35 mm (3 pcs.), adapter for turbine



**Ceramic milling set basic**  
Art. No. 2650/7

#### Scope of delivery:

Suction tub, separator, model table stainless steel




















### Ceramic milling set complete for S2 and S3

Art. No. 2650/05

**Scope of delivery:** Suction tub, separator, turbine T100, model table stainless steel, light head for turbine, diamond tool set for turbine 1.6 mm (8 pcs.), polishing set 2.35 mm (3 pcs.), adapter for turbine

Create your own set - according to your needs

 <p><b>S1 Basic unit</b> Art. No. 2805</p>	 <p><b>Collection tub without nozzle</b> Art. No. 2498</p>	 <p><b>Suction tub with nozzle</b> Art. No. 2470/5</p>	 <p><b>Separator</b> Art. No. 2655</p>
	 <p><b>Turbine T100</b> Art. No. 2640/1</p>	 <p><b>Model table stainless steel</b> Art. No. 2407/9</p>	 <p><b>Diamond tool set for turbine 1.6 mm (8 pcs.)</b> Art. No. 2660</p>
 <p><b>Light head for turbine</b> only possible in conjunction with 2480/1 when used in S1 2805 operation Art. No. 2510/1</p>	 <p><b>Power supply light head</b> Art. No. 2480/1</p>	 <p><b>Polishing set 2.35 mm (3 pcs.)</b> Art. No. 2665</p>	 <p><b>Adapter for turbine T100</b> Art. No. 2481</p>
 <p><b>Adapter for turbine NSK Aqua Presto</b> Art. No. 2481/1</p>	 <p><b>Adapter milling unit Degussa F1, F3, for Turbine T100</b> Art. No. 2481/6</p>	 <p><b>S W1 Watercooling for milling spindle</b> Art. No. 2990</p>	 <p><b>Connecting cable for S W1 together with S3 Premium</b> Art. No. 12548</p>

Adapter for turbines and cutters of other manufacturer on request.

■ Accessories for S1 Basic 
 ■ Accessories for S2 Profi 
 ■ Accessories for S3 Master and S3 Premium 
 ■ Accessories for S3 Premium 
 ■ For milling units of other manufacturers

# S W1

## Watercooling for milling spindle

### Optimal for zero-degree milling

The new S W1 water cooling in combination with all current Schick milling machines is the ideal alternative solution to wet milling with a turbine. For the first time, it offers the user the option of machining all-ceramic primary parts without the need for timeconsuming conversion of the milling machine. The water feed is controlled by means of a rotary knob integrated into the unit. Compared to a compressed air driven turbine, the motor milling spindle achieves absolutely precise results, especially in the case of zero-degree milling. Thus, the combination of a Schick motor milling spindle and S W1 water cooling fully meets the increased requirements in the field of digital fabrication of ceramic primary parts.

### Your benefits

- ▶ Cost-effective alternative to turbines
- ▶ Adjustable water feed
- ▶ No compressed air required
- ▶ Easy filling thanks to removable water tank
- ▶ Highest precision due to milling spindle without bearing clearance
- ▶ wear-free, low-maintenance spindle technology

### Important information:

- ▶ S W1 is suitable for Schick milling machines
- ▶ S 1 Basic, S2 Master, S2 Profi, S3 Master and S3 Premium
- ▶ Recommendation: Conversion of your existing Schick milling spindle to a high-grade steel version to prevent any possible corrosion
- ▶ Power supply for light head (art.-no. 2480/1) is needed if you want to use the S W1 with your S 1 Basic
- ▶ Additionally required for S3 Premium: connection cable Art. No. 12548

#### Technical data

Operating voltage	100-240 V
Operating frequency	50/60 Hz
Dimensions (WxHxD)	170 x 170 x 85 mm
Weight (empty) max.	1.4 kg
Pumping capacity	60 ml/min
capacity	max. 0.5 l

**S W1 Watercooling for Schick milling spindle**  
Art. No. 2990

#### Scope of delivery:

S W1 Station, water tank, S W1 light head with hose connection, reducer for chuck 3 mm to 1,6 mm, reducer for chuck 2,35 mm to 1,6 mm



## S W1 Ceramic milling set

**S W1 Ceramic milling set for S1 Basic milling unit**  
Art. No. 2990/1

**Scope of delivery:**

S W1 Watercooling for Schick milling spindle, collection tub without nozzle, model table stainless steel, power supply light head for S1 Basic

**S W1 Ceramic milling set for S2 und S3 milling units**  
Art. No. 2990/2

**Scope of delivery:**

S W1 Watercooling for Schick milling spindle, collection tub without nozzle, model table stainless steel



**Power supply**  
for automatic activation of  
dust extraction devices

# G2 Concept

Cut into the future

## Exact and time saving

Until now, many models have necessitated saw cutting by hand, due to the limited space between the preparation boundary or irregularities in the alignment of stumps. Superior to handsaws and diamond disc equipped handpieces, the G2's cutting power and precision dramatically shorten the previously time-consuming task of sectioning models. Within minutes, technicians can section models with smooth and straight die walls, without damaging margins. Unintentional start-ups and improper hand placement are avoided with the special touch sensors and laser hand-placement guides.

## Your benefits

- ▶ You can create models with smooth and straight die walls by cutting upwards – without damaging the preparation boundary
- ▶ No extra work needed thanks to straight die walls
- ▶ Suitable for model plaster and special model resin materials
- ▶ Within minutes, technicians can section models with smooth and straight die walls
- ▶ Suitable for nearly all model systems – with or without pins

### G2 Concept model saw cpl

Art. No. 1500/5

**Scope of delivery:** G2 Concept model saw  
with integrated line-marking laser, mains cable, diamond saw  
blade (Art. No. 1967/2)

unique


**Sensor I**  
for particularly safe  
two-hand operation

**Diamond saw blade**  
developed for high-precision  
and thin saw cuts

**Line-marking  
laser**  
for optimum user  
comfort and precise cut

**Connection  
for the dust  
extraction system**  
less air pollution and  
better working conditions  
protect your health

**Non-slip pad**  
for secure fixing  
of the models



**Ergonomic  
streamlined design**  
for an effortless operation

**Saw cut from below**  
prevents damage to  
the adjacent teeth or the  
preparation boundary

**Sensor II**  
prevents injuries with  
a special two-hand  
safety circuit

# G2 Concept

Functional and precise

## How it works

- The model is placed on the saw table and positioned with both hands. The integrated line marking laser enables exact cut positioning. When the saw table is lowered, the saw blade dips into the model from below, resulting in a clean cut parallel to the pins. The exact depth of a cut can now be determined. In the case of divergent tooth position, sawing is performed up to the preparation margin and manual fracture is performed at the generated predetermined fracture point.
- Since the saw cuts are parallel to each other, all model segments can be removed from the model without any problems.

## Maximum working safety

The model saw is equipped with two touch sensors, to prevent unintentional start-up of the saw. The process is only started when both hands are in a safe working position. In addition, the G2 Concept model saw can be connected to any commercially available laboratory suction.

## Technical data

Width	375 mm	Speed	7,500 rpm
Height	160 mm	Diamond saw blade	Ø 85 x 0,3 mm
Depth	390 mm		segmented
Weight	4.5 kg	Width of cut	0.3 mm
Operating voltage	230 V / 115 V / 100 V	Cutting depth	max. 26.5 mm
Operating frequency	50 / 60 Hz	Laser category	3R
Power consumption	130 Watt	Suction support	Ø 38 mm



See user video on  
[www.schick-dental.de/g2code-en](http://www.schick-dental.de/g2code-en)






You can save  
time and money

- ▶ G2 Concept model saw  
90 sec\*
- ▶ Saw cuts by hand  
180 sec\*\*

\*Determined in a practical test with  
a 11 segment model and 10 saw cuts



See user video on  
[www.schick-dental.de/g2comparison](http://www.schick-dental.de/g2comparison)



MADE  
IN GERMANY  
BY SCHICK

Development, production and  
assembly inhouse.

## Your benefits

- Personal point of contact
- Own service center
- Longtime availability of spare parts



[www.schick-dental.de](http://www.schick-dental.de)



Please scan here  
for Download

Our complete catalog in  
other languages can be  
found online!



**SCHICK**  
DENTAL

Schick GmbH  
Lehenkreuzweg 12  
D-88433 Schemmerhofen  
Telefon +49 7356 9500-0  
Telefax +49 7356 9500-95  
[info@schick-dental.de](mailto:info@schick-dental.de)  
[www.schick-dental.de](http://www.schick-dental.de)



Your local dealer

Subject to technical modification without prior notice