# THERMEO®

the smart mate by

pro **3d** ure  
medical

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## THERMEO®

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
### THERMEO® is flexible.

Who doesn't know it - „crooked impressions“ during conventional impression taking or imperfections during intraoral scanning that lead to „data holes“ on the surface of the generated data set? A conventional hard splint then wobbles on the plaster model and thus in the patient's mouth. The splint is subjected to greater stress at this point, the risk of fracture is significantly increased and the patient complains about a lack of wearing comfort. The unique THERMEO® splint material compensates for such inaccuracies with its patented thermomemory effect. The THERMEO® material becomes flexible at body temperature and adapts optimally to the tooth surface. This increases the wearing comfort of the splint and minimizes the risk of breakage. THERMEO® thus supports the satisfaction of the dental technician, dentist and patient through higher safety, exceptional wearing comfort and fewer „recalls“.

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## THERMEO® is smart.

Smart materials/responsive materials are materials that react independently to changing environmental conditions such as temperature. This reaction is reversible. The high-tech plastics of the THERMEO® family belong to this group of smart materials and are characterized by a unique patented thermomemory effect. Compared to hard dental splints, inaccuracies during impression taking can thus be compensated with a THERMEO® splint. THERMEO® products can also be used to produce functional therapeutic bite splints. THERMEO® splints are characterized by high wearing comfort and are self-adjusting due to the thermomemory effect. THERMEO®s special properties combine the advantages of hard as well as soft splint materials.

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Fig. 2: THERMEO® time depending stress relaxation at 73.4 °F (23 °C).\*

## THERMEO® is comfortable.

THERMEO® meets the highest medical technology standards. For example, THERMEO® does not use phthalate-based plasticizers, but instead employs a patented flexibilization technology that ensures unique biocompatibility in this class of material. Compared to other suppliers, this results in a significantly greater pressure drop on the tooth after insertion of the splint (Fig. 2 and 3), which in turn increases wearing comfort for the patient. Furthermore, due to the technology used, the thermomemory effect and thus the resilience of the THERMEO® material is greater than that of competitors (Fig. 4). The THERMEO® splint returns to its original position faster and more precisely due to body temperature. THERMEO® can be processed conventionally using the casting or tamping-pressing technique or integrated into the digital workflow in the form of THERMEO® blanks. The blanks are available in thicknesses of 16, 20 and 25 mm. This allows tooth-supported and also soft tissue-supported splints to be manufactured efficiently. Both the powder/liquid system and the milling rounds are additionally available as THERMEO® SO in shades A1, A2 and A3 for the fabrication of „snap-on splints“ (Fig. 1).

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Fig. 3: THERMEO® (orange) and competitor A (red) decrease in stress at 98.6 °F (37 °C) in comparison to 73.4 °F (23 °C).\*

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Fig. 4: THERMEO® (orange) and competitor A (red) elastic recovery after a 90° deformation depending on the temperature, 73.4 °F (23 °C) and 98.6 °F (37 °C).\*


Abb. 5: THERMEO® (orange) und Wettbewerber A (rot) volumetrischer Abbau nach 120.000 Zyklen Kausimulation, mm³.\*  
 Fig. 5: THERMEO® (orange) and competitor A (red) volumetric degradation after 120,000 cycles chewing simulation, mm³.\*



Abb. 6: THERMEO® (orange) und Wettbewerber A (rot) vertikaler Abbau nach 120.000 Zyklen Kausimulation, mm.\*  
 Fig. 6: THERMEO® (orange) and competitor A (red) vertical degradation after 120,000 cycles chewing simulation, mm.\*

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## THERMEO® is durable.

The materials of the THERMEO® family, characterized by their exceptional wearing comfort, have been developed for durability. The stability of the material was documented in an extensive chewing simulation study conducted by LMU Munich. Despite the higher flexibility and higher resilience of THERMEO® compared to its market competitors, comparable abrasion values were determined (Fig. 5 and 6). Thus, patient comfort is increased without reducing the clinical benefit. Furthermore, amine compounds and MMA were avoided in the development of THERMEO® in order to reduce the yellowing of the material and its potential for irritation and sensitization. Furthermore, tested cleaning recommendations for THERMEO® are available to the patient in order to sustainably ensure the special properties over the entire service life of the THERMEO® splint (see next page).

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- exceptional wearing comfort
- unique flexibility by thermomemory effect
- self adjusting
- no yellowing
- extreme durability
- long-lasting
- highest biocompatibility
- "snap on" splints available

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## THERMEO® splint Cleaning and Care

Comprehensive documentation and important tips on the use and care of THERMEO® splints help you and the patients use THERMEO® products for a long and sustainable time.

#### GERMANY

pro3dure medical GmbH  
Am Burgberg 13  
58642 Iserlohn  
phone +49 (0)2374 920050-0  
fax +49 (0)2374 920050-50

#### USA

pro3dure medical LLC  
Valley View Business Center  
9825 Valley View Road  
Eden Prairie, MN 55344  
phone +1 (952) 426 1928  
fax +1 (952) 681 7515

info@pro3dure.com  
www.pro3dure.com

pro **3d** ure  
medical  
idea to product.