Biphasic Calcium Sulfate - Chemical Structure and Regenerative Properties


Calcium sulfate (CS) features a unique position among all regenerative bone graft substitutes. With more than 100 years of documented clinical success, (Dreesmann 1892, Thomas and Puleo 2008) it has a longer history of clinical use than most currently available biomaterials.

In an in vivo scenario, proteins and other biological macromolecules may further retard the setting time to up to 200 min. This is impeding the application procedure as well as the clinical performance for a typical dental application significantly (Ricci 2000).

Although Orsini and co-workers have shown that there are no differences in the bone healing pattern between preset CS dihydrate granules and moldable CS-hemihydrate cement. (Orsini 2004), by using preset CS dihydrate the surgeon loses advantageous handling properties of the moldable cement derivative.

By making use of the Biphasic Calcium Sulfate, the previous setting issues can be solved: The CS hemihydrate component is controlling the cement consistency and moldability characteristics. The dihydrate component is regulating the setting properties.

Using this formulation, the setting process can be reduced from about 20 minutes to 3 minutes, also under in vivo conditions in presence of blood and saliva. The key performance criteria of CS-derived graft materials that was shown in numerous studies is the fast and effective vital bone formation as a result of the rapid dissolution in vivo. (Thomas and Puleo 2008, Toloue, Dahlin) These characteristics provoke advantageous clinical results in small and well-contained bony defects like extraction socket filling. (Toloue, Guarnieri).

Biphasic Calcium Sulfate Composition

Bond Apatite®

Bond Apatite® is composed BCS and synthetic hydroxyapatite granules. With this formulation, the concept of a graft cement can be successfully applied for a wider range of indications.

The BCS affects an injectable graft cement and features fast resorption and new bone formation. The hydroxyapatite granules provide a slow-resorbing osteoconductive scaffold to preserve the augmented volume also over a longer time period.

The BCS matrix is supposed to feature complete resorption. The slow-resorbing hydroxyapatite granule matrix serves as longer range space maintainer for optimized volume control.


Bond Apatite® is delivered in a specially designed ready-to-use syringe containing the granulated powder and physiological saline. Mixing the powder component with the liquid in the driver results in a viscous composite that is suitable for direct injection into the graft site.

Histological Evidence | Complete Regeneration 8 Months Post-Op

At 12 weeks after placing the Bond Apatite®, nearly 90% is replaced by the patient's own bone.

The HA particles occupying 1/3 of the composite graft are intended to slow down the overall resorption rate of the graft. The HA particles do not integrate with the bone, they are first encapsulated by connective tissue, slowly resorb, and later on, the remaining connective tissue under goes ossification.

3 Months Post-Op | Nearly 90% Bone Regeneration


Properties and Clinical Applications of Biphasic, Calcium Sulfate, Yahav et al, April 2020.

8 Months Post-Op | Complete Bone Regeneration

After eight months of healing time, the Bond Apatite® graft is almost entirely resorbed with no encapsulated materials remaining due to the complete transformation of the graft into the new vital bone.

The small and medium-sized HA particles resorb in 3-4 months. In contrast, the larger size particles, which comprise less than 10% of Bond Apatite®, will remain for a prolonged period, almost fully resorbed after eight months, and continues its resorption process up to a complete regeneration.

Properties and Clinical Applications of Biphasic, Calcium Sulfate, Yahav et al, April 2020.

More information and references can be found at: <https://www.augmabio.com/abca/e-book/>

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
- Thomas MV, Puleo DA, Calcium sulfate: A review. *J. Long Term Eff. Med. Implants* 2005;15(69):599-607.
- Thomas MV, Puleo DA, Calcium sulfate: Properties and clinical applications. *J. Biomed. Mater. Res. B Appl. Biomater.* 2009;88(2):597-610.
- Coetzee AS. Regeneration of bone in the presence of calcium sulfate. *Arch. Otolaryngol.* 1980;106:405-409.
- Turri A, Dahlin C. Comparative maxillary bone-defect healing by calcium-sulphate or deproteinized bovine bone particles and extra cellular matrix membranes in a guided bone regeneration setting: an experimental study in rabbits. *Clin. Oral Impl. Res.* 00, 2014, 1–6 doi: 10.1111/clr.12425.
- Toloue SM, Chesnoiu-Matei I, Blanchard SB. A clinical and histomorphometric study of calcium sulfate compared with freeze-dried bone allograft for alveolar ridge preservation. *J Periodontol.* 2012;83:847–855.
- Collins JR, Jiménez E, Martínez C, Polanco RT, Hirata R, Mousa R, Coelho PG, Bonfante EA, Tovar N. Clinical and Histological Evaluation of Socket Grafting Using Different Types of Bone Substitute in Adult Patients. *Implant Dent* 2014;23:489–495.
- Payne JM, Cobb CM, Rapley JW, Killoy WJ, Spencer P. Migration of human gingival fibroblasts over guided tissue regeneration barrier materials. *J Periodontol* 1996;67:236-244.
- Boden SD, Stevenson S. *Bone grafting and bone graft substitutes.* Philadelphia: Saunders, 1999.
- Pietrzak WS, Ronk R. Calcium sulfate bone void filler: A review and a look ahead. *J. Craniofac. Surg.* 2000;11(4):327-333.
- Ricci JL, Alexander H, Nadkarni P, Hawkins M, Turner J, Rosenblum SF, Brezenoff L, De Leonardis D, Pecora G. Biological mechanisms of calcium-sulfate replacement by bone. In: Davies JE , editor. *Bone Engineering.* Toronto: Em Squared Inc.; 2000:332-344.
- Ricci JL. Report on the use of calcium sulfate cement with dental cement. In: *Use of Calcium Sulfate Cement for Dental and maxillofacial Bone Repair*, edited by Pecora, G. and De Leonardis, D., Elite Publishing, Rome, 7 pages, 2001.
- Tay BKB, Patel VV, Bradford DS. Calcium sulfate- and calcium phosphate-based bone substitutes: Mimicry of the mineral phase of bone. *Orthop. Clin. North Am.* 1999;30(4):615-623.
- Anusavice KJ. Gypsum products. In: *Phillips' Science of Dental Materials.* St. Louis, edited by: Elsevier Health Sciences, 2003: 255-281.