

Premium Sterile Surgical Blade

# KIATO<sup>+</sup> PLUS




*The smart choice for sharp blades*


- Carbon Steel Blades made of Swiss steel grade 125Cr1.
- Stainless steel blades made of premium Swedish Steel.
- Individually packed in easy peel color coded foil and boxes.
- Fits all makes of handles. Conforms to ISO 7740.
- Sterilized by Gamma Radiation to a minimum dosage of 2.5 MRads.
- Blades conform to BS 2982 and Medical Device Directive 93/42/EEC.



100 blades per box, 10 box per case, 10 case per master shipping carton.



| Size            | 10    | 10A   | 11    | 11P   | 12    | 12D   | 13    | 15    | 15C   | 15T   |
|-----------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| Carbon Steel    | 18010 | 18009 | 18011 | 18031 | 18012 | 18025 | 18013 | 18015 | 18026 | 18045 |
| Stainless Steel | 18110 | 18109 | 18111 | 18131 | 18112 | 18125 | 18113 | 18115 | 18126 | 18145 |



| Size            | 18    | 19    | 20    | 21    | 22    | 22A   | 23    | 24    | 25    | 26    | 36    |
|-----------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| Carbon Steel    | 18027 | 18019 | 18020 | 18021 | 18022 | 18030 | 18023 | 18024 | 18028 | 18033 | 18029 |
| Stainless Steel | 18127 | 18119 | 18120 | 18121 | 18122 | -     | 18123 | 18124 | 18128 | 18133 | 18129 |

Blade size not to scale.

EC REP  
MDSS GmbH  
Schiffgraben 41  
30175 Hannover  
Germany



aditya  
DISPOMED

Aditya Dispomed Products Pvt. Ltd.  
Plot No. 19 Sector 6 IMT Manesar  
Gurgaon Haryana 122050 India  
Tel. +91 124 47 64 900  
Fax +91 124 47 64 911  
info@adityadispomed.com  
www.adityadispomed.com

A subsidiary of Kehr Surgical Pvt. Ltd.



**kiato** 

The smart choice for sharp blades



Kiato Plus scalpels are made with passion using a state-of-the-art robotized production line in compliance with all regulatory requirements.

Kiato Plus scalpels set the benchmark for all others to follow, from East to West. Rich on features and high on performance, Kiato Plus scalpels are user-friendly with a laser controlled durable sharp edge. In spite of its premium features, use of these scalpels enables hospitals to scale down their budgets. It is, therefore, not surprising that the Kiato Plus range of products is establishing a strong presence in markets worldwide.

Kiato Plus disposable scalpels are manufactured in an ISO 13485:2016 certified factory and are CE marked and gamma sterilised. They are individually packed, 10 to a box, in easy peel pouches. Each of these scalpels have a premium Swedish stainless steel surgical blade fitted to a plastic handle with a safety guard for safe handling.

User-friendly  
Ergonomic and safe  
Premium quality  
Superior performance



REF. 23110





**kiato** 




## User-friendly

The range of Kiato Plus scalpels come in easily differentiated color-coded boxes. Each of these fibre-free, easy peel packs display the blade shape for ease of identification. The blade size is also prominently printed on the handle with a graduated cm scale. Ease of use begins with the product packaging itself and makes for greater convenience.

REF. 23111

REF. 23112





## Ergonomic and safe

Made of ABS, the Kiato Plus scalpel handles are not just ergonomically designed but also stiffer and heavier, to facilitate ease of use. A serrated grip ensures greater comfort and control while the scale on the handle helps in determining incision length. For enhanced safety, each scalpel also has a unique stiff cap with lock tabs that prevents injuries and improves performance.



REF. 23115



REF. 23120

REF. 23121



REF. 23122

## Premium quality

The scalpel handles are molded in a full hot runner mold. The robotized handling system prevents any damage to the handle by gently placing it on the conveyor. Our agreement with Sandvik, Sweden ensures that we receive and use superior stainless steel quality. The blades are made of Sandvik grade 13C26 which conforms to the undermentioned specifications.

Stainless Steel Chemical Composition  
Ref: BS 5194: 1991, Type F

| Constituents    | Percentage (%) |
|-----------------|----------------|
| Carbon (C)      | 0.60 to 0.70   |
| Silicon (Si)    | 0.50 maximum   |
| Manganese (Mn)  | 1.00 maximum   |
| Chromium (Cr)   | 12.00 to 13.50 |
| Phosphorous (P) | 0.03 maximum   |
| Sulphur (S)     | 0.025 maximum  |
| Nickel (Ni)     | 0.50 maximum   |



**Kiato**

24



REF. 23123



REF. 23124



## Superior performance

Sharpness and durability are the true tests of a good quality scalpel. Both these factors are hallmarks of every Kiato Plus scalpel. The laser controlled edge in every scalpel guarantees benchmark sharpness and durability to ensure superior customer satisfaction and unflinching performance, even over time.

23






Aditya Dispomed Products Pvt Ltd  
Plot No. 19 Sector 6 I.M.T. Manesar  
Gurgaon Haryana 122 050 India  
[www.adityadispomed.com](http://www.adityadispomed.com)



MDSS GmbH  
Schiffgraben 41  
30175 Hannover  
Germany



KIATO Safety Scalpels help to minimize the risk of surgical blade injuries that occur before use, during passing, between steps of the procedure, after use and during disposal. The Safety Scalpel is designed to be used single handedly.

### 1. Ship/Pass position:




The retractable sheath securely covers the blade fully. The latch must be released to retract the sheath backwards to the Use Position. The sheath does not retract backwards without activation.

### 2. Use position:



The Surgical Blade is exposed and ready for use. Audible click confirms sheath is fully retracted.




Packed in easy peel rigid blister pouches. 10 scalpels per unit box. Boxes are colour coded for easy identification. 10 boxes per case. 10 case per Master Shipper.

## Directions for use


### To Use:

Press Release button marked On to depress latch on handle and retract sheath fully until it clicks. The blade is now exposed and the sheath is locked in the "USE" position




### To Pass:

Press Release button marked Off and slide sheath forward to initial position until it locks in the pass position.




### To Dispose

Slide Sheath to "pass position". Push dispose button inwards into handle until it clicks & locks.




### Step 1

Dispose in an approved sharps container.



### Step 2

- Premium Swedish Stainless Steel.
- Colour coded handle with Size No clearly printed
- Durable sharpness.
- Clearly printed instructions.
- Cap Safely locked in ship & pass position.
- Permanent Dispose button



### FITMENT 3

| Size | Reorder no |
|------|------------|
| 10   | 12910      |
| 11   | 12911      |
| 11P  | 12931      |
| 12   | 12912      |
| 12D  | 12925      |
| 13   | 12913      |
| 15   | 12915      |
| 15C  | 12926      |

### FITMENT 4

| Size | Reorder no |
|------|------------|
| 18   | 12927      |
| 19   | 12919      |
| 20   | 12920      |
| 21   | 12921      |
| 22   | 12922      |
| 23   | 12923      |
| 24   | 12924      |
| 25   | 12928      |



MDSS GmbH  
Schiffgraben 41  
30175 Hannover  
Germany



**Aditya Dispomed Products Pvt. Ltd.**  
Plot No. 19 Sector 6 IMT Manesar  
Gurgaon Haryana 122050 India  
Tel. +91 124 47 64 900  
Fax +91 124 47 64 911  
info@adityadispomed.com  
www.adityadispomed.com

The smart choice for sharp blades

A subsidiary of Kehr Surgical Pvt. Ltd.