



# PRODUCT CATALOGUE

# 2022

# → INDEX

|       |                                         |
|-------|-----------------------------------------|
| 6-7   | PIEZOSURGERY® <i>touch</i>              |
| 8-9   | PIEZOSURGERY® <i>white</i>              |
| 10    | enzymec                                 |
| 11    | universal cart                          |
| 12-13 | PIEZOSURGERY® spare parts               |
| 14    | insert sets                             |
| 15    | basic set/osteotomy set                 |
| 16    | sinus lift lateral set/piezo lift set   |
| 17    | sinus lift set/ implant prep set pro    |
| 18    | extraction set/ explantation set        |
| 19    | resective perio set/retro surgical set  |
| 20-21 | indications                             |
| 22    | trays                                   |
| 23-25 | implant prep inserts                    |
| 26-29 | osteotomy inserts                       |
| 30-31 | osteoplasty inserts                     |
| 32    | elevation inserts/periosteum inserts    |
| 33    | endo retro inserts                      |
| 34    | extraction inserts/explantation inserts |
| 35    | periodontal inserts                     |
| 36    | implant cleaning inserts                |
| 37    | crown prep tips                         |

38-39 MULTIPIEZO

40-51 MICROPIEZO

42 COMPACT PIEZO LED/COMPACT PIEZO P2K

43-45 ultrasound accessories

46 insert sets

47 set maintenance/set scaling

48 set perio maintenance/set perio treatment

49 set implant cleaning/set endo

50 set endo revision

52 scaling inserts

53 perio universal inserts

54 perio anatomic inserts

55 implant cleaning inserts

56 endo inserts

57 endo revision inserts

58 endo retro inserts/ extraction inserts

59 crown prep inserts

→ ULTRASOUND

## → AIR-POLISHING

60-61 COMBI *touch*

62-63 STARJET

64-65 prophylaxis powders

66 universal cart

68-71 spare parts

## → DISINFECTION

72-73 OzoActive

## → LED-CURING

74-75 starlight uno

76-77 starlight pro

78-79 starlight s<sup>+</sup>/starlight sx5

80 spare parts

82-83 appropriate settings for PIEZOSURGERY® inserts

84-85 appropriate settings for ultrasound inserts

86-87 mectron online education

88-91 distributors

We reserve the right to make technical changes and changes in the scope of the supply. No liability can be accepted for any errors. The reference numbers of the units are meant for standard packagings; they may differ in some countries and according to different voltages.

# PIEZOSURGERY® *touch*

The future of bone surgery? Get in touch.



## PIEZOSURGERY® *touch*

### Equipment set:

- 1 PIEZOSURGERY® *touch* unit
- 1 LED-handpiece, complete with cord
- 1 kit of 8 tubes for pump
- 1 torque wrench K8
- 1 foot pedal
- power cord and suitcase

Ref. no.

05120061



## → TECHNICAL DATA

### **sizes (l x w x h in mm):**

300 x 250 x 95 mm

### **weight:**

3.2 kg

### **operating frequency:**

24.000 – 36.000 Hz

### **power types:**

6 functions  
IMPLANT  
CORTICAL  
CANCELLOUS  
SPECIAL  
PERIO  
ENDO

### **peristaltic pump capacity:**

adjustable on the touch screen:

ENDO / PERIO - 7 levels: from 0 to 6  
(from 0 to approximately 75 ml/min)  
SPECIAL / CANCELLOUS / CORTICAL /  
IMPLANT - 6 capacity levels: from 1 to 6  
(from 8 to approximately 75 ml/min)

### **power supply voltage:**

digital 100-240 VAC, 50/60 Hz

**max. power absorbed:** 120 VA

### **warranty:**

2 years, handpiece 1 year

## PIEZOSURGERY® *white*

### Equipment set:

1 PIEZOSURGERY® *white* unit

1 kit of 8 tubes for pump

1 torque wrench K8

1 foot pedal

power cord and suitcase  
(handpiece and inserts not included)

Ref. no. 05120063

### LED-HANDPIECE COMPLETE WITH CORD AND TRAY

- swivel-type LED light can be directed at the insert edge
- handpiece and handpiece cord (including the irrigation line) are fully sterilizable all at once

Ref. no. 03120134



### HANDPIECE COMPLETE WITH CORD AND TRAY

handpiece and handpiece cord (including the irrigation line) are fully sterilizable all at once

Ref. no. 03120217



# PIEZOSURGERY® *white*

## Get started in bone surgery!







## → TECHNICAL DATA

**sizes (l x w x h in mm):**

300 x 250 x 95 mm

**weight:**

3.2 kg

**operating frequency:**

24.000 – 36.000 Hz

**power types:**

6 functions

IMPLANT

CORTICAL

CANCELLOUS

SPECIAL

PERIO

ENDO

**peristaltic pump capacity:**

adjustable on the touch screen:

ENDO / PERIO - 7 levels: from 0 to 6

(from 0 to approximately 75 ml/min)

SPECIAL / CANCELLOUS / CORTICAL /

PERIO - 6 capacity levels: from 1 to 6

(from 8 to approximately 75 ml/min)

**power supply voltage:**

digital 100-240 VAC, 50/60 Hz

**max. power absorbed:** 120 VA

**warranty:**

2 years, handpiece 1 year

ENZYMEC is intended for use in cleaning of dental and general surgical instruments, and can be used during pre-cleaning in automatic equipments, ultrasonic cleaners and washer disinfectors according to the instructions and recommendations of the equipment manufacturer.

#### → ENZYMEC

| Product                          | Ref. no. |
|----------------------------------|----------|
| Enzymec bottle (1 liter)         | 03900001 |
| Enzymec 4-bottle pack (4 liters) | 03900002 |

# → ENZYMEC

## Enzymatic detergent

- enzymatic solution for efficient removal of organic residues
- specifically dedicated to the “clean” function of all PIEZOSURGERY® devices
- with practical dosage cap



# → UNIVERSAL CART

For daily use

**universal cart** – maximum mobility for all your table-top devices. All supply lines are hidden in the easily accessible cable room. The power supply is brought at the tray's height, for an easy connection to your devices.



## → SIZES:

**total (l x w x h in mm):**

480 x 460 x 850 mm

**upper tray:**

480 x 370 in 770 mm height

**middle tray:**

420 x 310 in 440 mm height

**lower tray:**

480 x 370 in 65 mm height

**weight:**

29 kg

## → UNIVERSAL CART

Product

Ref. no.

**universal cart**

03540011

**universal cart with  
water & air connections**

03540012



surgical tray



torque wrench K8



irrigation set



PIEZOSURGERY® insert box

# SPARE PARTS

High quality spare parts for professionals



DPS adaptor



PIEZOSURGERY®  
insert adaptor

## ADAPTORS FOR THERMODISINFECTION

Product

Ref. no.

**DPS Adaptor** (for PIEZOSURGERY® dental irrigation line) \*

04610006

**spare filter for DPS Adaptor** \*

04590006

**PIEZOSURGERY® insert Adaptor** \*

04610010

## PIEZOSURGERY® COMMON SPARE PARTS

**surgical tray**

04610004

**dynamometric torque wrench K8**

02900109

**dynamometric torque wrench K7** (for tipholder DB2)

02900081

**irrigation set** (50 pieces)

03230006

**PIEZOSURGERY® insert box**

02900111

**insert tray**

02900037

**insert box with bone expander tray**

02900107

**insert box with implant insert tray**

02900108

**bag support rod**

01380002

\* not disposable

## PIEZOSURGERY® *touch* AND PIEZOSURGERY® *white*

### Product

**LED-handpiece complete with cord** (PS *touch* and PS *white*)

**PIEZOSURGERY® *white* handpiece complete with cord**

**handpiece LED-cone**

**handpiece metal cone**

**handpiece holder**

**O-ring set for handpiece holder**

**mobile handpiece holder**

**kit of 8 tubes for pump**

**bone grafting kit** (liquid flow rate reducer, equipped with 8 silicone tubes and 1 tube-tube connector)

**protection foils**

(pack of 10 pcs. - each in a single sterile packaging)

**footswitch** (PS *touch*)

**footswitch** (PS *white*)

**PIEZOSURGERY® *touch* and *white* suitcase**

Ref. no.

03120134

03120217

03020158

03020171

02900105

03020160

03110015

02900104

02900106

03230012

02900102

02900114

04440006

## PIEZOSURGERY® I + II + 3

**handpiece complete with cord**

**handpiece holder** (PSII + 3)

**kit of 8 tubes for pump** (PSI)

**kit of 8 tubes for pump** (PSII + 3,

for handpiece cords with metal connector)

**kit of 8 tubes for pump** (PSII + 3,

for handpiece cords with plastic connector)

**peristaltic pump rolling** (PSII + 3)

**tube-tube connector**

**footswitch** (PSI + II)

**footswitch** (PS3)

**PIEZOSURGERY® 3 suitcase**

03120023

02900040

02900022

02900039

02900104

03210006

00680027

02900020

02900078

04440007

kit of 8 tubes for pump

protection foils

tube-tube  
connector

LED-handpiece


mobile  
handpiece  
holder

handpiece complete  
with cord



# INSERT SETS

## Set of inserts for different applications



## → BASIC SET



## → OSTEOTOMY SET



### → INSERT SETS

Product

Ref. no.

#### **basic set**

01520001

equipped with inserts

OT2, OT7, OP1, OP3 and EX1

description of the single inserts can  
be found on page 26, 27, 30 and 34

#### **osteotomy set**

01520010

equipped with inserts

OT7, OT7S-4, OT7S-3, OT8R and OT8L

description of the single inserts can  
be found on page 27

## → INSERT SETS

Product

Ref. no.

### **sinus lift lateral set\***

01520023

equipped with inserts

SLC, SLO-H, SLS, SLE1 and SLE2

description of the single inserts can  
be found on page 28, 31 and 32

### **piezo lift set**

01520024

equipped with inserts

PL1, PL2 and PL3

description of the single inserts can  
be found on page 29

\* developed in collaboration with Prof. Tomaso Vercellotti  
and Dr. Philippe Russe

## → SINUS LIFT LATERAL SET



## → PIEZO LIFT SET





## → SINUS LIFT SET



## → IMPLANT PREP SET PRO



### → INSERT SETS

| Product                                                                                                                                                                                         | Ref. no. |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>sinus lift set</b>                                                                                                                                                                           | 01520002 |
| equipped with inserts<br>OT1, OT5, EL1, EL2 and EL3<br>description of the single inserts can<br>be found on page 26 and 32                                                                      |          |
| <b>implant prep set pro</b>                                                                                                                                                                     | 01520011 |
| equipped with inserts<br>IM1S, IM2A, IM2P, IM3A, IM3P, OT4, IM4A,<br>IM4P, P2-3, P3-4, 3 PINS IM1S, 3 PINS 2-2,4<br>description of the single inserts can<br>be found on page 23, 24, 25 and 26 |          |

## → INSERT SETS

Product

Ref. no.

### **extraction set**

equipped with inserts

EX1, EX2, EX3, PS2 and PS6

description of the single inserts can  
be found on page 34 and 35

01520003

### **explantation set**

equipped with inserts

EXP3-R, EXP3-L, EXP4-R, EXP4-L

description of the single inserts can  
be found on page 34

01520021

## → EXTRACTION SET



## → EXPLANTATION SET



## → RESECTIVE PERIO SET



## → RETRO SURGICAL SET



### → INSERT SETS

| Product                                                                 | Ref. no. |
|-------------------------------------------------------------------------|----------|
| <b>resective perio set *</b>                                            | 01520022 |
| equipped with inserts<br>OT13, OT14, OP5A, OP8** and OP9                |          |
| description of the single inserts can<br>be found on page 28, 30 and 31 |          |
| <b>retro surgical set</b>                                               | 01520005 |
| equipped with inserts<br>OP7, EN1, EN3, EN5R, EN5L                      |          |
| description of the single inserts can<br>be found on page 31 and 33     |          |

\* developed in collaboration with Professor Leonardo Trombelli  
from the University of Ferrara

\*\* only available for PIEZOSURGERY® 3, PIEZOSURGERY® touch  
and PIEZOSURGERY white

# .....> INDICATIONS

Surgical choice: all

→ SINUS LIFT  
CRESTAL APPROACH



→ SINUS LIFT  
LATERAL APPROACH



→ IMPLANT SITE  
PREPARATION



→ RIDGE EXPANSION



→ PERIOSTEUM  
PREPARATION



→ EXTRACTIONS



→ PIEZO-LIFT

→ PL1 → PL2 → PL3

→ STANDARD

→ OPTIONAL

→ SLC → SLO-H → SLS → SLE1 → SLE2 → OP3 → OT1 → OT1A → OT5 → OT5A → OT5B → EL1 → EL2 → EL3

→ STANDARD

→ IM1S → IM2A → IM3A → IM4A → IM2P → IM3P → IM4P → OT4 → P2-3 → P3-4

→ OPTIONAL

→ IM1 AL → IM2A-15 → IM2.8A → IM3A-15 → IM3.4A → IM2P-15 → IM2.8P → IM3P-15 → IM3.4P

→ STANDARD

→ OPTIONAL

→ OT7 → OT4 → OP5 → OT2 → OT7A → OT7S-4 → OT7S-3 → OT12 → OT12S → OT7-20

→ STANDARD

→ PR1 → PR2

→ STANDARD

→ EX1 → EX2 → EX3 → PS2

→ STANDARD

→ EXP3-R → EXP3-L → EXP4-R → EXP4-L

→ STANDARD

→ OPTIONAL

→ OT7 → OP5 → OT8R → OT8L → OT6 → OT7A → OT7S-4 → OT7S-3 → OT12 → OT12S → OT7-20

→ STANDARD → OPTIONAL

→ OP3 → OP1 → OP2 → OP3A

→ STANDARD

→ OPTIONAL

→ OP7 → PS2 → EN1 → EN2 → EN3 → EN4 → EN5R → EN5L → EN6R → EN6L → OP3

→ STANDARD

→ OT1 → OT5

→ STANDARD → OPTIONAL

OT7S-4 → OT7S-3 → OT7 → OT7A

→ STANDARD

→ OP5A → OP8\* → OP9 → OT13 → OT14 → OP5 → OP3 → PS2 → PP1 → ICP + IC1

→ OPTIONAL

→ PS1 → PS6 → PP10 → PP11 → PP12 → OP2 → OP3A → OP4 → OP6 → OP6A

→ STANDARD

→ DB2 → CROWN PREP TIPS

\* only available for PIEZOSURGERY® 3, PIEZOSURGERY® touch  
and PIEZOSURGERY white



→ EXPLANTATION



→ BONE BLOCK GRAFTING



→ BONE CHIP GRAFTING/  
BONE MODELING



→ ENDODONTICS



→ OSTEOTOMY  
CLOSE TO NERVES



→ CORTICOTOMY  
TECHNIQUE



→ PERIODONTAL SURGERY  
AND IMPLANT CLEANING



→ PROSTHETICS

# → TRAYS

Do you like boxing? Don't miss them then.

## → PIEZOSURGERY® INSERT BOX

- fully autoclavable
- perfect storage of all PIEZOSURGERY® inserts
- light, easy to use, compact and roomy

Ref. no. 02900111

## → INSERT TRAY

- stainless steel tray with depth markings
- for storage

Ref. no. 02900037



PIEZOSURGERY® insert box

# → INSERTS

For every clinical needs, the proper insert.

## → IMPLANT PREP



IM1S



PIN IM1S



PIN 2-2.4



PIN 2-3



IM1-AL



PIN IM1-AL



PIN 2-2.4-AL



PIN 2-3-AL

### → IMPLANT PREP INSERTS


| Product                                                                                                                       | Ref. no. |
|-------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>insert IM1S</b><br>initial pilot osteotomy<br>diamond grain size 30 µm                                                     | 03510019 |
| <b>PIN IM1S</b>                                                                                                               | 03740007 |
| <b>PIN 2-2.4</b>                                                                                                              | 03740003 |
| <b>PIN 2-3</b><br>alignment pins, dedicated to<br>IM1S, IM2 (Ø 2), IM3 (Ø 3), OT4 (Ø 2.4) , IM1 SP                            | 03740001 |
| <b>insert IM1 AL</b><br>initial pilot osteotomy<br>diamond grain size 30 µm                                                   | 03510010 |
| <b>PIN IM1 AL</b>                                                                                                             | 03740004 |
| <b>PIN 2-2.4 AL</b>                                                                                                           | 03740005 |
| <b>PIN 2-3 AL</b><br>30 mm alignment pins, for<br>extraction sites dedicated to<br>M1 AL, IM2A (Ø 2), OT4 (Ø 2.4), IM3A (Ø 3) | 03740006 |

## → IMPLANT PREP INSERTS

| Product                                          | Ref. no. |
|--------------------------------------------------|----------|
| <b>insert IM2A</b>                               | 03510001 |
| <b>insert IM2A-15</b>                            | 03510015 |
| <b>insert IM2P</b>                               | 03510002 |
| <b>insert IM2P-15</b>                            | 03510016 |
| pilot osteotomy in anterior and posterior region |          |
| <b>insert IM2.8A</b>                             | 03510025 |
| <b>insert IM2.8P</b>                             | 03510026 |
| <b>insert IM3A</b>                               | 03510003 |
| <b>insert IM3A-15</b>                            | 03510017 |



to enlarge or to finalize the implant site preparation;  
insert with double irrigation to avoid overheating

## → IMPLANT PREP





## → IMPLANT PREP



### → IMPLANT PREP INSERTS

| Product                                                                                                       | Ref. no. |
|---------------------------------------------------------------------------------------------------------------|----------|
| <b>insert IM3P</b>                                                                                            | 03510004 |
| <b>insert IM3P-15</b>                                                                                         | 03510018 |
| <b>insert IM3.4A</b>                                                                                          | 03510027 |
| <b>insert IM3.4P</b>                                                                                          | 03510028 |
| to enlarge or to finalize the implant site preparation;<br>insert with double irrigation to avoid overheating |          |
| <b>insert IM4A</b>                                                                                            | 03510005 |
| <b>insert IM4P</b>                                                                                            | 03510006 |
| to finalize the implant site preparation;<br>insert with double irrigation to avoid overheating               |          |
| <b>insert P2-3</b>                                                                                            | 03510020 |
| <b>insert P3-4</b>                                                                                            | 03510021 |
| to optimize concentricity of implant site preparation<br>between Ø 2 and Ø 3 mm<br>diamond grain size 150 µm  |          |

## → OSTEOTOMY INSERTS

| Product            | Ref. no. |
|--------------------|----------|
| <b>insert OT1</b>  | 03370001 |
| <b>insert OT1A</b> | 03370008 |

to finalize the osteotomy in proximity of soft tissues  
(for example: sinus membrane, vessel, alveolar nerve)

OT1 diamond grain size 91 µm

OT1A diamond grain size 126 µm

|                                                                                                                                                         |          |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>insert OT2</b>                                                                                                                                       | 03370002 |
| osteotome of great precision in anatomically thin<br>structures (for example, ridge expansion, interdental<br>corticotomies, non traumatic nasal spine) |          |

|                                                                                                                                               |          |
|-----------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>insert OT4</b>                                                                                                                             | 03370004 |
| to correct pilot osteotomy axis, to finalize the implant<br>site preparation close to the alveolar nerve; sinus crestal<br>approach technique |          |
| diamond grain size 151 µm                                                                                                                     |          |

|                    |          |
|--------------------|----------|
| <b>insert OT5</b>  | 03370005 |
| <b>insert OT5A</b> | 03370009 |
| <b>insert OT5B</b> | 03370010 |

non-traumatic, to finalize the osteotomy or osteoplasty  
on thin bone and/or near delicate anatomic structures

OT5 Ø 1,7 - diamond grain size 91 µm

OT5A Ø 1,9 - diamond grain size 91 µm

OT5B Ø 1,7 - diamond grain size 126 µm

## → OSTEOTOMY



# → OSTEOTOMY



OT6

5 teeth,  
thickness 0.75 mm



OT7A

5 teeth,  
thickness 0.75 mm



OT7

5 teeth,  
thickness 0.55 mm



OT7-20

5 teeth,  
thickness 0.6 mm



OT7S-4

4 teeth,  
thickness 0.35 mm



OT7S-3

3 teeth,  
thickness 0.35 mm



OT8R

thickness 0.6 mm



OT8L

thickness 0.6 mm

## → OSTEOTOMY INSERTS

| Product                                                                                                                                         | Ref. no. |
|-------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>insert OT6</b>                                                                                                                               | 03370006 |
| osteotomy of large bone sections during maxillofacial surgery                                                                                   |          |
| <b>insert OT7</b>                                                                                                                               | 03370007 |
| <b>insert OT7A</b>                                                                                                                              | 03370011 |
| all the osteotomy technique in maxilla and mandible                                                                                             |          |
| – ridge expansion                                                                                                                               |          |
| – corticotomy technique                                                                                                                         |          |
| – bone block grafting                                                                                                                           |          |
| <b>insert OT7-20*</b>                                                                                                                           | 03370017 |
| ridge expansion technique, bone block grafting (from chin/mandible ramus), Le Fort I osteotomy, bilateral sagittal split osteotomy              |          |
| <b>insert OT7S-4</b>                                                                                                                            | 03370014 |
| <b>insert OT7S-3</b>                                                                                                                            | 03370015 |
| very thin osteotomy, corticotomy for orthodontic microsurgery technique, root separation in dental extraction technique and periodontal surgery |          |
| <b>insert OT8R</b>                                                                                                                              | 03370012 |
| <b>insert OT8L</b>                                                                                                                              | 03370013 |
| all the osteotomy technique in maxilla and mandible                                                                                             |          |
| – bone block grafting                                                                                                                           |          |

\* tighten up this insert only with torque wrench K8

## OSTEOTOMY INSERTS

Product

Ref. no.

### insert OT9

03370016

to finalize the osteotomy or osteoplasty near delicate anatomic structures; removal of the sinus basal cortex; insert with double irrigation to avoid overheating

### insert OT11

03370019

removal of the sinus basal cortex and elevation of the sinus membrane during crestal sinus lift procedures

### insert OT12

03370020

### insert OT12S

03370021

dedicated to the osteotomy techniques where the standard piezoelectric bone saws may have a difficult surgical approach

### insert OT13

03370022

### insert OT14

03370023

non-traumatic bone defect preparation on buccal and/or palatal cortical plates in resective periodontal surgery

OT13 Ø 1.7 - diamond grain size 150 µm

OT14 Ø 2.3 - diamond grain size 150 µm

### insert SLO-H

03370025

bone window antrostomy in sinus lift procedure by vestibular approach

## OSTEOTOMY



## → OSTEOTOMY



### → OSTEOTOMY INSERTS

| Product                                                                                                                | Ref. no.        |
|------------------------------------------------------------------------------------------------------------------------|-----------------|
| <b>insert PL1</b>                                                                                                      | <b>03370026</b> |
| sinus floor reaching and safe bony ring removal during crestal sinus lift procedure (PIEZO-LIFT)                       |                 |
| <b>insert PL2</b>                                                                                                      | <b>03370027</b> |
| sinus floor consumption and initial membrane elevation during crestal sinus lift procedure (PIEZO-LIFT)                |                 |
| <b>insert PL3</b>                                                                                                      | <b>03370028</b> |
| removal of the sinus basal cortex and elevation of the sinus membrane during crestal sinus lift procedure (PIEZO-LIFT) |                 |

## OSTEOPLASTY

### OSTEOPLASTY INSERTS

| Product                                                                                                                                                  | Ref. no. |
|----------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>insert OP1</b>                                                                                                                                        | 03380001 |
| bone remodelling and harvesting of bone chips                                                                                                            |          |
| <b>insert OP2</b>                                                                                                                                        | 03380002 |
| bone osteotomy, remodelling and removal of inflammatory tissue                                                                                           |          |
| <b>insert OP3</b>                                                                                                                                        | 03380003 |
| <b>insert OP3A</b>                                                                                                                                       | 03380011 |
| periodontal osteotomy, crown lengthening, bone chips harvesting, inflammatory tissue removal (cyst, etc.)                                                |          |
| <b>insert OP4</b>                                                                                                                                        | 03380004 |
| interproximal osteoplasty and root planing<br>diamond grain size 126 µm                                                                                  |          |
| <b>insert OP5</b>                                                                                                                                        | 03380005 |
| <b>insert OP5A</b>                                                                                                                                       | 03380012 |
| root debridement and root planing during resective and regenerative periodontal surgery<br>OP5 diamond grain size 30 µm<br>OP5A diamond grain size 90 µm |          |
| <b>insert OP6</b>                                                                                                                                        | 03380006 |
| root preparation in periodontal surgery<br>OP6 diamond grain size 126 µm                                                                                 |          |



## → OSTEOPLASTY



### → OSTEOPLASTY INSERTS

| Product                                                                     | Ref. no. |
|-----------------------------------------------------------------------------|----------|
| <b>insert OP7</b>                                                           | 03380007 |
| peri-apical maxillary bone osteotomy access,<br>inflammatory tissue removal |          |
| <b>insert OP8*</b>                                                          | 03380008 |
| <b>insert OP9</b>                                                           | 03380009 |
| Interproximal osteoplasty without damaging<br>adjacent root surfaces        |          |
| OP8* cutting edge from 0.7 - 1.3 mm                                         |          |
| OP9 cutting edge from 1 - 2 mm                                              |          |
| <b>insert SLC</b>                                                           | 03380013 |
| bone consuming during lateral sinus lift procedures                         |          |

\* only available for PIEZOSURGERY® 3, PIEZOSURGERY® touch  
and PIEZOSURGERY white

### → ELEVATION INSERTS

| Product                                                                   | Ref. no. |
|---------------------------------------------------------------------------|----------|
| <b>insert EL1</b>                                                         | 03390001 |
| separation of the sinus membrane,<br>2 mm around the frame of bony window |          |
| <b>insert EL2</b>                                                         | 03390002 |
| <b>insert EL3</b>                                                         | 03390003 |
| separation of the sinus membrane in internal zones                        |          |
| <b>insert SLS</b>                                                         | 03390006 |
| separation of the sinus membrane<br>during lateral sinus lift procedures  |          |
| <b>insert SLE1</b>                                                        | 03390007 |
| elevation of the sinus membrane from the sinus floor                      |          |
| <b>insert SLE2</b>                                                        | 03390008 |
| elevation of the sinus membrane<br>from the palatal sinus wall            |          |

### → PERIOSTEUM INSERTS

|                           |          |
|---------------------------|----------|
| <b>insert PR1</b>         | 03390004 |
| <b>insert PR2</b>         | 03390005 |
| subperiosteal preparation |          |

## → ELEVATION



## → PERIOSTEUM





## → ENDO RETRO



### → ENDO RETRO INSERTS

| Product                                                                                                                                         | Ref. no. |
|-------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <b>insert EN1</b><br>diamond-coated endo apical debrider 3 mm,<br>apical root debridement<br>diamond grain size 30 µm                           | 02170001 |
| <b>insert EN2</b><br>smooth endo apical debrider 3 mm,<br>gentle apical root debridement                                                        | 02170002 |
| <b>insert EN3</b><br>diamond-coated endo apical debrider 2.2 mm, apical<br>root debridement<br>diamond grain size 30 µm                         | 02170005 |
| <b>insert EN4</b><br>smooth endo apical debrider 2.2 mm,<br>gentle apical root debridement                                                      | 02170006 |
| <b>insert EN5R</b>                                                                                                                              | 02170007 |
| <b>insert EN5L</b><br>right and left angled, diamond-coated endo apical<br>debrider 2.2 mm, apical root debridement<br>diamond grain size 30 µm | 02170008 |
| <b>insert EN6R</b>                                                                                                                              | 02170009 |
| <b>insert EN6L</b><br>right and left angled, smooth endo apical<br>debrider 2.2 mm, apical root debridement                                     | 02170010 |

## → EXTRACTION INSERTS

| Product                                            | Ref. no. |
|----------------------------------------------------|----------|
| <b>insert EX1</b>                                  | 03400001 |
| to cut off the ankylosis, root fraction techniques |          |
| <b>insert EX2</b>                                  | 03400002 |
| analogue to EX1 in the posterior regions           |          |
| <b>insert EX3</b>                                  | 03400003 |
| analogue to EX1 in the posterior regions           |          |

## → EXPLANTATION INSERTS

|                      |          |
|----------------------|----------|
| <b>insert EXP3-R</b> | 03400004 |
| <b>insert EXP3-L</b> | 03400005 |
| <b>insert EXP4-R</b> | 03400006 |
| <b>insert EXP4-L</b> | 03400007 |
| explantations        |          |

## → EXTRACTION



## → EXPLANTATION





PS1



PS2



PS6



PP1



PP10



PP11



PP12

## → PERIODONTAL INSERTS

| Product                                                                                        | Ref. no. |
|------------------------------------------------------------------------------------------------|----------|
| <b>insert PS1</b><br>root scaling                                                              | 03180001 |
| <b>insert PS2</b><br>scaling and inflammatory tissue removal<br>fractured root apex extraction | 03180003 |
| <b>insert PS6</b><br>root scaling                                                              | 03180005 |
| <b>insert PP1</b><br>root planing                                                              | 03180002 |
| <b>insert PP10</b><br>gentle subgingival concretions removal                                   | 03180006 |
| <b>insert PP11</b><br>gentle subgingival concretions removal                                   | 03180007 |
| <b>insert PP12</b><br>gentle subgingival concretions removal                                   | 03180008 |

## → IMPLANT CLEANING



### → IMPLANT CLEANING TIPS

Product

Ref. no.

#### **Implant cleaning set P**

01520020

(consists of a tipholder ICP and 5 tips IC1)

#### **tipholder ICP**

03570003

#### **tip IC1 (5 pieces)**

02900112

efficient and soft cleaning of implants and restorations; implant cleaning tip IC1 can be easily tightened without any special tools and is fully autoclavable


## → CROWN PREP



key AB1



dynamometric  
wrench K7



tipholder DB2  
with crown prep tips

### → Available lengths and diameters of crown prep tips

Preparation with chamfer technique

→ Length 10 mm → Diamond coating

| diameter | D120     | D90     | D60     |
|----------|----------|---------|---------|
| 1.2 mm   |          | TA12D90 | TA12D60 |
| 1.4 mm   | TA14D120 | TA14D90 | TA14D60 |
| 1.6 mm   | TA16D120 | TA16D90 | TA16D60 |

Preparation with vertical preparation technique

→ Length 10 mm → Diamond coating

| diameter | D90     | D60     |
|----------|---------|---------|
| 1.2 mm   | TF12D90 | TF12D60 |
| 1.6 mm   | TF16D90 | TF16D60 |

### → CROWN PREP TIPHOLDER FOR PIEZOSURGERY® HANDPIECES

Product

Ref. no.

**crown prep set PS**  
equipped with DB2, AB1 and K7

02160010

**tipholder DB2**

03570002

### → CROWN PREP TIPHOLDER ACCESSORIES

**key AB1**

03580001

**dynamometric wrench K7**

02900081

### → CROWN PREP TIPS FOR CHAMFER T.

|                                |          |
|--------------------------------|----------|
| <b>crown prep tip TA12D60</b>  | 03590001 |
| <b>crown prep tip TA14D60</b>  | 03590002 |
| <b>crown prep tip TA16D60</b>  | 03590003 |
| <b>crown prep tip TA12D90</b>  | 03590004 |
| <b>crown prep tip TA14D90</b>  | 03590005 |
| <b>crown prep tip TA16D90</b>  | 03590006 |
| <b>crown prep tip TA14D120</b> | 03590007 |
| <b>crown prep tip TA16D120</b> | 03590008 |

### → CROWN PREP TIPS FOR VERTICAL T.

|                               |          |
|-------------------------------|----------|
| <b>crown prep tip TF12D60</b> | 03590010 |
| <b>crown prep tip TF16D60</b> | 03590011 |
| <b>crown prep tip TF12D90</b> | 03590012 |
| <b>crown prep tip TF16D90</b> | 03590013 |



## MULTIPIEZO

### Equipment set

1 MULTIPIEZO

1 LED ultrasonic handpiece slim

2 inserts: scaling insert S1-S and perio insert P3

2 bottles, 500 ml

2 safety bottle caps


Ref. no.

05060023

# MULTIPIEZO

Looking beautiful is looking really beautiful.


## RESTORATIVE MODE PULSE MODE



### Pulse Mode

Pulse Mode (level 6 of restorative function) transforms sine wave-shaped ultrasound oscillations to a new profile, characterized by specific power peaks, which results in an optimized performance for extractions and prosthetics.

## STANDARD MODE SOFT MODE



### SOFT MODE

- special modulation for gentle but still powerful insert movement
- maximum comfort for sensitive patients



## TECHNICAL DATA

**sizes (l x w x h in mm):**

320 x 230 x 145

**weight:**

2.4 kg

**operating frequency:**

24.000 – 36.000 Hz

**power types:**

3 functions, each with seven power levels, SOFT MODE

**liquid supply:**

via removable bottle

**liquid quantity:**

adjustable in 7 levels

**liquid flow rate:**

between 0 and 28 ml

**warranty:**

2 years, handpiece 1 year

# ..... MICROPIEZO

A new era in ultrasound devices.

The new advanced entry level ultrasound scaler with three functions:

- scaling and perio treatments
- endodontic treatments
- prosthetic treatments

Patients benefit from gentle ultrasound scaling in SOFT MODE.



## ..... MICROPIEZO

1 MICROPIEZO

1 LED ultrasonic handpiece slim

2 inserts: scaling insert S1-S and perio insert P3

1 liquid container, 500 ml

1 footpedal

Ref. no.

05080013







## → TECHNICAL DATA

**sizes (l x w x h in mm):**

161 x 214 x 171

**weight:**

0.94 kg

**operating frequency:**

24.000 – 36.000 Hz

**power types:**

3 functions (endo, perio-scaling, prosthetics),  
SOFT MODE, each with seven power levels

**liquid supply:**

via removable bottle

**liquid quantity:**

adjustable in 7 levels

**liquid flow rate:**

between 0 and 33 ml

## → TECHNICAL DATA

### operating frequency:

24.000 – 36.000 Hz

### power types:

adjustable via the dental chair

### power supply:

24 VAC, 50-60 Hz / 32 VDC

### water supply:

via the dental unit

### warranty:

2 years, handpiece 1 year

## → COMPACT PIEZO LED

### Equipment set

1 electronic module

1 LED ultrasonic handpiece slim

2 inserts: scaling insert S1-S and perio insert P3

Ref. no. 05090008

## → COMPACT PIEZO P2K

### Equipment set

1 electronic module

1 ultrasonic handpiece

2 inserts: scaling insert S1-S and perio insert P3

Ref. no. 05090003

# COMPACT PIEZO LED COMPACT PIEZO P2K

Compact efficiency or directly built-in.



Including feedback-system  
for constant power.

# SPARE PARTS

High quality spare parts for your ultrasound device.



DSC adaptor for ultrasound inserts



scaler insert adaptor

## ADAPTORS FOR THERMODISINFECTION

Product

**DSC Adaptor (for scaler handpiece irrigation line)\***  
**spare filter for DSC Adaptor\***  
**scaler insert Adaptor\***

Ref. no.

04610007  
 04590006  
 04610009

## COMMON SPARE PARTS

**dynamometric torque wrench K10**  
**dynamometric torque wrench K6**  
**dynamometric wrench K7 (for tipholder DB1)**  
**water filter**  
**insert tray**  
**universal cart**  
**universal cart with water & air connections**

02900137-001  
 02900074  
 02900081  
 00420004  
 02900072  
 03540011  
 03540012

\* not disposable



dynamometric torque wrench K10



dynamometric torque wrenches K6 and K7



water filter



safety  
bottle cap  
(blue)



liquid container (MICROPIEZO),  
incl. cap, 500ml



plastic bottle  
(MULTIPIEZO)  
with complete  
bottle cap

#### ➔ SPARE PARTS ULTRASOUND UNITS WITH LED HANDPIECE

Product

Ref. no.

**LED ultrasonic handpiece slim**

03120238

**kit light guide, metal cone and blue O-ring**

for LED ultrasonic handpiece slim

02900146

**LED ultrasonic handpiece**

03120142

**handpiece LED-cone**

03020158

**handpiece metal cone**

03020171

#### ➔ SPARE PARTS ULTRASOUND UNITS WITHOUT LED HANDPIECE

**ultrasonic handpiece (grey)**

03120016

**ultrasonic handpiece (blue)**

03120132

**handpiece metal cone**

00710034

#### ➔ BOTTLES AND BOTTLE CAPS FOR PIEZO SMART AND MULTIPIEZO RANGE

**plastic bottle (piezo smart), 300 ml,**

with blue complete bottle cap

03020066

**plastic bottle (piezo smart), 300 ml,**

with grey complete bottle cap

03020067

**plastic bottle (piezo smart), 300 ml, without cap**

03230007

**plastic bottle (MULTIPIEZO), 500 ml, without cap**

03150102

**complete bottle cap, blue**

03020081

**complete bottle cap, grey**

03020082

**safety bottle cap, grey**

01950015

**safety bottle cap, blue**

01950014

## → SPARE PARTS MULTIEPIEZO RANGE

| Product                                                     | Ref. no. |
|-------------------------------------------------------------|----------|
| <b>O-ring set for bottle holder</b> (6 pieces)              | 03020175 |
| <b>peristaltic pump</b> (white)                             | 03210001 |
| <b>peristaltic pump</b> (black - from serial no. 106005993) | 03210009 |
| <b>footswitch</b>                                           | 02900114 |

## → SPARE PARTS MICROPIEZO

|                                                                                         |          |
|-----------------------------------------------------------------------------------------|----------|
| <b>LED ultrasonic handpiece slim</b>                                                    | 03120238 |
| <b>kit light guide, metal cone and blue O-ring</b><br>for LED ultrasonic handpiece slim | 02900146 |
| <b>liquid container</b> (incl. cap, 500ml)                                              | 02900145 |
| <b>footpedal</b>                                                                        | 02900144 |
| <b>pump</b>                                                                             | 03210009 |
| <b>power adapter</b>                                                                    | 03160024 |
| <b>kit O-ring + spring for liquid container</b>                                         | 02900147 |

## → SPARE PARTS PIEZO SMART

|                                                |          |
|------------------------------------------------|----------|
| <b>O-ring set for bottle holder</b> (3 pieces) | 03020068 |
| <b>peristaltic pump</b>                        | 03210007 |
| <b>footswitch</b>                              | 02900055 |
| <b>power adapter piezo smart</b>               | 03020105 |

## → SPARE PARTS MICROPIEZO S

|                     |          |
|---------------------|----------|
| <b>footswitch</b>   | 02900001 |
| <b>water filter</b> | 00420004 |



LED ultrasonic handpiece slim



LED ultrasonic handpiece



ultrasonic handpiece (grey)



ultrasonic handpiece (blue)

# → INSERT SETS

An insert set for every indication.



→ K10 – protects your handpiece and inserts

- sterilisable dynamometric torque wrench
- supplied with inserts



## → SET MAINTENANCE



## → SET SCALING



### → INSERT SETS

Product

Ref. no.

#### **set maintenance**

02160035-001

equipped with inserts

S1, S2, S5, P2, P3 and

5 torque wrenches K10

description of the single inserts

can be found on page 52 and 53

#### **set scaling**

02160001-001

S1, S2, S3, S4, S5 and

5 torque wrenches K10

description of the single inserts

can be found on page 52

## → INSERT SETS

Product

### **set perio maintenance**

equipped with inserts  
P10, P11, P12, P13, P14 and  
5 torque wrenches K10  
description of the single inserts  
can be found on page 54

### **set perio treatment**

equipped with inserts  
P2, P3, P15, P16R, P16L and  
5 torque wrenches K10  
description of the single inserts  
can be found on page 53 and 54

Ref. no.

02160006-001

02160036-001

## → SET PERIO MAINTENANCE



## → SET PERIO TREATMENT





## → SET IMPLANT CLEANING



## → SET ENDO



### → INSERT SETS

| Product                           | Ref. no.     |
|-----------------------------------|--------------|
| <b>implant cleaning set S</b>     | 02160033-001 |
| equipped with inserts             |              |
| 1 insert ICS                      |              |
| 1 torque wrench K10               |              |
| 5 Implant cleaning tips IC1       |              |
| description of the single inserts |              |
| can be found on page 55           |              |

### **set endo** 02980020

|                                                        |
|--------------------------------------------------------|
| equipped with inserts                                  |
| E1 (file holder 120°)                                  |
| E2 (file holder 90°)                                   |
| 6 NiTi-files ISO 15, 27 mm, 6 NiTi-files ISO 20, 27 mm |
| 6 NiTi-files ISO 25, 27 mm, 6 NiTi-files ISO 15, 31 mm |
| 6 NiTi-files ISO 20, 31 mm, 6 NiTi-files ISO 25, 31 mm |
| 1 wrench K1                                            |
| 1 wrench for endo files                                |
| description of the single inserts                      |
| can be found on page 56                                |

## → SET ENDO REVISION



### → INSERT SETS

Product

Ref. no.

#### **set endo revision**

02160008-001

equipped with inserts

ER1, ER2, ER3, ER4, ER5 and

5 torque wrenches K10

description of the single inserts

can be found on page 57





→ every insert  
is supplied with its  
own torque wrench

## → SCALING INSERTS

Product

Ref. no.

### **insert S1**

02960001-001

universal curette with semicircular diameter for considerable tartar removal

### **insert S1-S**

02960009-001

slim universal curette with triangled slightly curved surface for removal of supragingival tartar, subgingival concretions and biofilm

### **insert S2**

02960002-001

universal curette with triangled, slightly curved surface, efficient in the interdental spaces and posterior surfaces

### **insert S3**

02960003-001

flat, with rounded edges for considerable supragingival tartar removal

### **insert S4**

02960004-001

universal curette with 45° angled, triangled and slightly curved surface for tartar removal on mesial and distal interdental surfaces in the posterior area

### **insert S5**

02960005-001

similar shape to S1 insert but longer and thinner for gentle supra- and subgingival tartar removal and for gingivitis

### **insert S6**


02960006-001

contra-angled universal curette with semicircular diameter, powerful insert (twice the power of S1 insert) for considerable tartar removal

# → INSERTS

For a wide range of applications

## → SCALING



## → PERIO UNIVERSAL



→ every insert  
is supplied with its  
own torque wrench

### → PERIO UNIVERSAL INSERTS

| Product                                                                      | Ref. no.     |
|------------------------------------------------------------------------------|--------------|
| <b>insert P1</b><br>long, straight working tip with a circular surface       | 02970001-001 |
| <b>insert P2</b><br>curved working tip with a circular surface               | 02970002-001 |
| <b>insert P3</b><br>slim insert smooth like a probe                          | 02970003-001 |
| <b>insert P4</b><br>extra long, straight working tip with a circular surface | 02970004-001 |



→ every insert  
is supplied with its  
own torque wrench

## → PERIO ANATOMIC INSERTS

| Product                                                                        | Ref. no.     |
|--------------------------------------------------------------------------------|--------------|
| <b>insert P10</b>                                                              | 03080001-001 |
| extra long, curved working tip with a circular surface                         |              |
| <b>insert P11</b>                                                              | 03080002-001 |
| right (inclination of 15°) curved insert with round tip                        |              |
| <b>insert P12</b>                                                              | 03080003-001 |
| left (inclination of 15°) curved insert with round tip                         |              |
| <b>insert P13</b>                                                              | 03080004-001 |
| right (inclination of 15°) curved insert with spherical (Ø = 0.8 mm) round tip |              |
| <b>insert P14</b>                                                              | 03080005-001 |
| left (inclination of 15°) curved insert with spherical (Ø = 0.8 mm) round tip  |              |
| <b>insert P15</b>                                                              | 03080006-001 |
| universal curette shape with round surface                                     |              |
| <b>insert P16R</b>                                                             | 03080007-001 |
| right angled periodontal curette                                               |              |
| <b>insert P16L</b>                                                             | 03080008-001 |
| left angled periodontal curette                                                |              |

## → PERIO ANATOMIC



## ➔ IMPLANT CLEANING



➔ implant cleaning tip IC1 can be easily tightened without any special tools and is fully autoclavable



➔ every insert is supplied with its own torque wrench

### ➔ IMPLANT CLEANING INSERTS

| Product                                                                               | Ref. no.     |
|---------------------------------------------------------------------------------------|--------------|
| <b>Implant cleaning set S</b><br>(consists of a tipholder ICS and 5 tips IC1)         | 02160033-001 |
| <b>tipholder ICS</b>                                                                  | 03570004-001 |
| <b>tip IC1 (5 pieces)</b><br>efficient and soft cleaning of implants and restorations | 02900112     |



→ every insert  
is supplied with its  
own torque wrench

## → ENDO INSERTS

Product

Ref. no.

### insert E1

02350001-001

120° file holder for endo files and inserts  
with a shaft of 0,8 mm

### insert E2

02350002-001

90° file holder for endo files and inserts  
with a shaft of 0,8 mm

## → NITI-FILES

highly flexible NiTi-files for mectron file holders,  
two lengths, packages of 12 pieces

12 x endo file, ISO 15, 31 mm 02730001

12 x endo file, ISO 20, 31 mm 02730002

12 x endo file, ISO 25, 31 mm 02730003

12 x endo file, ISO 15, 27 mm 02740001

12 x endo file, ISO 20, 27 mm 02740002

12 x endo file, ISO 25, 27 mm 02740003

## → SPARE PARTS

key K1 for file holders

E1 and E2

02270003

key for screw of file holders

E1 and E2

02270004


screw for file holder E1

01260005

screw for file holder E2

01260003

## → ENDO



ISO 15



ISO 20



ISO 25



ISO 15



ISO 20



ISO 25

length 27 mm

length 31 mm



## → ENDO REVISION



→ every insert  
is supplied with its  
own torque wrench

### → ENDO REVISION INSERTS

| Product                                                                                   | Ref. no.     |
|-------------------------------------------------------------------------------------------|--------------|
| <b>insert ER1</b>                                                                         | 03450001-001 |
| angled shaft with small, spherical (Ø 1.7 mm) and D30 diamond coated tip                  |              |
| <b>insert ER2</b>                                                                         | 03450002-001 |
| angled insert with 0.7 mm Ø, working length up to 10 mm, the last 5 mm D30 diamond coated |              |
| <b>insert ER3</b>                                                                         | 03450003-001 |
| angled insert with 0.6 mm Ø, working length up to 10 mm, the last 5 mm D30 diamond coated |              |
| <b>insert ER4</b>                                                                         | 03450004-001 |
| angled insert with 0.6 mm Ø, working length up to 20 mm, the last 5 mm D30 diamond coated |              |
| <b>insert ER5</b>                                                                         | 03450005-001 |
| angled insert with 0.5 mm Ø, working length up to 24 mm                                   |              |



→ every insert  
is supplied with its  
own torque wrench

## → ENDO RETRO INSERTS

Product

Ref. no.

### insert RD3

03050006-001

angled shaft, with 90° angled working tip  
and fine diamond coating (D30)

### insert RD4

03050007-001

angled shaft, with 90° right angled working tip  
and fine diamond coating (D30)

### insert RD5

03050008-001

angled shaft, with 90° left angled working tip  
and fine diamond coating (D30)



→ equipped with  
K6 torque wrench

## → EXTRACTION INSERTS

Product

Ref. no.

### insert ME1\*

03130001-001

### insert ME2\*


03130002-001

### insert ME3\*

03130003-001

The specially shaped inserts ME1, ME2 and ME3 – with a blade thickness of only 0.35 mm – perfectly fulfill the anatomical demands of non surgical extractions. They offer maximum intra-operative sensitivity, allowing a great preservation of alveolar bone and periodontal tissues. Furthermore, the physiological solution guarantees maximum visibility thanks to the cavitation effect which leads to a temporary hemostasis of blood vessels.

## → ENDO RETRO



RD3

RD4

RD5

## → EXTRACTION\*



ME1

ME2

ME3

→ EXCLUSIVE  
FOR PULSE  
MODE\*


## → CROWN PREP



key AB1



dynamometric wrench K7



tipholder DB1  
with crown prep tips

### → Available lengths and diameters of crown prep tips

Preparation with chamfer technique

→ Length 10 mm → Diamond coating

| diameter | D120     | D90     | D60     |
|----------|----------|---------|---------|
| 1.2 mm   |          | TA12D90 | TA12D60 |
| 1.4 mm   | TA14D120 | TA14D90 | TA14D60 |
| 1.6 mm   | TA16D120 | TA16D90 | TA16D60 |

Preparation with vertical preparation technique

→ Length 10 mm → Diamond coating

| diameter | D90     | D60     |
|----------|---------|---------|
| 1.2 mm   | TF12D90 | TF12D60 |
| 1.6 mm   | TF16D90 | TF16D60 |

### → CROWN PREP TIPHOLDER FOR ULTRASOUND SCALER HANDPIECES

| Product                                                       | Ref. no. |
|---------------------------------------------------------------|----------|
| <b>crown prep set scaler</b><br>equipped with DB1, AB1 and K7 | 02160009 |
| <b>tipholder DB1</b>                                          | 03570001 |

### → CROWN PREP TIPHOLDER ACCESSORIES

|                               |          |
|-------------------------------|----------|
| <b>key AB1</b>                | 03580001 |
| <b>dynamometric wrench K7</b> | 02900081 |

### → CROWN PREP TIPS FOR CHAMFER T.

|                                |          |
|--------------------------------|----------|
| <b>crown prep tip TA12D60</b>  | 03590001 |
| <b>crown prep tip TA14D60</b>  | 03590002 |
| <b>crown prep tip TA16D60</b>  | 03590003 |
| <b>crown prep tip TA12D90</b>  | 03590004 |
| <b>crown prep tip TA14D90</b>  | 03590005 |
| <b>crown prep tip TA16D90</b>  | 03590006 |
| <b>crown prep tip TA14D120</b> | 03590007 |
| <b>crown prep tip TA16D120</b> | 03590008 |

### → CROWN PREP TIPS FOR VERTICAL T.

|                               |          |
|-------------------------------|----------|
| <b>crown prep tip TF12D60</b> | 03590010 |
| <b>crown prep tip TF16D60</b> | 03590011 |
| <b>crown prep tip TF12D90</b> | 03590012 |
| <b>crown prep tip TF16D90</b> | 03590013 |

### → COMBI touch STANDARD

#### Equipment set

- 1 LED ultrasonic handpiece slim
- 2 inserts: scaling insert S1-S and perio insert P3
- 1 bottle (500ml)
- 1 grey safety bottle cap
- 1 air-polishing handpiece 120°
- 1 bottle of prophylaxis powder soft+ (250g)
- 1 bottle of prophylaxis powder sensitive+ - glycine powder (160g)

Ref. no.

05020015

### → COMBI touch PERIO

#### Equipment set

- 1 LED ultrasonic handpiece slim
- 2 inserts: scaling insert S1-S and perio insert P3
- 1 bottle (500ml)
- 1 grey safety bottle cap
- 1 air-polishing handpiece 120°
- 1 air-polishing handpiece PERIO
- 1 pack Subgingival Perio Tips (40 pcs)
- 1 wrench k9
- 1 bottle prophylaxis powder soft+ (250g)
- 1 bottle prophylaxis powder sensitive+ - glycine powder (160g)

Ref. no.


05020015

## → COMBI touch

Ultrasound and air-polishing in one unit.

- flexible supra- and subgingival air-polishing
- glass touch screen, easy to clean
- removable powder containers
- SOFT MODE for ultra-gentle scaling

→ NEW  
SLIM  
HANDPIECES



- widest range of ultrasound indications in the market



## → TECHNICAL DATA

**sizes (l x w x h in mm):** 410 x 260 x 145

**weight:** 4.8 kg

**water supply:**

working pressure from 1 to 6 bars

**air supply:**

input pressure between 4 and 8 bars

**ultrasound part:**

**operating frequency:**

24 KHz to 36 KHz

**power mode:**

"endo", "perio/scaler", "restorative", "SOFT MODE"

power levels from 1 to 5

**irrigation:**

tap water line – stepless adjustment

**peristaltic pump – adjustable on the**

**touch screen:** 7 flow levels from 1

(approx. 5 ml/min) to 7 (approx. 30 ml/min)

possibility to switch off the irrigation in the

power modes "endo", "perio/scaler" and

"restorative" at power levels from 1 to 4

**polishing part:**

**polisher function:**

selectable via touch screen: "Prophy" or "Perio"

**irrigation:**

continuous regulation via knob,

water heated via heating element.

**warranty:**

2 years, handpiece 1 year

### 7 different turbine adapters

- Bien-Air UNIFIX® and UNIFIX® L
- KaVo MULTIFLEX® and MULTIFLEX® lux
- SIRONA® R and B
- W&H® ROTO QUICK® and ROTO QUICK® lux
- NSK® MACHLITE® and PHATELUS®
- Midwest
- Borden

- UNIFIX® is a registered trademark of Bien-Air (BIEN-AIR®) Holding SA
- MULTIFLEX® is a registered trademark of Kaltenbach & Voigt (KaVo®) GmbH
- SIRONA® is a registered trademark of SIRONA DENTAL SYSTEMS GMBH
- ROTO QUICK® is a registered trademark of W&H® Deutschland GmbH
- MACHLITE® and PHATELUS® are registered trademarks of Nakanishi Corporation and Nakanishi Shika Kikai (NSK®) Ltd

### → STARJET BASIC

#### Equipment set

1 STARJET

1 spray head 120°

1 bottle prophylaxis powder soft+ (250g)


1 bottle prophylaxis powder sensitive+ - glycine powder (160g)

#### Product

|                         | Ref. no. |
|-------------------------|----------|
| STARJET basic (BienAir) | 05150074 |
| STARJET basic (KaVo)    | 05150075 |
| STARJET basic (Sirona)  | 05150076 |
| STARJET basic (W&H)     | 05150077 |
| STARJET basic (Midwest) | 05150078 |
| STARJET basic (Borden)  | 05150079 |
| STARJET basic (NSK)     | 05150080 |

# → STARJET

## The 2 face air-polisher.



**One device, two indications – easy switch from supra- to subgingival air-polishing with a PROPHY to PERIO turn.**

- maximum flexibility as it can be used with powders with different granulometries
- optimal access thanks to 3 available spray heads – 120°, PERIO and 90° optional
- permanent control of the powder level thanks to the transparent powder chamber
- easy powder filling thanks to lateral opening



## → TECHNICAL DATA

**sizes (l x w x h in mm):**

200 x 90 x 50

**water supply:**

0 – 3.0 bar

**air supply:**

2.0 – 3.0 bar

**warranty:**

1 year

## → STARJET

### Equipment set

1 STARJET

1 spray head 120°

1 spray head PERIO;

Subgingival Perio Tips (40 pcs)

1 wrench K9

1 bottle prophylaxis powder soft+ (250g)

1 bottle prophylaxis powder sensitive+ - glycine powder (160g)

Product

Ref. no.

STARJET (BienAir)

05150074

STARJET (KaVo)

05150075

STARJET (Sirona)

05150076

STARJET (W&H)

05150077

STARJET (Midwest)

05150078

STARJET (Borden)

05150079

STARJET (NSK)

05150080

## → PROPHYLAXIS POWDER

| Product                                                                       | Ref. no. |
|-------------------------------------------------------------------------------|----------|
| <b>Prophylaxis powder smooth</b><br>pack with four bottles<br>(250 g each)    | 03140010 |
| <b>Prophylaxis powder soft+</b><br>pack with four bottles<br>(250 g each)     | 03140027 |
| <b>Prophylaxis powder sensitive+</b><br>pack with two bottles<br>(160 g each) | 03140021 |

# → MECTRON PROPHYLAXIS POWDERS.

Higher performance for bright results.

## → PROPHYLAXIS POWDER SMOOTH \*

spherical  
particles gently  
remove supra-  
gingival biofilm  
while respecting  
the surrounding  
tissues



- Calcium carbonate-based
- Average particle size between 45 µm and 75 µm
- Neutral taste

## → PROPHYLAXIS POWDER SOFT+

gently and  
effectively  
removes  
supragingival  
biofilm and  
stains on the  
enamel surface



- Sodium bicarbonate-based
- Average particle size 40 µm
- Mint flavour

\* available until stock lasts





## → PROPHYLAXIS POWDER SENSITIVE<sup>+</sup>

gently removes  
supra- and sub-  
gingival biofilm



- Glycine-based
- Average particle size 25 µm
- Sweet taste

| INDICATIONS                                      | SMOOTH                                                                              | SOFT <sup>+</sup>                                                                   | SENSITIVE <sup>+</sup>                                                              |
|--------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Supragingival biofilm removal                    |  |  |  |
| Removal of stains and pigmentations              |                                                                                     |  |                                                                                     |
| Subgingival biofilm removal                      |                                                                                     |                                                                                     |  |
| Periodontal pocket detoxification                |                                                                                     |                                                                                     |  |
| Implant maintenance                              |                                                                                     |                                                                                     |  |
| Peri-implantitis treatment                       |                                                                                     |                                                                                     |  |
| Fixed orthodontic appliances                     |  |  |  |
| Cleaning of prosthesis and restoration materials |  |                                                                                     |  |
| Restoration maintenance                          |  |  |  |

#### → SIZES:

**total (l x w x h in mm):**

480 x 460 x 850 mm

**upper tray:**

480 x 370 in 770 mm height

**middle tray:**

420 x 310 in 440 mm height

**lower tray:**

480 x 370 in 65 mm height

**weight:**

29 kg

# → UNIVERSAL CART

## Let's move the devices!

The **universal cart** – maximum mobility for all your table-top devices. All supply lines are hidden in the easily accessible cable room. The power supply (and the water & air connection – where available) is brought at the tray's height, for an easy connection to your devices.



#### → UNIVERSAL CART

Product

Ref. no.

**universal cart**

03540011

**universal cart with  
water & air connections**

03540012





water filter



plastic bottle (scaler) –  
COMBI touch



LED ultrasonic handpiece slim



LED ultrasonic handpiece



dynamometric  
torque wrench K6



dynamometric  
torque wrench K10

# SPARE PARTS

## COMMON SPARE PARTS COMBI touch

Product

Ref. no.

dynamometric torque wrench K10

02900137-001

dynamometric torque wrench K6

02900074

dynamometric torque wrench K7

02900081

plastic bottle (scaler), 500 ml, without cap

03150102

complete bottle cap, blue

03020081

complete bottle cap, grey

03020082

safety bottle cap, grey

01950014

safety bottle cap, blue

01950015

Subgingival Perio Tip (40 pieces)

02900120

wrench K9

02900122

cleaning needles kit 0,4 x 57 mm and 0,8 x 60 mm

02900151

O-ring set COMBI touch

03020206

powder chamber PROPHY

03020194

powder chamber PERIO

03020200

powder chamber cap

03020195

water filter

00420004

water supply hose

02900119

air supply hose

02900012

suitcase COMBI touch

04440019

footswitch

02900114

## → SPARE PARTS COMBI *touch* WITH SLIM HANDPIECES

| Product                                                                         | Ref. no. |
|---------------------------------------------------------------------------------|----------|
| LED ultrasonic handpiece slim                                                   | 03120238 |
| kit light guide metal cone and blue O-ring<br>for LED ultrasonic handpiece slim | 02900146 |
| air-polishing handpiece 120°                                                    | 03120240 |
| air-polishing handpiece PERIO                                                   | 03120241 |
| air-polishing handpiece 90°                                                     | 03120239 |
| cleaning needles kit 0,4 x 57 mm and 0,8 x 60 mm                                | 02900151 |

## → SPARE PARTS COMBI *touch* WITH STANDARD HANDPIECES up to serial no. 205010618

|                                                  |          |
|--------------------------------------------------|----------|
| LED ultrasonic handpiece                         | 03120224 |
| handpiece LED-cone                               | 03020205 |
| handpiece metal cone                             | 03020171 |
| 90° spray head                                   | 03340007 |
| 120° spray head                                  | 03340008 |
| PERIO spray head                                 | 03340009 |
| O-ring set COMBI <i>touch</i>                    | 03020206 |
| cleaning needles kit 0,4 x 57 mm and 0,8 x 57 mm | 02900149 |

## → COMMON SPARE PARTS COMBI AND TURBODENT

|                                |              |
|--------------------------------|--------------|
| dynamometric torque wrench K10 | 02900137-001 |
| torque wrench K6               | 02900074     |
| cleaning needle                | 02900003     |
| O-ring for powder chamber cap  | 00410018     |
| footswitch                     | 02900001     |
| water filter                   | 00420004     |
| O-ring water coupling          | 00410024     |



air-polishing handpiece 120°



air-polishing handpiece PERIO



air-polishing handpiece 90°



90° spray head



120° spray head



PERIO spray head



ultrasonic handpiece (blue)



ultrasonic handpiece (grey)



spray head set turbodent p/s



spray head, to easyjet serial no. 080XXXXY



spray head, from easyjet serial no. 080XXXXY



O-ring Set

## ADAPTORS FOR THERMODISINFECTION

Product

Ref. no.

**DSC Adaptor** (for scaler handpiece irrigation line) \*  
**spare filter for DSC Adaptor** \*  
**scaler insert Adaptor** \*

04610007  
 04590006  
 04610009

## SPARE PARTS COMBI, up to serial no. 046XXXXY / TURBODENT, up to serial no. 047XXXXY

**ultrasonic handpiece (grey)**  
**ultrasonic handpiece cord (grey)**  
**spray head set turbodent p/s**  
**air polishing handpiece cord (grey)**  
**O-ring spray head, small**  
**O-ring spray head, large**

03120016  
 02830022  
 03340003  
 03120009  
 00410013  
 00410014

## SPARE PARTS COMBI, from serial no. 203XXXXY / TURBODENT, from serial no. 202XXXXY

**ultrasonic handpiece (blue)**  
**ultrasonic handpiece cord (white)**  
**spray head**  
**spray head set**  
 (spray head, silicon protector, long cleaning needle and short cleaning needle)  
**air polishing handpiece cord**  
**silicon protector**

03120132  
 02830042  
 03340006  
 03020121  
 03120133  
 03150033

## SPARE PARTS STARJET

**90° spray head**  
**120° spray head**  
**PERIO spray head**  
**Subgingival Perio Tip (40 pieces)**

03340015  
 03340014  
 03340013  
 02900120

## SPARE PARTS STARJET

| Product                                                                    | Ref. no.     |
|----------------------------------------------------------------------------|--------------|
| <b>wrench K9</b>                                                           | 02900122     |
| <b>O-ring set STARJET (2 black + blue)</b>                                 | 02900134-001 |
| <b>set powder chamber cap + 1 O-ring blue</b>                              | 02900133-001 |
| <b>cleaning kit STARJET</b>                                                | 02900136-001 |
| (cleaning needle spray head, cleaning needle STARJET body, cleaning brush) |              |

## COMMON SPARE PARTS EASYJET PRO

|                                         |          |
|-----------------------------------------|----------|
| <b>cleaning needle (long)</b>           | 02900050 |
| <b>silicon protector for spray head</b> | 03150033 |
| <b>O-ring Borden turbine adapter</b>    | 01970005 |
| <b>O-ring Midwest turbine adapter</b>   | 01970013 |

## SPARE PARTS EASYJET PRO, up to serial no. 080XXXXXX

|                                                                                            |          |
|--------------------------------------------------------------------------------------------|----------|
| <b>spray head set easyjet pro</b>                                                          | 03020061 |
| (spray head, silicon protector, cleaning needle, small and large O-ring)                   |          |
| <b>O-ring set easyjet pro</b>                                                              | 02900051 |
| (1 x small and 1 x large O-ring for the spray head, 1 x O-ring for the powder chamber cap) |          |
| <b>powder chamber cap, transparent</b>                                                     | 03020059 |
| <b>powder chamber cap, blue</b>                                                            | 03020060 |

## SPARE PARTS EASYJET PRO (from serial no. 201XXXXXX) AND EASYJET PERIO

|                                                                                      |          |
|--------------------------------------------------------------------------------------|----------|
| <b>spray head</b>                                                                    | 03340006 |
| <b>spray head set</b>                                                                | 03020121 |
| (spray head, silicon protector, long cleaning needle and short cleaning needle)      |          |
| <b>powder chamber cap, blue</b>                                                      | 03020123 |
| <b>easyjet O-ring set</b>                                                            | 02900087 |
| (1 small and 1 large O-ring for the spray head, 1 O-ring for the powder chamber cap) |          |



PERIO spray head



120° spray head



90° spray head



powder chamber cap

# .....> OZOACTIVE

## The Ozone's benefits for Dentistry.

The **OzoActive** device is specially designed to generate ozonated water at a concentration range that is very safe for both patients and operators.

- minimize the risk of aerosol for a safer clinical practice during the daily activity
- Ozone gas has a high oxidation potential when used as an antimicrobial agent against bacteria and viruses
- effective and faster sterilization and healing of the surgical site
- postoperative pain reduction thanks to its anti-inflammatory action
- perio pockets cleansing and disinfection and periimplantitis management
- hygiene treatment of fixed and mobile prosthesis and orthodontic patients

### .....> OZOACTIVE

#### Equipment set

- 1 device body
- 1 multivoltage power adapter
- 1 footswitch
- Instructions for use
- 1 setting sheet
- 3 disposable irrigation set for PIEZOSURGERY®
- 1 irrigation set for prophylaxis ultrasound units
- 1 pump cover for prophylaxis ultrasound units

Ref. no.

05270001







OzoActive generates ozone gas for distribution into the irrigation lines of mectron PIEZOSURGERY® and ultrasound units.



## → TECHNICAL DATA

**sizes (l x w x h in mm):**

150 x 140 x 150

**weight:**

1.6 kg

**power types:**

2 programs: Oral Surgery and Prophylaxis

**liquid and ozone quantity:**

adjustable from 0.011 mg to 0.079 mg/liter

synchronized use: level 3 on PIEZOSURGERY®

corresponds to level 3 on OzoActive

**power supply:**

100-240 V~ 50/60 Hz

**max. power consumption:**

36 W

**warranty:**

2 years

## → OZOACTIVE ACCESSORIES

Product

Ref. no.

**Aqualab OzoActive -  
Prophylaxis Irrigation Kit -  
reusable (1 pc.)**

02900157

**PIEZOSURGERY®  
Irrigation Kit - disposable  
(1 pack of 30 pcs.)**

03230018

**Aqualab OzoActive footpedal**

02900156

# → STARLIGHT UNO

Colorful and even easiest handling

LED curing light available in 6 different colours.  
One touch handling, outstanding curing performance, exchangeable battery ... persuasive features, persuasive price!

## → STARLIGHT UNO

1 handpiece, 1 charging station

Color

- ☐ white
- ☐ grey
- ☐ turquoise
- ☐ yellow
- ☐ light blue
- ☐ violet

Ref. no.

- 05100092
- 05100093
- 05100094
- 05100095
- 05100096
- 05100097





## → TECHNICAL DATA

**sizes (l x Ø in mm):**

190 x 21

**weight:**

77 g

**light source:**

high power LED

**wavelength:**

dominant wavelength 440 – 465 nm

**light intensity\*:**

> 1.500 mW/cm<sup>2</sup>

**average life span LED:**

1.800.000 cycles of 20 sec

**battery:**

rechargeable Lithium-ion battery  
(no memory-effect)

**capacity of the accumulator:**

240 cycles of 10 sec

**charging time:**

4 hours (240 min)

**warranty:**

3 years, battery excluded

\* measured with CURE RITE® radiometer by Dentsply

# → STARLIGHT PRO

Lightweight, powerful, versatile



## → STARLIGHT PRO

1 handpiece, 1 charging station

Color

□ white, shiny

■ yellow, shiny

■ silver

■ black, matt  
carbon fiber look

Ref. no.

05100083

05100085

05100087

05100086



The lithium-ion battery of the cordless **starlight pro** provides 320 working cycles of 10 seconds – after only 120 minutes of charging. The power of more than 1.400 mW/cm<sup>2</sup> hardens a 2-mm thick layer in only 10 seconds. All this with only 105 g.



#### → TECHNICAL DATA

**sizes (l x Ø in mm):**

190 x 23

**weight:**

105 g

**light source:**

high power LED

**wavelength:**

dominant wavelength 440 – 465 nm

**optical fiber:**

Ø 8 mm black, click-clack-connection, sterilisable at 135° C

**light intensity\*:**

Ø 8 mm optical fiber, black > 1.400 mW/cm<sup>2</sup>

**average life span LED:**

1.800.000 cycles of 20 sec

**accumulator:**

lithium-ion-accumulator (no memory-effect)

**capacity of the accumulator:**

320 cycles of 10 sec

**charging time:**

about 120 min

**warranty:**

3 years (light guide and battery excluded)

\* measured with CURE RITE® radiometer by Dentsply

## → TECHNICAL DATA

**sizes (l x Ø in mm):**

141 x 23

**weight:**

102 g

**wavelength:**

440 – 465 nm

**optical fiber:**

Ø 8 mm, sterilisable at 135° C

**light intensity\*:**

Ø 8 mm optical fiber > 1.400 mW/cm<sup>2</sup>

**average life span LED:**

1.800.000 cycles of 10 sec

**warranty:**

3 years (light guide excluded)

## → STARLIGHT S+

Product

starlight s+

Ref. no.

05100098

\* measured with CURE RITE® radiometer by Dentsply

# → STARLIGHT S+

flexibility integrated - 3 operating modes  
easy handling

starlight s+ can be easily built into your dental unit and is simply placed in the handpiece holder – exactly where you need it.



# **STARLIGHT SX5**

Built-in kit – the integrated solution



The alternative for dental units. Easily integrated and simply handled.

## **TECHNICAL DATA**

**sizes (l x Ø in mm):**

148 x 22

**weight**

65 g

**wavelength:**

440 – 465 nm

**optical fiber:**

Ø 8 mm, click-clack-connection,  
sterilizable at 135° C

**light intensity\*:**

Ø 8 mm optical fiber > 1.400 mW/cm<sup>2</sup>

**average life span LED:**

1.800.000 cycles of 20 sec

**warranty:**

3 years (light guide excluded)

## **STARLIGHT SX5**

Product

starlight sx5

Ref. no.

05100003

\* measured with CURE RITE® radiometer by Dentsply



optical protection, 8mm



black optical fiber with metal click-clack connection Ø 8 mm



transparent optical fiber Ø 8 mm



black optical fiber Ø 8 mm



transparent optical fiber Ø 4,5 mm

## SPARE PARTS

### SPARE PARTS STARLIGHT UNO

Product

Ref. no.

**battery starlight uno**

03020222

**protective sleeves** (500 pieces)

03230015

**optical protection** (package with 2 pieces)

02900130

### COMMON SPARE PARTS STARLIGHT PRO, S, P AND SLER AND BLUELIGHT

**optical protection, sterilizable Ø 8 mm** (package of 2 pcs.)

02900018

**optical protection, sterilizable Ø 4.5 mm** (package of 2 pcs.)

02900019

### SPARE PARTS STARLIGHT PRO AND ORTHO with metal click-clack connection of the optical fiber

**black optical fiber with metal click-clack connection Ø 8 mm**

02900086

### SPARE PARTS STARLIGHT PRO, S+ AND SX5

**transparent optical fiber Ø 8 mm**

02750009

**black optical fiber Ø 8 mm**

02750012

**transparent optical fiber Ø 4.5 mm**

02750013

**starlight sx5 cord**

02830029

### SPARE PARTS STARLIGHT SLER®

**transparent optical fiber with metal click-clack connection Ø 8 mm**

02900085

**black optical fiber with metal click-clack connection Ø 8mm**

02900086







# APPROPRIATE SETTINGS

| INSERT | PIEZOSURGERY® <i>touch</i><br>PIEZOSURGERY® <i>white</i> | PIEZOSURGERY® 3 |
|--------|----------------------------------------------------------|-----------------|
|        | FUNCTION                                                 | MODE POWER      |

| IMPLANT PREP |          |      |         |
|--------------|----------|------|---------|
| IM1S         | special* | BONE | SPECIAL |
| IM1 AL       | implant  | IMPL | IMPLANT |
| IM2A         | implant  | IMPL | IMPLANT |
| IM2A-15      | implant  | IMPL | IMPLANT |
| IM2P         | implant  | IMPL | IMPLANT |
| IM2P-15      | implant  | IMPL | IMPLANT |
| IM2.8A       | implant  | IMPL | IMPLANT |
| IM2.8P       | implant  | IMPL | IMPLANT |
| IM3A         | implant  | IMPL | IMPLANT |
| IM3A -15     | implant  | IMPL | IMPLANT |
| IM3P         | implant  | IMPL | IMPLANT |
| IM3P-15      | implant  | IMPL | IMPLANT |
| IM3.4A       | implant  | IMPL | IMPLANT |
| IM3.4P       | implant  | IMPL | IMPLANT |
| IM4A         | implant  | IMPL | IMPLANT |
| IM4P         | implant  | IMPL | IMPLANT |
| P2-3         | implant  | IMPL | IMPLANT |
| P3-4         | implant  | IMPL | IMPLANT |

| OSTEOTOMY |                     |      |                    |
|-----------|---------------------|------|--------------------|
| OT1       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT1A      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT2       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT4       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT5       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT5A      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT5B      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT6       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT7       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT7A      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |

| INSERT | PIEZOSURGERY® <i>touch</i><br>PIEZOSURGERY® <i>white</i> | PIEZOSURGERY® 3 |
|--------|----------------------------------------------------------|-----------------|
|        | FUNCTION                                                 | MODE POWER      |

| OSTEOTOMY |                     |      |                    |
|-----------|---------------------|------|--------------------|
| OT7S-3    | special*            | BONE | SPECIAL*           |
| OT7S-4    | special*            | BONE | SPECIAL*           |
| OT7-20    | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT8L      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT8R      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT9       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT11      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT12      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT12S     | special*            | BONE | SPECIAL*           |
| OT13      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OT14      | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| SLO-H     | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| PL1       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| PL2       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| PL3       | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |

| OSTEOPLASTY |                     |      |                    |
|-------------|---------------------|------|--------------------|
| OP1         | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OP2         | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OP3         | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OP3A        | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OP4         | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OP5         | perio               | ROOT | PERIO              |
| OP5A        | perio               | ROOT | PERIO              |
| OP6         | perio               | ROOT | PERIO              |
| OP7         | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| OP8         | special/perio       | ROOT | PERIO              |
| OP9         | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |
| SLC         | cortical/cancellous | BONE | CORTICAL/SPONGIOUS |

# FOR PIEZOSURGERY® INSERTS

| INSERT | PIEZOSURGERY® <i>touch</i><br>PIEZOSURGERY® <i>white</i> | PIEZOSURGERY® 3 |
|--------|----------------------------------------------------------|-----------------|
|        | FUNCTION                                                 | MODE POWER      |

## SINUS MEMBRANE ELEVATION

|      |          |               |
|------|----------|---------------|
| EL1  | endo*    | ROOT ENDO*    |
| EL2  | endo*    | ROOT ENDO*    |
| EL3  | endo*    | ROOT ENDO*    |
| SLS  | special* | BONE SPECIAL* |
| SLE1 | special* | BONE SPECIAL* |
| SLE2 | special* | BONE SPECIAL* |

## PERIOSTEUM PREPARATION

|     |                     |                         |
|-----|---------------------|-------------------------|
| PR1 | cortical/cancellous | BONE CORTICAL/SPONGIOUS |
| PR2 | cortical/cancellous | BONE CORTICAL/SPONGIOUS |

## ENDO RETRO

|      |       |            |
|------|-------|------------|
| EN1  | endo* | ROOT ENDO* |
| EN2  | endo* | ROOT ENDO* |
| EN3  | endo* | ROOT ENDO* |
| EN4  | endo* | ROOT ENDO* |
| EN5R | endo* | ROOT ENDO* |
| EN5L | endo* | ROOT ENDO* |
| EN6R | endo* | ROOT ENDO* |
| EN6L | endo* | ROOT ENDO* |

## EXTRACTION

|     |                     |                         |
|-----|---------------------|-------------------------|
| EX1 | cortical/cancellous | BONE CORTICAL/SPONGIOUS |
| EX2 | cortical/cancellous | BONE CORTICAL/SPONGIOUS |
| EX3 | cortical/cancellous | BONE CORTICAL/SPONGIOUS |

## EXPLANTATION

|        |                             |                                 |
|--------|-----------------------------|---------------------------------|
| EXP3-R | special/cancellous/cortical | BONE SPECIAL/SPONGIOUS/CORTICAL |
| EXP3-L | special/cancellous/cortical | BONE SPECIAL/SPONGIOUS/CORTICAL |
| EXP4-R | special/cancellous/cortical | BONE SPECIAL/SPONGIOUS/CORTICAL |
| EXP4-L | special/cancellous/cortical | BONE SPECIAL/SPONGIOUS/CORTICAL |

| INSERT | PIEZOSURGERY® <i>touch</i><br>PIEZOSURGERY® <i>white</i> | PIEZOSURGERY® 3 |
|--------|----------------------------------------------------------|-----------------|
|        | FUNCTION                                                 | MODE POWER      |

## PERIO

|      |       |            |
|------|-------|------------|
| PP1  | perio | ROOT PERIO |
| PS1  | perio | ROOT PERIO |
| PS2  | perio | ROOT PERIO |
| PS6  | perio | ROOT PERIO |
| PP10 | perio | ROOT PERIO |
| PP11 | perio | ROOT PERIO |
| PP12 | perio | ROOT PERIO |

## IMPLANT CLEANING

|           |            |                 |
|-----------|------------|-----------------|
| ICP + IC1 | endo/perio | ROOT ENDO/PERIO |
|-----------|------------|-----------------|

## RESTORATIVE

|     |               |              |
|-----|---------------|--------------|
| DB2 | special/perio | BONE SPECIAL |
|     |               | ROOT PERIO   |

\* Maximum power permitted

# APPROPRIATE SETTINGS



| function | endo | perio / scaler | restorative | endo | perio / scaler | restorative | endo | perio / scaler | prosthetics |  |
|----------|------|----------------|-------------|------|----------------|-------------|------|----------------|-------------|--|
|----------|------|----------------|-------------|------|----------------|-------------|------|----------------|-------------|--|

## SCALING

|                        |  |       |  |       |  |       |  |       |  |        |
|------------------------|--|-------|--|-------|--|-------|--|-------|--|--------|
| S1-S                   |  | 1-2** |  | 1-2** |  | 1-2** |  | 1-2** |  | 0-30%  |
| S1, S2, S3, S4, S5, S6 |  | 1-6** |  | 1-6** |  | 1-3** |  |       |  | 0-100% |

## PERIO

|                                         |  |       |  |       |  |       |  |       |  |       |
|-----------------------------------------|--|-------|--|-------|--|-------|--|-------|--|-------|
| P1, P2, P3, P4, P10, P11, P12, P13, P14 |  | 1-3** |  | 1-3** |  | 1-3** |  | 1-3** |  | 0-50% |
| P15, P16R, P16L                         |  | 1-3** |  | 1-3** |  | 1-3** |  |       |  | 0-50% |

## IMPLANT CLEANING

|           |  |       |  |       |  |       |  |  |  |       |
|-----------|--|-------|--|-------|--|-------|--|--|--|-------|
| ICS + IC1 |  | 1-2** |  | 1-2** |  | 1-2** |  |  |  | 0-30% |
|-----------|--|-------|--|-------|--|-------|--|--|--|-------|

## ENDO

|        |     |  |     |  |       |  |  |  |  |       |
|--------|-----|--|-----|--|-------|--|--|--|--|-------|
| E1, E2 | 1-6 |  | 1-5 |  | 1-2** |  |  |  |  | 0-10% |
|--------|-----|--|-----|--|-------|--|--|--|--|-------|

# FOR ULTRASOUND INSERTS



| function | endo | perio / scaler | restorative | endo | perio / scaler | restorative | endo | perio / scaler | prosthetics |
|----------|------|----------------|-------------|------|----------------|-------------|------|----------------|-------------|
|----------|------|----------------|-------------|------|----------------|-------------|------|----------------|-------------|

## ENDO REVISION

|                    |  |     |  |     |  |     |  |        |
|--------------------|--|-----|--|-----|--|-----|--|--------|
| ER1                |  | 1-6 |  | 1-5 |  | 1-3 |  | 0-100% |
| ER2, ER3, ER4, ER5 |  | 1-3 |  | 1-3 |  | 1-3 |  | 0-50%  |

## ENDO RETRO

|               |  |     |  |     |  |     |  |       |
|---------------|--|-----|--|-----|--|-----|--|-------|
| RD3, RD4, RD5 |  | 1-3 |  | 1-3 |  | 1-3 |  | 0-50% |
|---------------|--|-----|--|-----|--|-----|--|-------|

## CROWN PREP

|     |  |     |  |     |  |     |  |         |
|-----|--|-----|--|-----|--|-----|--|---------|
| DB1 |  | 4-6 |  | 4-5 |  | ○** |  | 50-100% |
|-----|--|-----|--|-----|--|-----|--|---------|

## EXTRACTION

|                   |  |   |  |   |  |  |  |  |
|-------------------|--|---|--|---|--|--|--|--|
| ME1, ME2, ME3 *** |  | 6 |  | 5 |  |  |  |  |
|-------------------|--|---|--|---|--|--|--|--|

\* These devices are not provided with a power level indicator. The percentage value indicated in the table refers to the adjustment percentage of the scale of the potentiometer supplied with the device.

\*\* Can be used with or without "SOFT MODE"-function.

\*\*\* not available in the United States and in Italy market.


# → MECTRON ONLINE EDUCATION

We are happy to share with you our online education platform.

It will allow you to attend free of charge our live and on-demand educational webinars and view our clinical video library.

A fantastic possibility to take a deeper view on clinical opportunities helpful for your daily practice.

Visit our mectron education platform and register to receive notifications about upcoming webinars.





Free of charge registration  
<https://education.mectron.com/>

# **DISTRIBUTORS**

Please check also our permanently updated distributor list under [www.mectron.com](http://www.mectron.com).

| COUNTRY              | DEALER                                          | PHONE             | E-MAIL                     | MECTRON PRODUCTS | PIEZOSURGERY® PRODUCTS |
|----------------------|-------------------------------------------------|-------------------|----------------------------|------------------|------------------------|
| → ALBANIA            | Rodio Pharma                                    | +355 68 2066757   | info@rodioharma.al         | ●                | ●                      |
| → ALBANIA            | Dental Leader Corsi                             | +355 69 2097616   | mema_ed@yahoo.it           | ●                | ●                      |
| → ALGERIA            | Dentina                                         | +213 02 3772515   | dentinadental@hotmail.com  | ●                | ●                      |
| → ARGENTINA          | Biopolimeros SA                                 | +54 351 424 5528  | ventas@biopolimeros.com.ar | ●                | ●                      |
| → AUSTRALIA          | BioMeDent Pty Ltd                               | +61 1300 792624   | info@biomedent.com.au      | ●                | ●                      |
| → AUSTRIA            | mectron Deutschland Vertriebs GmbH              | +49 221 492015 0  | info@mectron.de            | ●                | ●                      |
| → AZERBAIJAN         | ADENTAA Dental Supplier                         | +994 12 5251117   | info@adentta.az            | ●                | ●                      |
| → BELGIUM            | Dental Promotion & Innovation SA                | +32 25 581700     | info@dentalpromotion.be    | ●                | ●                      |
| → BOLIVIA            | GENTE DENTAL                                    | +32 67 844284     | vegadent23@hotmail.com     | ●                | ●                      |
| → BOSNIA-HERZEGOVINA | Mikodental BH d.o.o.                            | +381 15 345123    | mikodental@gmail.com       | ●                | ●                      |
| → BRAZIL             | Implantec - Comércio e Assistência Técnica      | +55 31 34983326   | jarbas@implantecmg.com.br  |                  | ●                      |
| → BULGARIA           | Medical Depot Rik Ood                           | +359 2 8052500    | office@mdrk.bg             | ●                | ●                      |
| → BULGARIA           | Medicalpark Ltd                                 | +359 878 55 09 04 | import@medicalpark.bg      | ●                | ●                      |
| → CANADA             | Piezosurgery® Inc. / Mectron US                 | +1 614 459 4922   | info@piezosurgery.us       | ●                | ●                      |
| → CHILE              | Berghem                                         | +56 2 2144643     | info@berghem.cl            | ●                | ●                      |
| → CHINA              | Shanghai Tian Ying Medical Instrument Co., Ltd. | +86 21 51060620   | sales@tydental.com         | ●                | ●                      |
| → COLOMBIA           | REPRESENTACIONES A&D LTDA.                      | +57 9156636       | info@ayd.com.co            | ●                | ●                      |
| → COSTA RICA         | ODONTOWEB PROFESIONAL SA                        | +506 22 34 31 33  | info@odontoweb.com         | ●                | ●                      |
| → CROATIA            | Dental Grupa D.O.O.                             | +385 52 752106    | info@dentalgrupa.hr        | ●                | ●                      |
| → CZECH REPUBLIC     | mectron.cz s.r.o.                               | +420 777 761107   | info@mectron.cz            | ●                | ●                      |
| → DENMARK            | Eltident APS                                    | +45 70271747      | tina@eltident.com          | ●                | ●                      |
| → ECUADOR            | PRODENTEC CIA. LTDA.                            | +593 72831232     | ventas@prodentec.com.ec    | ●                | ●                      |



| COUNTRY     | DEALER                              | PHONE            | E-MAIL                      | MECTRON PRODUCTS | PIEZOSURGERY® PRODUCTS |
|-------------|-------------------------------------|------------------|-----------------------------|------------------|------------------------|
| → EGYPT     | Zaident for import dental supply    | +20 111 222 5076 | zaident@link.net            | ●                | ●                      |
| → ESTONIA   | Prodent OÜ                          | +372 656 5982    | info@prodent.ee             | ●                | ●                      |
| → ESTONIA   | Retent                              | +372 7 301 400   | retent@retent.ee            | ●                | ●                      |
| → FINLAND   | Hammasväline Oy                     | +358 10 588 6011 | leila.aalto@hammasvaline.fi | ●                |                        |
| → FINLAND   | Plandent Oy                         | +358 20 77952210 | markku.nappari@plandent.com | ●                |                        |
| → FRANCE    | mectron France                      | +39 0185 1759069 | mectronfrance@mectron.fr    | ●                | ●                      |
| → GEORGIA   | UNIMEDI LTD                         | +995 595900344   | info@unimedi.ge             | ●                | ●                      |
| → GERMANY   | mectron Deutschland Vertriebs GmbH  | +49 221 492015 0 | info@mectron.de             | ●                | ●                      |
| → GREECE    | Green Dental Ike                    | +30 210 3803055  | siakavaras.g@gmail.com      | ●                | ●                      |
| → GUATEMALA | IMFOHSA                             | +502 2317 2100   | info@imfohsa.com            | ●                | ●                      |
| → HONG KONG | Denco (HK) Medical Supplies Ltd.    | +852 2493 7112   | hugo@denco.com.hk           | ●                | ●                      |
| → HUNGARY   | Dental Unio Kft                     | +36 96 518047    | info@dentalunio.hu          | ●                |                        |
| → HUNGARY   | Medi-cont Kft.                      | +36 96 516360    | info@medi-cont.hu           | ●                | ●                      |
| → ICELAND   | Dentalvarehuset Spesialistsenter AS | +47 22 333900    | thomas@dentalvarehuset.no   | ●                | ●                      |
| → INDONESIA | PT. Andini Sarana                   | +62 21 4610515   | hadiwira@andini.co.id       | ●                |                        |
| → IRAN      | Nik Nam Teb Ayandeh                 | +98 21 66467007  | info@niknamteb.com          | ●                | ●                      |
| → IRELAND   | General Medical                     | +44 1380 734990  | info@generalmedical.co.uk   | ●                | ●                      |
| → ISRAEL    | Bicon Isreal                        | +972 86450479    | biconisrael@gmail.com       |                  | ●                      |
| → ISRAEL    | E. Dental Ltd.                      | +972 36915507    | elih@e-dental.co.il         | ●                |                        |
| → ITALY     | mectron s.p.a.                      | +39 0185 35361   | mectron@mectron.com         | ●                | ●                      |
| → JAPAN     | Implatex Co., Ltd                   | +81 3 58508555   | itx@itx.co.jp               |                  | ●                      |
| → JAPAN     | Tokyo Dental Industrial Co. Ltd.    | +81 3 38310176   | orii@tokyodental.co.jp      | ●                |                        |

# **DISTRIBUTORS**

| COUNTRY        | DEALER                                | PHONE              | E-MAIL                       | MECTRON<br>PRODUCTS | PIEZOSURGERY®<br>PRODUCTS |
|----------------|---------------------------------------|--------------------|------------------------------|---------------------|---------------------------|
| → JORDAN       | Ibn Rushd Dental Supplies (IRDS)      | +962 655 26162     | tareq@irds.jo                | ●                   | ●                         |
| → KUWAIT       | Al-Sayafe Medical & Pharmacy Supplies | +965 2245 4924     | alsayafe@alsayafemedical.com | ●                   | ●                         |
| → LATVIA       | AS "Recipe Plus"                      | +371 678 1584242   | info@recipe.lv               | ●                   |                           |
| → LEBANON      | Interbio SAL                          | +961 1491891       | chadi.maroun@interbio.me     | ●                   | ●                         |
| → LITHUANIA    | Jsc Med Grupe                         | +370 5 2761377     | giedre@medgrupe.lt           | ●                   | ●                         |
| → LITHUANIA    | Vilimekso Servisas JSC                | +370 5 2796709     | info@vilimeksas.lt           | ●                   |                           |
| → LUXEMBOURG   | T.D.L.s.a.                            | +352 446881        | info@tdlsa.net               | ●                   | ●                         |
| → MALAYSIA     | Galla Dental Cooperation Sdn Bhd      | +603 6259 3800     | info@galladental.com         |                     | ●                         |
| → MEXICO       | DENTADEC S.A. DE C.V.                 | +55 55454520       | daniel@dentadec.com          |                     | ●                         |
| → MOLDOVA      | VIVODENT SRL CENTRUL STOMATOLOCIC     | +373 22 226522     | office@vivodent.md           | ●                   | ●                         |
| → MONTENEGRO   | MIKODENTAL d.o.o.                     | +381 15 345123     | mikodental@gmail.com         | ●                   | ●                         |
| → MOROCCO      | OrthoZenith                           | +212 5 37 67 45 12 | orthozenith@gmail.com        | ●                   | ●                         |
| → NETHERLANDS  | Hofmeester Dental Bv                  | +31 884 100318     | verkoop@hofmeester.nl        | ●                   | ●                         |
| → NICARAGUA    | INDENT S.A.                           | +505 2223 3517     | info@indentsa.com            | ●                   | ●                         |
| → NORWAY       | Dentalvarehuset Spesialistsenter AS   | +47 22 333900      | thomas@dentalvarehuset.no    | ●                   | ●                         |
| → NORWAY       | Technomedics Norge AS                 | +47 6988 7920      | ek@technomedics.no           | ●                   | ●                         |
| → PAKISTAN     | Glow Pak International                | +92 42 58460902    | info@glowpak.com             | ●                   | ●                         |
| → POLAND       | FM Produkty Dla Stomatologii          | +48 12 632 19 27   | biuro@fmdental.pl            | ●                   | ●                         |
| → PORTUGAL     | Incotrading S.A.                      | +34 91 3807490     | incotrading@incotrading.net  | ●                   | ●                         |
| → ROMANIA      | S.I.T.E.A. Romania Srl                | +40 215 280320     | fabrizio@sitea.ro            | ●                   | ●                         |
| → RUSSIA       | Anas Ltd                              | +7 495 626 5455    | clients@anas.ru              | ●                   | ●                         |
| → SAUDI ARABIA | Cigalah healthcare company            | +966 12 231 6194   | ifarra@cigalah.com.sa        | ●                   | ●                         |

| COUNTRY                | DEALER                                                | PHONE             | E-MAIL                             | MECTRON<br>PRODUCTS | PIEZOSURGERY®<br>PRODUCTS |
|------------------------|-------------------------------------------------------|-------------------|------------------------------------|---------------------|---------------------------|
| → SERBIA               | MIKODENTAL Depot                                      | +381 15 345123    | mikodental@gmail.com               | ●                   | ●                         |
| → SINGAPORE            | Raydent Supplies (S) Pte Ltd                          | +65 9665 8987     | jeanchng@raydent.com.sg            | ●                   | ●                         |
| → SLOVAKIA             | mectron.cz s.r.o.                                     | +420 777 761107   | info@mectron.cz                    | ●                   | ●                         |
| → SLOVENIA             | Korum D.O.O.                                          | +386 1 5612533    | korum@siol.net                     | ●                   | ●                         |
| → SLOVENIA             | Prodent Int. D.O.O.                                   | +386 1 5204800    | info@prodent.si                    | ●                   | ●                         |
| → SOUTH KOREA          | MD Safe Co. Ltd.                                      | +82 2 334 2815    | jclosa@naver.com                   | ●                   | ●                         |
| → SPAIN                | Incotrading S.A.                                      | +34 91 3807490    | incotrading@incotrading.net        | ●                   | ●                         |
| → SWEDEN               | Dentalmind                                            | +46 346 48800     | marcus@dentalmind.com              | ●                   | ●                         |
| → SWITZERLAND          | Bien-Air Dental SA                                    | +41 032 344 64 64 | sales-ch@bienair.com               | ●                   | ●                         |
| → SYRIA                | EBDA Engineering & Development                        | +963 1157 34444   | management@ebda-eng.net            | ●                   | ●                         |
| → TAIWAN               | Sudent Corporation                                    | +886 2822 80 861  | sudent.co@msa.hinet.net            | ●                   | ●                         |
| → THAILAND             | Udom Medical Equipment Co. Ltd.                       | +662 322 7124     | wisit@udom.co.th                   | ●                   | ●                         |
| → TUNISIA              | New Med Plus                                          | +216 31 561 563   | newmedplus@gmail.com               | ●                   | ●                         |
| → TURKEY               | Tekmil Tibbi Arac                                     | +90 21 2337837    | info@tekmil.com.tr                 | ●                   | ●                         |
| → UKRAINE              | Galit                                                 | +380 352 433807   | office@galit.te.ua                 | ●                   | ●                         |
| → UNITED ARAB EMIRATES | New Al Farwaniya Surgicals & Medical Equipment L.L.C. | +971 2 6775447    | info@alfarwaniya.com               | ●                   | ●                         |
| → UNITED KINGDOM       | General Medical                                       | +44 1380 734990   | info@generalmedical.co.uk          | ●                   | ●                         |
| → USA                  | Piezosurgery® Inc./Mectron US                         | +1 614 459 4922   | info@piezosurgery.us               | ●                   | ●                         |
| → VIETNAM              | Pacific Dental Supply                                 | +84 8 391121270   | customercahcm@dentalpacific.com.vn | ●                   | ●                         |

mectron s.p.a.,  
via Loreto 15/A, 16042 Carasco (Ge), Italia,  
tel +39 0185 35361, fax +39 0185 351374

.....> [www.mectron.com](http://www.mectron.com) or [mectron@mectron.com](mailto:mectron@mectron.com)



© Copyright mectron S.p.A., Carasco, Italy  
All rights reserved. Texts, pictures and graphics of mectron brochures are protected by copyright and other protection laws. Without written approval of mectron S.p.A. the contents may not be copied, distributed, changed or made available to third parties for commercial purposes.

The following products are registered under MDR 2017/745. For graphic reasons the styling of the product names might be different. The registered product names correspond as following:

- |                      |                     |                          |                        |
|----------------------|---------------------|--------------------------|------------------------|
| • PIEZOSURGERY touch | PIEZOSURGERY® touch | • COMPACT PIEZO LED      | COMPACT PIEZO LED      |
| • PIEZOSURGERY touch | PIEZOSURGERY® touch | • COMPACT PIEZO P2K      | COMPACT PIEZO P2K      |
| • combi touch        | COMBI touch         | • starjet                | STARJET                |
| • multipiezo         | MULTIPIEZO          | • SUBGINGIVAL PERIO TIPS | Subgingival Perio Tips |
| • MICROPIEZO         | MICROPIEZO          |                          |                        |