

ULTIMATE CHAIRSIDE CORRECTION APPLIANCE



Class II & III Corrector

Class II Setup



Class III Setup



RESULTS MATTER

BEFORE



AFTER



Class III X-Rays: Case photos compliments of Dr. Andy Hayes of St. Louis, MO


CS[®]5 Components



TWIST-LOCK PIVOT
Slot Size
.022



SELF-LIGATING PIVOT
Slot Size
.022



TRADITIONAL PIVOT
Sizes Available
.018 x .025 & .022 x .025



CS[®]5 SPRING
Sizes Available
7mm & 10mm



HYBRID SCREW WITH NYLOK™



MAGNETIC CHARGER

Charges the Magnetic Allen Pencil.
Please remove after charge.

MAGNETIC ALLEN PENCIL

CS[®]5 5-Patient Kit



Important Note: Any previous CS[®] Systems (CS[®]200 and CS[®]3) must not be intermixed with new CS[®]5 parts.



www.dynaflex.com | 866.346.5665



SCAN TO
LEARN MORE