


# The Place Where Technology Meets Art.

STARES FAREE PROCESS!





Krausz Remount-Guide System.

# SIMPLIFY YOUR IMPLANT PROSTHETIC WORKFLOW.



### Remount Guide™ - for Peace of Mind:

- Smile design for optimal aesthetics.
- Prosthetic driven planning.
- Team communication tools.
- Strong durable temporary restoration.
- High-end aesthetic.
- Perfect balanced occlusion.
- Reduce chair-time by over 50% by using Remount locator pins.
- Increase patient comfort.
- Achieve patient satisfaction.

### Maximize the guide:

- Remount (guide A) seats securely on remaining teeth facilitating positioning of location pins.
- Remount (guide B) now seats precisely following full extraction of the teeth utilizing location pin sites.

### **Advantages:**

- Secure, stable accurate surgical guide.
- Prevents the need for free-hand surgery in extraction sites.
- Ensures predictable implant placement and positioning in accordance with plan.
- Excellent prosthetic outcome.

### At Krausz Dental Solutions, we take an holistic approach:

With the use of advanced digital tools and diagnostics we evaluate aesthetic requirements and bone structure.

Following a complete digital mock-up we determine:

- Ideal tooth position.
- Optimal implant planning.
- Choice of prosthetic components.
- Ultimate aesthetic solution.

No matter how complex the case, with Remount we provide you with streamlined procedures and an aesthetic outcome.

### Our method is simple, productive and cost-effective:

- Step-by-step streamlined procedure.
- Eliminates unnecessary visits.
- Easy and simple procedure even when the patient is under sedation.
- Cut treatment time by half.
- Eliminates unwanted laboratory procedures.

With years of professional experience in computer guided planning and prosthetic planning, Krausz Dental Solutions is a "one stop shop".



Important Notice: Suitable for use with all implant systems.





# Benefits of Our Solution.






## Planning -

The Remount-Guide combines full arch design with implant planning.

Prosthetic driven procedure ensures optimal aesthetics and function



### Control -

Stay in control!

A predictable and safe guided system offering the clinician full control throughout the planning and surgical stages. Intuitive software and hardware.



### Time efficient -

Save valuable chair-time when fitting the immediate temporary bridge using our advanced Remount system and prosthetic location pins.



### Accuracy -

Perfect occlusion, aesthetics and functionality.

Predictable, efficient and safe.

## Pre-Plan: Total Digital Workflow.









## Step by Step.



Pin-Guide





Step 2

**Implant Remount-Guide** 





Step 3

**Remount-Bridge** 







\*Available for pilot drill or fully guided



## The Place Where Technology

## **Meets Art.**

#### **Krausz Dental Ltd.**

Cu

+972-3-574-7263



solutions@krausz.co.il



www.krausz.co.il

9

138 Jabotinsky St., Ramat Gan, Israel


Scan the QR code and register our app for easy ordering!

Please feel free to contact us with any questions or concerns.





Premium Prosthetics





10 YEAR WARRANTY!

# ZIRCONIA BRIDGE WITH TITANIUM CORE.



## **Prisma Benefits.**



### Unique titanium infrastructure


with Zirconia overlay offers long-term reliability and prevents fractures.





#### Passive Fit -


Controlled digital workflow and high end milling procedures guarantee passive fit.





### Improved Aesthetics -

Our advanced CAD/CAM manufacturing eliminates unwanted screw access channels in the aesthetic zone.





### 10 Year Warranty\* -

Confidence in our product.
Your precision Prisma bridge is
delivered with 10 year warranty.



## Advantages of Our Solution.

### Prisma Bridge - Zirconia Bridge with Titanium Core



- A monolithic zirconia overlay on precision milled titanium bar delivered with a 10 year warranty\*.
- A fully digital workflow offering a predictable result. 3D design file
   + smile design are shared with you for viewing and approval. A
   diagnostic try-in will be sent for verification of aesthetic + seating
   prior to the final Prisma bridge.
- Simple, streamlined and user friendly process.
- Precision milled titanium bar with unique design gives optimal support for the zirconia overlay.
- The titanium bar offers excellent biocompatibility with the implant / abutment interface and the surrounding soft tissue.
- Precision manufacturing, high accuracy, excellent passive fit.
- Reduce lateral forces due to splinting of implants.
- Prisma bridge offers angled screw access channels up to 30 degrees.
- The restoration may be at implant or abutment level.
- Suitable for both partial or full arch restoration.

# Simple Digital Workflow.

01

Acquisition of digital data - Intra oral scan & Clinical images



Preview of Prisma bridge and smile design.





Try-in of verification 
An aluminum bar with resin overlay.

Final Prisma bridge with zirconia monolithic overlay or milled hybrid PMMA.
















**INCREASE PRODUCTIVITY** 

## **PMMA Bridge**

Acrylic overlay suitable for use against full arch acrylic restoration.





## **Final Bridge**

Zirconia Bridge - Recommended for all types of restorations, with a 10 year warranty\*. High stability restoration.





\*The warranty is given according to the company's terms and conditions





## The Place Where Technology

## **Meets Art.**

#### **Krausz Dental Ltd.**

Cu

+972-3-574-7263



solutions@krausz.co.il



www.krausz.co.il



138 Jabotinsky St., Ramat Gan, Israel

Scan the QR code and register our app for easy ordering!

Please feel free to contact us with any questions or concerns.



