

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
Imes-Icore	0.6	3	NC	Ball end with 1 flute	PMMA/WAX
Imes-Icore	1.0	3	NC	Ball end with 1 flute	PMMA/WAX
Imes-Icore	2.5	3	NC	Ball end with 1 flute	PMMA/WAX
Imes-Icore	0.3	3	DLC	Ball end	ZrO2
Imes-Icore	0.6	3	DLC	Ball end	ZrO2
Imes-Icore	1.0	3	DLC	Ball end	ZrO2
Imes-Icore	1.5	3	DLC	Ball end	ZrO2
Imes-Icore	2.0	3	DLC	Ball end	ZrO2
Imes-Icore	2.5	3	DLC	Ball end	ZrO2
Imes-Icore	0.5	3	DLC	Flat end	ZrO2
Imes-Icore	1.5	3	DLC	Flat end	ZrO2
Imes-Icore	0.6	3	DC	Ball end	ZrO2
Imes-Icore	1.0	3	DC	Ball end	ZrO2
Imes-Icore	1.5	3	DC	Ball end	ZrO2
Imes-Icore	2.0	3	DC	Ball end	ZrO2
Imes-Icore	2.5	3	DC	Ball end	ZrO2
Imes-Icore	0.5	3	DC	Flat end	ZrO2
Imes-Icore	1.5	3	DC	Flat end	ZrO2
Imes-Icore	0.6	3	GD	Ball end	ZrO2
Imes-Icore	1.0	3	GD	Ball end	ZrO2
Imes-Icore	1.5	3	GD	Ball end	ZrO2
Imes-Icore	2.0	3	GD	Ball end	ZrO2
Imes-Icore	2.5	3	GD	Ball end	ZrO2
Imes-Icore	0.5	3	GD	Flat end	ZrO2
Imes-Icore	1.5	3	GD	Flat end	ZrO2
Imes-Icore	0.6	3	Diamond	Ball end	E.max
Imes-Icore	1.0	3	Diamond	Ball end	E.max
Imes-Icore	2.5	3	Diamond	Ball end	E.max
Imes-Icore	0.6	6	NC	Ball end with 1 flute	PMMA/WAX
Imes-Icore	1.0	6	NC	Ball end with 1 flute	PMMA/WAX
Imes-Icore	2.5	6	NC	Ball end with 1 flute	PMMA/WAX
Imes-Icore	0.3	6	DLC	Ball end	ZrO2
Imes-Icore	0.6	6	DLC	Ball end	ZrO2
Imes-Icore	1.0	6	DLC	Ball end	ZrO2
Imes-Icore	2.0	6	DLC	Ball end	ZrO2
Imes-Icore	2.5	6	DLC	Ball end	ZrO2
Imes-Icore	1.5	6	DLC	Flat end	ZrO2
Imes-Icore	0.6	6	DC	Ball end	ZrO2
Imes-Icore	1.0	6	DC	Ball end	ZrO2
Imes-Icore	2.0	6	DC	Ball end	ZrO2
Imes-Icore	2.5	6	DC	Ball end	ZrO2
Imes-Icore	0.6	6	ALN	Ball end	Metal
Imes-Icore	1.0	6	ALN	Ball end	Metal
Imes-Icore	1.5	6	ALN	Ball end	Metal
Imes-Icore	2.0	6	ALN	Ball end	Metal
Imes-Icore	3.0	6	ALN	Ball end	Metal
Imes-Icore	1.0	6	ALN	Flat end	Metal

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
Imes-Icore	1.5	6	ALN	Flat end	Metal
Imes-Icore	2.0	6	ALN	Flat end	Metal
Imes-Icore	1.5	6	ALN	Torus end	Metal
Imes-Icore	2.0	6	ALN	Torus end	Metal
Imes-Icore	1.5	6	ALN	Drill Mill	Metal
Imes-Icore	2.2	6	ALN	Drill Mill	Metal
Imes-Icore	2.5	6	ALN	Drill Mill	Metal
Imes-Icore	0.6	6	Diamond	Ball end	E.max
Imes-Icore	1.0	6	Diamond	Ball end	E.max
Imes-Icore	2.5	6	Diamond	Ball end	E.max
Roland	0.5	4	NC	Ball end	PMMA/WAX
Roland	0.6	4	NC	Ball end	PMMA/WAX
Roland	0.8	4	NC	Ball end	PMMA/WAX
Roland	1.0	4	NC	Ball end	PMMA/WAX
Roland	2.0	4	NC	Ball end	PMMA/WAX
Roland	0.3	4	DLC	Ball end	ZrO2
Roland	0.5	4	DLC	Ball end	ZrO2
Roland	0.6	4	DLC	Ball end	ZrO2
Roland	0.8	4	DLC	Ball end	ZrO2
Roland	1.0	4	DLC	Ball end	ZrO2
Roland	1.5	4	DLC	Ball end	ZrO2
Roland	2.0	4	DLC	Ball end	ZrO2
Roland	2.5	4	DLC	Ball end	ZrO2
Roland	0.6	4	DLC	Flat end	ZrO2
Roland	0.5	4	DLC	Flat end	ZrO2
Roland	1.0	4	DLC	Flat end	ZrO2
Roland	1.5	4	DLC	Flat end	ZrO2
Roland	2.0	4	DLC	Flat end	ZrO2
Roland	0.3	4	DC	Ball end	ZrO2
Roland	0.5	4	DC	Ball end	ZrO2
Roland	0.6	4	DC	Ball end	ZrO2
Roland	0.8	4	DC	Ball end	ZrO2
Roland	1.0	4	DC	Ball end	ZrO2
Roland	1.5	4	DC	Ball end	ZrO2
Roland	2.0	4	DC	Ball end	ZrO2
Roland	2.5	4	DC	Ball end	ZrO2
Roland	0.6	4	DC	Flat end	ZrO2
Roland	0.5	4	DC	Flat end	ZrO2
Roland	1.0	4	DC	Flat end	ZrO2
Roland	1.5	4	DC	Flat end	ZrO2
Roland	2.0	4	DC	Flat end	ZrO2
Roland	0.3	4	GD	Ball end	ZrO2
Roland	0.5	4	GD	Ball end	ZrO2
Roland	0.6	4	GD	Ball end	ZrO2
Roland	0.8	4	GD	Ball end	ZrO2
Roland	1.0	4	GD	Ball end	ZrO2

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
Roland	2.0	4	GD	Ball end	ZrO2
Roland	0.6	4	Diamond	Ball end	E.max
Roland	1.0	4	Diamond	Ball end	E.max
Roland	2.5	4	Diamond	Ball end	E.max
Roland	1.0	4	ALN	Ball end	Metal
Roland	2.0	4	ALN	Ball end	Metal
Roland	3.0	4	ALN	Ball end	Metal
Zirkonzahn	0.5	3	NC	Ball end	PMMA/WAX
Zirkonzahn	1.0	3	NC	Ball end	PMMA/WAX
Zirkonzahn	2.0	3	NC	Ball end	PMMA/WAX
Zirkonzahn	0.3	3	NC	Ball end	ZrO2
Zirkonzahn	0.5	3	NC	Ball end	ZrO2
Zirkonzahn	0.6	3	NC	Ball end	ZrO2
Zirkonzahn	1.0	3	NC	Ball end	ZrO2
Zirkonzahn	2.0	3	NC	Ball end	ZrO2
Zirkonzahn	0.6	3	NC	Flat end	ZrO2
Zirkonzahn	1.5	3	NC	Flat end	ZrO2
Zirkonzahn	0.3	3	DLC	Ball end	ZrO2
Zirkonzahn	0.5	3	DLC	Ball end	ZrO2
Zirkonzahn	0.6	3	DLC	Ball end	ZrO2
Zirkonzahn	1.0	3	DLC	Ball end	ZrO2
Zirkonzahn	2.0	3	DLC	Ball end	ZrO2
Zirkonzahn	0.6	3	DLC	Flat end	ZrO2
Zirkonzahn	1.5	3	DLC	Flat end	ZrO2
Zirkonzahn	0.3	3	DC	Ball end	ZrO2
Zirkonzahn	0.5	3	DC	Ball end	ZrO2
Zirkonzahn	0.6	3	DC	Ball end	ZrO2
Zirkonzahn	1.0	3	DC	Ball end	ZrO2
Zirkonzahn	2.0	3	DC	Ball end	ZrO2
Zirkonzahn	0.6	3	DC	Flat end	ZrO2
Zirkonzahn	1.5	3	DC	Flat end	ZrO2
Zirkonzahn	0.5	6	NC	Ball end	PMMA/WAX
Zirkonzahn	1.0	6	NC	Ball end	PMMA/WAX
Zirkonzahn	2.0	6	NC	Ball end	PMMA/WAX
Zirkonzahn	0.3	6	NC	Ball end	ZrO2
Zirkonzahn	0.5	6	NC	Ball end	ZrO2
Zirkonzahn	0.6	6	NC	Ball end	ZrO2
Zirkonzahn	1.0	6	NC	Ball end	ZrO2
Zirkonzahn	2.0	6	NC	Ball end	ZrO2
Zirkonzahn	0.6	6	NC	Flat end	ZrO2
Zirkonzahn	1.5	6	NC	Flat end	ZrO2
Zirkonzahn	0.3	6	DLC	Ball end	ZrO2
Zirkonzahn	0.5	6	DLC	Ball end	ZrO2
Zirkonzahn	0.6	6	DLC	Ball end	ZrO2
Zirkonzahn	1.0	6	DLC	Ball end	ZrO2
Zirkonzahn	2.0	6	DLC	Ball end	ZrO2
Zirkonzahn	0.6	6	DLC	Flat end	ZrO2

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
Zirkonzahn	1.5	6	DLC	Flat end	ZrO2
Zirkonzahn	0.3	6	DC	Ball end	ZrO2
Zirkonzahn	0.5	6	DC	Ball end	ZrO2
Zirkonzahn	0.6	6	DC	Ball end	ZrO2
Zirkonzahn	1.0	6	DC	Ball end	ZrO2
Zirkonzahn	2.0	6	DC	Ball end	ZrO2
Zirkonzahn	0.6	6	DC	Flat end	ZrO2
Zirkonzahn	1.5	6	DC	Flat end	ZrO2
Zirkonzahn	0.6	6	Diamond	Ball end	E.max
Zirkonzahn	1.25	6	Diamond	Ball end	E.max
Zirkonzahn	2.5	6	Diamond	Ball end	E.max
Zirkonzahn	1.0	6	ALN	Ball end	Metal
Zirkonzahn	1.5	6	ALN	Ball end	Metal
Zirkonzahn	2.0	6	ALN	Ball end	Metal
Zirkonzahn	3.0	6	ALN	Ball end	Metal
Zirkonzahn	1.5	6	ALN	Torus end	Metal
Zirkonzahn	2.0	6	ALN	Torus end	Metal
Zirkonzahn	3.0	6	ALN	Torus end	Metal
Zirkonzahn	1.5	6	ALN	Flat end	Metal
Zirkonzahn	2.0	6	ALN	Flat end	Metal
Zirkonzahn	3.0	6	ALN	Flat end	Metal
VHF 40	0.6	3	NC	Ball end	PMMA/WAX
VHF 40	1.0	3	NC	Ball end	PMMA/WAX
VHF 40	2.0	3	NC	Ball end	PMMA/WAX
VHF 40	2.5	3	NC	Ball end with 1 flute, and German raw material	PMMA/WAX
VHF 35/40	0.3	3	NC	Ball end	ZrO2
VHF 35/40	0.6	3	NC	Ball end	ZrO2
VHF 35/40	1.0	3	NC	Ball end	ZrO2
VHF 35/40	2.0	3	NC	Ball end	ZrO2
VHF 35/40	0.3	3	DLC	Ball end	ZrO2
VHF 35/40	0.6	3	DLC	Ball end	ZrO2
VHF 35/40	1.0	3	DLC	Ball end	ZrO2
VHF 35/40	2.0	3	DLC	Ball end	ZrO2
VHF 35/40	0.6	3	DLC	Flat end	ZrO2
VHF 35/40	1.2	3	DLC	Flat end	ZrO2
VHF 35/40	0.3	3	DC	Ball end	ZrO2
VHF 35/40	0.6	3	DC	Ball end	ZrO2
VHF 35/40	1.0	3	DC	Ball end	ZrO2
VHF 35/40	2.0	3	DC	Ball end	ZrO2
VHF 35/40	1.2	3	DC	Flat end	ZrO2
VHF 35/40	0.6	3	DC	Flat end	ZrO2

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
VHF 35/40	0.3	3	GD	Ball end	ZrO2
VHF 35/40	0.6	3	GD	Ball end	ZrO2
VHF 35/40	1.0	3	GD	Ball end	ZrO2
VHF 35/40	2.0	3	GD	Ball end	ZrO2
VHF 35	0.6	3	Diamond	Ball end	E.max
VHF 35	1.0	3	Diamond	Ball end	E.max
VHF 35	2.4/2.6	3	Diamond	2.4 ball/2.6 flat	E.max
VHF 32/35	0.6	3	ALN	Ball end	Metal
VHF 32/35	1.0	3	ALN	Ball end	Metal
VHF 32/35	2.0	3	ALN	Ball end with 2/4 flutes	Metal
VHF 32/35	1.2	3	ALN	Flat end	Metal
Wieland Select	0.3	3	NC	Ball end	PMMA/WAX
Wieland Select	0.7	3	NC	Ball end	PMMA/WAX
Wieland Select	1.0	3	NC	Ball end	PMMA/WAX
Wieland Select	2.5	3	NC	Ball end	PMMA/WAX
Wieland Mini	0.3	3	NC	Ball end	ZrO2
Wieland Mini	0.7	3	NC	Ball end	ZrO2
Wieland Mini	1.0	3	NC	Ball end	ZrO2
Wieland Mini	2.5	3	NC	Ball end	ZrO2
Wieland Select	0.3	3	NC	Ball end	ZrO2
Wieland Select	0.7	3	NC	Ball end	ZrO2
Wieland Select	1.0	3	NC	Ball end	ZrO2
Wieland Select	2.5	3	NC	Ball end	ZrO2
Wieland Select	1.2	3	NC	Flat end	ZrO2
Wieland Select	0.6	3	NC	Flat end	ZrO2
Wieland Mini/Select	0.3	3	DLC	Ball end	ZrO2
Wieland Mini/Select	0.7	3	DLC	Ball end	ZrO2
Wieland Mini/Select	1.0	3	DLC	Ball end	ZrO2
Wieland Mini/Select	2.5	3	DLC	Ball end	ZrO2
Wieland Mini/Select	1.2	3	DLC	Flat end	ZrO2
Wieland Mini/Select	0.6	3	DLC	Flat end	ZrO2
Wieland Select	0.3	3	DC	Ball end	ZrO2
Wieland Select	0.7	3	DC	Ball end	ZrO2
Wieland Select	1.0	3	DC	Ball end	ZrO2
Wieland Select	2.5	3	DC	Ball end	ZrO2
Wieland Select	1.2	3	DC	Flat end	ZrO2
Wieland Select	0.6	3	DC	Flat end	ZrO2
Wieland Select	0.3	3	GD	Ball end	ZrO2
Wieland Select	0.7	3	GD	Ball end	ZrO2
Wieland Select	1.0	3	GD	Ball end	ZrO2
Wieland Select	2.5	3	GD	Ball end	ZrO2
Wieland Mini	0.6	3	Diamond	Ball end	E.max
Wieland Mini	1.0	3	Diamond	Ball end	E.max
Wieland Mini	2.0	3	Diamond	Ball end	E.max
Wieland Mini	3.0	3	Diamond	Ball end	E.max

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
Amann Girrbach	0.3	3	NC	Ball end	PMMA/WAX
Amann Girrbach	0.6	3	NC	Ball end	PMMA/WAX
Amann Girrbach	1.0	3	NC	Ball end	PMMA/WAX
Amann Girrbach	2.5	3	NC	Ball end	PMMA/WAX
Amann Girrbach	0.3	3	DLC	Ball end	ZrO2
Amann Girrbach	0.6	3	DLC	Ball end	ZrO2
Amann Girrbach	1.0	3	DLC	Ball end	ZrO2
Amann Girrbach	2.5	3	DLC	Ball end	ZrO2
Amann Girrbach	0.6	3	DLC	Flat end	ZrO2
Amann Girrbach	1.2	3	DLC	Flat end	ZrO2
Amann Girrbach	0.3	3	DC	Ball end	ZrO2
Amann Girrbach	0.6	3	DC	Ball end	ZrO2
Amann Girrbach	1.0	3	DC	Ball end	ZrO2
Amann Girrbach	2.5	3	DC	Ball end	ZrO2
Amann Girrbach	0.6	3	DC	Flat end	ZrO2
Amann Girrbach	1.2	3	DC	Flat end	ZrO2
Amann Girrbach	0.3	3	GD	Ball end	ZrO2
Amann Girrbach	0.6	3	GD	Ball end	ZrO2
Amann Girrbach	1.0	3	GD	Ball end	ZrO2
Amann Girrbach	2.5	3	GD	Ball end	ZrO2
Amann Girrbach	0.6	3	GD	Flat end	ZrO2
Amann Girrbach	1.2	3	GD	Flat end	ZrO2
Amann Girrbach	0.4	3	Diamond	Ball end	E.max
Amann Girrbach	1.0	3	Diamond	Ball end	E.max
Amann Girrbach	1.4	3	Diamond	Flat end	E.max
Amann Girrbach	1.8	3	Diamond	Ball end	E.max
Amann Girrbach	1.0	3	ALN	Ball end	Metal
Amann Girrbach	2.0	3	ALN	Ball end	Metal
XTCERA	0.6	3	NC	Ball end	PMMA/WAX
XTCERA	1.0	3	NC	Ball end	PMMA/WAX
XTCERA	2.0	3	NC	Ball end	PMMA/WAX
XTCERA	0.6	4	NC	Ball end	PMMA/WAX
XTCERA	1.0	4	NC	Ball end	PMMA/WAX
XTCERA	2.0	4	NC	Ball end	PMMA/WAX
XTCERA	0.6	3	DLC	Ball end	ZrO2
XTCERA	1.0	3	DLC	Ball end	ZrO2
XTCERA	2.0	3	DLC	Ball end	ZrO2
XTCERA	0.6	4	DLC	Ball end	ZrO2
XTCERA	1.0	4	DLC	Ball end	ZrO2
XTCERA	2.0	4	DLC	Ball end	ZrO2
XTCERA	0.6	3	DC	Ball end	ZrO2
XTCERA	1.0	3	DC	Ball end	ZrO2
XTCERA	2.0	3	DC	Ball end	ZrO2
XTCERA	0.6	4	DC	Ball end	ZrO2
XTCERA	1.0	4	DC	Ball end	ZrO2
XTCERA	2.0	4	DC	Ball end	ZrO2
XTCERA	0.6	3	GD	Ball end	ZrO2

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
XTCERA	1.0	3	GD	Ball end	ZrO2
XTCERA	2.0	3	GD	Ball end	ZrO2
XTCERA	0.6	4	GD	Ball end	ZrO2
XTCERA	1.0	4	GD	Ball end	ZrO2
XTCERA	2.0	4	GD	Ball end	ZrO2
XTCERA	0.6	3	Diamond	Ball end	E.max
XTCERA	1.0	3	Diamond	Ball end	E.max
XTCERA	2.5	3	Diamond	Ball end	E.max
XTCERA	0.6	4	Diamond	Ball end	E.max
XTCERA	1.0	4	Diamond	Ball end	E.max
XTCERA	2.5	4	Diamond	Ball end	E.max
XTCERA	0.6	6	Diamond	Ball end	E.max
XTCERA	1.0	6	Diamond	Ball end	E.max
XTCERA	2.5	6	Diamond	Ball end	E.max
XTCERA	M1.4	6	ALN	Screw end	Metal
XTCERA	M1.6	6	ALN	Screw end	Metal
XTCERA	M1.8	6	ALN	Screw end	Metal
XTCERA	M2.0	6	ALN	Screw end	Metal
XTCERA	1.5	4	ALN	Ball end	Metal
XTCERA	2.0	4	ALN	Ball end	Metal
XTCERA	3.0	4	ALN	Ball end	Metal
XTCERA	1.0	6	ALN	Ball end	Metal
XTCERA	2.0	6	ALN	Ball end	Metal
XTCERA	3.0	6	ALN	Ball end	Metal
XTCERA	1.5x14	6	ALN	Torus end	Metal
XTCERA	1.5x18	6	ALN	Torus end	Metal
XTCERA	2.0	6	ALN	Torus end	Metal
XTCERA	0.5	6	ALN	Flat end	Metal
XTCERA	1.0	6	ALN	Flat end	Metal
XTCERA	1.5	6	ALN	Flat end	Metal
XTCERA	2.0	6	ALN	Flat end	Metal
XTCERA	2.2	6	ALN	Drill Mill	Metal
XTCERA	1.5	6	ALN	Drill Mill	Metal
SIRONA MCX5	0.5	3	NC	Ball end	PMMA/WAX
SIRONA MCX5	1.0	3	NC	Ball end	PMMA/WAX
SIRONA MCX5	2.5	3	NC	Ball end	PMMA/WAX
SIRONA MCX5	0.5	3	DLC	Ball end	ZrO2
SIRONA MCX5	1.0	3	DLC	Ball end	ZrO2
SIRONA MCX5	2.5	3	DLC	Ball end	ZrO2
SIRONA MCX5	0.5	3	DC	Ball end	ZrO2
SIRONA MCX5	1.0	3	DC	Ball end	ZrO2
SIRONA MCX5	2.5	3	DC	Ball end	ZrO2
SIRONA MCX5	0.5	3	GD	Ball end	ZrO2
SIRONA MCX5	1.0	3	GD	Ball end	ZrO2
SIRONA MCX5	2.5	3	GD	Ball end	ZrO2
SIRONA MCXL	1.0	N/A	Diamond	N/A	E.max
SIRONA MCXL	1.4	N/A	Diamond	N/A	E.max

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
SIRONA MCXL	1.8	N/A	Diamond	N/A	E.max
SIRONA MCX5	0.6	3	Diamond	Ball end	E.max
SIRONA MCX5	1.2	3	Diamond	Ball end	E.max
SIRONA MCX5	1.4	3	Diamond	Flat end	E.max
SIRONA MCX5	2.2	3	Diamond	Flat end	E.max
ARUM with ring	0.6	4	NC	Ball end	PMMA/WAX
ARUM with ring	1.0	4	NC	Ball end	PMMA/WAX
ARUM with ring	1.5	4	NC	Ball end	PMMA/WAX
ARUM with ring	2.0	4	NC	Ball end	PMMA/WAX
ARUM without ring	0.6	6	NC	Ball end	PMMA/WAX
ARUM without ring	1.0	6	NC	Ball end	PMMA/WAX
ARUM without ring	2.0	6	NC	Ball end	PMMA/WAX
ARUM without/with ring	0.3	4	DLC	Ball end	ZrO2
ARUM without/with ring	0.6	4	DLC	Ball end	ZrO2
ARUM without/with ring	1.0	4	DLC	Ball end	ZrO2
ARUM without/with ring	2.0	4	DLC	Ball end	ZrO2
ARUM without ring	0.6	6	DLC	Ball end	ZrO2
ARUM without ring	1.0	6	DLC	Ball end	ZrO2
ARUM without ring	2.0	6	DLC	Ball end	ZrO2
ARUM without/with ring	0.3	4	DC	Ball end	ZrO2
ARUM without/with ring	0.6	4	DC	Ball end	ZrO2
ARUM without/with ring	1.0	4	DC	Ball end	ZrO2
ARUM without/with ring	2.0	4	DC	Ball end	ZrO2
ARUM without ring	0.6	6	DC	Ball end	ZrO2
ARUM without ring	1.0	6	DC	Ball end	ZrO2
ARUM without ring	2.0	6	DC	Ball end	ZrO2
ARUM without/with ring	0.3	4	GD	Ball end	ZrO2
ARUM without/with ring	0.6	4	GD	Ball end	ZrO2
ARUM without/with ring	1.0	4	GD	Ball end	ZrO2
ARUM without/with ring	2.0	4	GD	Ball end	ZrO2
ARUM without ring	1.5	6	ALN	Ball end	Metal
ARUM without ring	2.0	6	ALN	Ball end	Metal
ARUM without ring	3.0	6	ALN	Ball end	Metal
ARUM without ring	1.0	6	ALN	Ball end	Metal
ARUM without ring	0.6	6	ALN	Ball end	Metal
ARUM with ring	1.5	4	ALN	Ball end	Metal

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
ARUM with ring	2.0	4	ALN	Ball end	Metal
ARUM with ring	1.0	4	ALN	Ball end	Metal
ARUM without ring	1.0	6	ALN	Flat end	Metal
ARUM without ring	1.5	6	ALN	Flat end	Metal
ARUM without ring	2.0	6	ALN	Flat end	Metal
ARUM with ring	0.6	4	Diamond	Ball end	E.max
ARUM with ring	1.0	4	Diamond	Ball end	E.max
ARUM with ring	1.5	4	Diamond	Ball end	E.max
ARUM with ring	2.5	4	Diamond	Ball end	E.max
ARUM without ring	0.6	6	Diamond	Ball end	E.max
ARUM without ring	1.0	6	Diamond	Ball end	E.max
ARUM without ring	1.5	6	Diamond	Ball end	E.max
ARUM without ring	2.5	6	Diamond	Ball end	E.max
YENADENT	0.6	4	NC	Ball end	PMMA/WAX
YENADENT	1.0	4	NC	Ball end	PMMA/WAX
YENADENT	2.0	4	NC	Ball end	PMMA/WAX
YENADENT	0.6	4	DLC	Ball end	ZrO2
YENADENT	1.0	4	DLC	Ball end	ZrO2
YENADENT	2.0	4	DLC	Ball end	ZrO2
YENADENT	0.6	4	DC	Ball end	ZrO2
YENADENT	1.0	4	DC	Ball end	ZrO2
YENADENT	2.0	4	DC	Ball end	ZrO2
YENADENT	0.6	4	GD	Ball end	ZrO2
YENADENT	1.0	4	GD	Ball end	ZrO2
YENADENT	2.0	4	GD	Ball end	ZrO2
YENADENT	1.5	4	ALN	Ball end	Metal
YENADENT	2.0	4	ALN	Ball end	Metal
YENADENT	3.0	4	ALN	Ball end	Metal
YENADENT	1.0	4	ALN	Ball end	Metal
YENADENT	1.4	4	ALN	Flat end	Metal
YENADENT	1.8	4	ALN	Flat end	Metal
YENADENT	0.6	4	Diamond	Ball end	E.max
YENADENT	1.0	4	Diamond	Ball end	E.max
YENADENT	2.0	4	Diamond	Ball end	E.max
Dentium	0.5	3	NC	Ball end	PMMA/WAX
Dentium	1.0	3	NC	Ball end	PMMA/WAX
Dentium	2.0	3	NC	Ball end	PMMA/WAX
Dentium	0.5	6	NC	Ball end	PMMA/WAX
Dentium	1.0	6	NC	Ball end	PMMA/WAX
Dentium	2.0	6	NC	Ball end	PMMA/WAX
Dentium	0.5	3	DLC	Ball end	ZrO2
Dentium	1.0	3	DLC	Ball end	ZrO2
Dentium	2.0	3	DLC	Ball end	ZrO2
Dentium	0.5	4	DLC	Ball end	ZrO2
Dentium	1.0	4	DLC	Ball end	ZrO2
Dentium	2.0	4	DLC	Ball end	ZrO2
Dentium	0.5	3	DC	Ball end	ZrO2

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
Dentium	1.0	3	DC	Ball end	ZrO2
Dentium	2.0	3	DC	Ball end	ZrO2
Dentium	0.5	4	DC	Ball end	ZrO2
Dentium	1.0	4	DC	Ball end	ZrO2
Dentium	2.0	4	DC	Ball end	ZrO2
Dentium	0.5	3	GD	Ball end	ZrO2
Dentium	1.0	3	GD	Ball end	ZrO2
Dentium	2.0	3	GD	Ball end	ZrO2
Dentium	0.5	4	GD	Ball end	ZrO2
Dentium	1.0	4	GD	Ball end	ZrO2
Dentium	2.0	4	GD	Ball end	ZrO2
Dentium	0.6	3	Diamond	Ball end	E.max
Dentium	1.0	3	Diamond	Ball end	E.max
Dentium	2.0	3	Diamond	Ball end	E.max
Dentium	0.6	6	Diamond	Ball end	E.max
Dentium	1.0	6	Diamond	Ball end	E.max
Dentium	2.0	6	Diamond	Ball end	E.max
Dentium	1.5	6	ALN	Ball end	Metal
Dentium	2.0	6	ALN	Ball end	Metal
Dentium	3.0	6	ALN	Ball end	Metal
Dent Mill	0.3	3.175	DLC	Ball end	ZrO2
Dent Mill	0.6	3.175	DLC	Ball end	ZrO2
Dent Mill	0.8	3.175	DLC	Ball end	ZrO2
Dent Mill	1.0	3.175	DLC	Ball end	ZrO2
Dent Mill	2.0	3.175	DLC	Ball end	ZrO2
Dent Mill	3.0	3.175	DLC	Ball end	ZrO2
Dent Mill	0.3	3.175	DC	Ball end	ZrO2
Dent Mill	0.6	3.175	DC	Ball end	ZrO2
Dent Mill	0.8	3.175	DC	Ball end	ZrO2
Dent Mill	1.0	3.175	DC	Ball end	ZrO2
Dent Mill	2.0	3.175	DC	Ball end	ZrO2
Dent Mill	3.0	3.175	DC	Ball end	ZrO2
Dent Mill	0.3	3.175	GD	Ball end	ZrO2
Dent Mill	0.6	3.175	GD	Ball end	ZrO2
Dent Mill	0.8	3.175	GD	Ball end	ZrO2
Dent Mill	1.0	3.175	GD	Ball end	ZrO2
Dent Mill	2.0	3.175	GD	Ball end	ZrO2
Dent Mill	3.0	3.175	GD	Ball end	ZrO2
3Shape Lava 38/50	0.6	3	NC	Ball end	PMMA/WAX
3Shape Lava 38/50	1.0	3	NC	Ball end	PMMA/WAX
3Shape Lava 38/50	2.0	3	NC	Ball end	PMMA/WAX
3Shape Lava 38/50	0.6	3	DLC	Ball end	ZrO2
3Shape Lava 38/50	1.0	3	DLC	Ball end	ZrO2
3Shape Lava 38/50	2.0	3	DLC	Ball end	ZrO2
3Shape Lava 38/50	0.6	3	DC	Ball end	ZrO2
3Shape Lava 38/50	1.0	3	DC	Ball end	ZrO2
3Shape Lava 38/50	2.0	3	DC	Ball end	ZrO2

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
3Shape Lava 38/50	0.6	3	GD	Ball end	ZrO2
3Shape Lava 38/50	1.0	3	GD	Ball end	ZrO2
3Shape Lava 38/50	2.0	3	GD	Ball end	ZrO2
Kavo	1.0	6	DC	Ball end	ZrO2
Kavo	2.0	6	DC	Ball end	ZrO2
Upcera 38/40	0.3	3	NC	Ball end	PMMA/WAX
Upcera 38/40	0.6	3	NC	Ball end	PMMA/WAX
Upcera 38/40	1.0	3	NC	Ball end	PMMA/WAX
Upcera 38/40	2.0	3	NC	Ball end	PMMA/WAX
Upcera 40	1.0	3	NC	Flat end	PMMA/WAX
Upcera 40	2.0	3	NC	Flat end	PMMA/WAX
Upcera 38/40	0.3	3	DLC	Ball end	ZrO2
Upcera 38/40	0.6	3	DLC	Ball end	ZrO2
Upcera 38/40	1.0	3	DLC	Ball end	ZrO2
Upcera 38/40	2.0	3	DLC	Ball end	ZrO2
Upcera 38	1.5	3	DLC	Ball end	ZrO2
Upcera 38/40	1.5	3	DLC	Flat end	ZrO2
Upcera 38/40	1.0	3	DLC	Flat end	ZrO2
Upcera 38/40	2.0	3	DLC	Flat end	ZrO2
Upcera 38/40	0.3	3	DC	Ball end	ZrO2
Upcera 38/40	0.6	3	DC	Ball end	ZrO2
Upcera 38/40	1.0	3	DC	Ball end	ZrO2
Upcera 38/40	2.0	3	DC	Ball end	ZrO2
Upcera 38	1.5	3	DC	Ball end	ZrO2
Upcera 38/40	1.5	3	DC	Flat end	ZrO2
Upcera 38/40	1.0	3	DC	Flat end	ZrO2
Upcera 38/40	2.0	3	DC	Flat end	ZrO2
Upcera 38/40	0.3	3	GD	Ball end	ZrO2
Upcera 38/40	0.6	3	GD	Ball end	ZrO2
Upcera 38/40	1.0	3	GD	Ball end	ZrO2
Upcera 38/40	2.0	3	GD	Ball end	ZrO2
Upcera	0.6	6	Diamond	Ball end	E.max
Upcera	1.0	6	Diamond	Ball end	E.max
Upcera	2.5	6	Diamond	Ball end	E.max
Upcera	0.6	3	Diamond	Ball end	E.max
Upcera	1.0	3	Diamond	Ball end	E.max
Upcera	2.5	3	Diamond	Ball end	E.max
Upcera	1.0	6	ALN	Ball end	Metal
Upcera	1.5	6	ALN	Ball end	Metal
Upcera	2.0	6	ALN	Ball end	Metal
Upcera	3.0	6	ALN	Ball end	Metal
Upcera	1.0	6	ALN	Flat end	Metal
Upcera	2.0	6	ALN	Flat end	Metal
Upcera	1.5	6	ALN	Torus end	Metal
Upcera	2.0	6	ALN	Torus end	Metal
Upcera	1.5	6	ALN	Drill Mill	Metal
Upcera	2.5	6	ALN	Drill Mill	Metal

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
Adite AMD	0.6	4	NC	Ball end with 1 flute	PMMA/WAX
Adite AMD	1.0	4	NC	Ball end with 1 flute	PMMA/WAX
Adite AMD	2.0	4	NC	Ball end with 1 flute	PMMA/WAX
Adite AMD	0.3	4	DLC	Ball end	ZrO2
Adite AMD	0.6	4	DLC	Ball end	ZrO2
Adite AMD	1.0	4	DLC	Ball end	ZrO2
Adite AMD	2.0	4	DLC	Ball end	ZrO2
Adite AMD	1.0-Long	4	DLC	Ball end	ZrO2
Adite AMD	2.0-Long	4	DLC	Ball end	ZrO2
Adite AMD	1.5	4	DLC	Flat end	ZrO2
Adite AMD	0.3	4	DC	Ball end	ZrO2
Adite AMD	0.6	4	DC	Ball end	ZrO2
Adite AMD	1.0	4	DC	Ball end	ZrO2
Adite AMD	2.0	4	DC	Ball end	ZrO2
Adite AMD	1.0-Long	4	DC	Ball end	ZrO2
Adite AMD	2.0-Long	4	DC	Ball end	ZrO2
Adite AMD	1.5	4	DC	Flat end	ZrO2
Adite AMD	0.3	4	GD	Ball end	ZrO2
Adite AMD	0.6	4	GD	Ball end	ZrO2
Adite AMD	1.0	4	GD	Ball end	ZrO2
Adite AMD	2.0	4	GD	Ball end	ZrO2
Adite AMD	1.5	6	ALN	Ball end	Metal
Adite AMD	1.0	6	ALN	Ball end	Metal
Adite AMD	3.0	6	ALN	Ball end	Metal
Adite AMD	0.6	6	Diamond	Ball end	E.max
Adite AMD	1.0	6	Diamond	Ball end	E.max
Adite AMD	2.5	6	Diamond	Ball end	E.max
Dental Plus	0.5	4	DLC	Ball end	ZrO2
Dental Plus	1.0	4	DLC	Ball end	ZrO2
Dental Plus	2.0	4	DLC	Ball end	ZrO2
Dental Plus	0.5	4	DC	Ball end	ZrO2
Dental Plus	1.0	4	DC	Ball end	ZrO2
Dental Plus	2.0	4	DC	Ball end	ZrO2
Dental Plus	0.5	4	GD	Ball end	ZrO2
Dental Plus	1.0	4	GD	Ball end	ZrO2
Dental Plus	2.0	4	GD	Ball end	ZrO2
Dental Plus	1.0	4	ALN	Ball end	Metal
Dental Plus	1.5	4	ALN	Ball end	Metal
Dental Plus	2.0	4	ALN	Ball end	Metal
Dental Plus	3.0	4	ALN	Ball end	Metal
Dental Plus	1	4	ALN	Flat end	Metal
Dental Plus	1.5	4	ALN	Flat end	Metal
Dental Plus	1.5	4	ALN	Torus end	Metal
Dental Plus	1	4	ALN	Torus end	Metal
Dental Plus	0.6	4	Diamond	Ball end	E.max
Dental Plus	1.0	4	Diamond	Ball end	E.max
Dental Plus	2.0	4	Diamond	Ball end	E.max

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
Dental Plus	2.5	4	Diamond	Ball end	E.max
PM 3/5	0.5	3	NC	Ball end	PMMA/WAX
PM 3/5	1.0	3	NC	Ball end	PMMA/WAX
PM 3/5	2.5	3	NC	Ball end	PMMA/WAX
PM 7	0.6	6	NC	Ball end with 1 flute	PMMA/WAX
PM 7	1.0	6	NC	Ball end with 1 flute	PMMA/WAX
PM 7	1.5	6	NC	Ball end with 1 flute	PMMA/WAX
PM 7	2.5	6	NC	Ball end with 1 flute	PMMA/WAX
PM 7	5.0	6	NC	Ball end with 1 flute	PMMA/WAX
PM 3/5	0.5	3	DLC	Ball end	ZrO2
PM 3/5	1.0	3	DLC	Ball end	ZrO2
PM 3/5	2.5	3	DLC	Ball end	ZrO2
PM 7	0.5	6	DLC	Ball end	ZrO2
PM 7	1.0	6	DLC	Ball end	ZrO2
PM 7	2.5	6	DLC	Ball end	ZrO2
PM 3/5	0.5	3	DC	Ball end	ZrO2
PM 3/5	1.0	3	DC	Ball end	ZrO2
PM 3/5	2.5	3	DC	Ball end	ZrO2
PM 7	0.5	6	DC	Ball end	ZrO2
PM 7	1.0	6	DC	Ball end	ZrO2
PM 7	2.5	6	DC	Ball end	ZrO2
PM 3/5	0.5	3	GD	Ball end	ZrO2
PM 3/5	1.0	3	GD	Ball end	ZrO2
PM 3/5	2.5	3	GD	Ball end	ZrO2
PM 7	0.5	6	Diamond	Ball end	E.max
PM 7	1.0	6	Diamond	Ball end	E.max
PM 7	2.0	6	Diamond	Ball end	E.max
PM 7	3.0	6	Diamond	Ball end	E.max
PM 7	0.6	6	ALN	Ball end	Metal
PM 7	1.0	6	ALN	Ball end	Metal
PM 7	1.5	6	ALN	Ball end	Metal
PM 7	2.0	6	ALN	Ball end	Metal
PM 7	3.0	6	ALN	Ball end	Metal
PM 7	3.0	6	ALN	Flat end	Metal
UP3D	0.3	4	NC	Ball end	PMMA/WAX
UP3D	0.6	4	NC	Ball end	PMMA/WAX
UP3D	1.0	4	NC	Ball end	PMMA/WAX
UP3D	2.0	4	NC	Ball end	PMMA/WAX
UP3D	0.3	4	DLC	Ball end	ZrO2
UP3D	0.6	4	DLC	Ball end	ZrO2
UP3D	1.0	4	DLC	Ball end	ZrO2
UP3D	2.0	4	DLC	Ball end	ZrO2
UP3D	0.3	4	DC	Ball end	ZrO2
UP3D	0.6	4	DC	Ball end	ZrO2
UP3D	1.0	4	DC	Ball end	ZrO2
UP3D	2.0	4	DC	Ball end	ZrO2
UP3D	0.3	4	GD	Ball end	ZrO2

2022 YUTE TOOLS Burs List9.15

info@yutedent.com

whatsapp:+86 18680320038

System	Diameter, mm	Shank, mm	Coating	Type	Materials
UP3D	0.6	4	GD	Ball end	ZrO2
UP3D	1.0	4	GD	Ball end	ZrO2
UP3D	2.0	4	GD	Ball end	ZrO2
UP3D	0.6	4	Diamond	Ball end	E.max
UP3D	1.0	4	Diamond	Ball end	E.max
UP3D	2.0	4	Diamond	Ball end	E.max
DOF	1.0	4	DLC	Ball end	ZrO2
DOF	2.0	4	DLC	Ball end	ZrO2
DOF	1.0	4	DC	Ball end	ZrO2
DOF	2.0	4	DC	Ball end	ZrO2
DOF	1.0	4	Diamond	Ball end	E.max
DOF	2.5	4	Diamond	Ball end	E.max
Origin	0.6	4	DLC	Ball end	ZrO2
Origin	1.0	4	DLC	Ball end	ZrO2
Origin	2.0	4	DLC	Ball end	ZrO2
Origin	0.6	4	DC	Ball end	ZrO2
Origin	1.0	4	DC	Ball end	ZrO2
Origin	2.0	4	DC	Ball end	ZrO2