


Controlled Memory Spring NiTi File

# Spring Endo File Hi Plus

Compared to the existing spring endo file, this specially designed endo file has maximized the flexibility, performance, and safety to allow access inside the root canal by maintaining the bent shape of the shank.



## ● THE BEST FLEXIBILITY



Excellent flexibility complete with both heat treatment and a spring structure

## ● THE BEST PERFORMANCE



Curved root canal



Patients who cannot open their mouth wide




Treating the mesio-buccal canal of the maxillary posterior teeth



Minimizing the canal deformation during canal shaping

# WHY Denflex Spring Endo File Hi Plus?



**Spring Endo File Hi Plus**  
The best file for a curved root canal and posterior teeth



SAFE

## The spring structure minimizes the risk of a fracture inside the canal

- While the load is dispersed towards the spring, the micro vibration of the spring when rotating provides a much easier fragment removal with the file because the force of the fragment being stuck inside the canal is not strong compared to the existing file without spring structure.



EFFICIENT

## The spring structure and heat treatment allow efficient procedure

- **Sufficient cutting force** due to the micro vibration of the spring.
- The spring structure provides 100% rotatory power even when bent at 50~60 degrees.
- The advanced durability provides greater numbers of uses compared to general endo files.
- The maximized flexibility can preserve the original canal shape as much as possible during canal shaping.



COMFORTABLE

## The spring structure provides a highly convenient procedure for both the patient and operator

- The flexible structure facilitates the operation of the patients who cannot open their mouth wide.
- Excessive preparation is unnecessary even when the canal has curved too much curved from the tooth axis.
- **Optimized for a curved root canal.**

## Product Information

Packaging: 6 files



### Spring Endo File Hi Plus of 6 files

CW: Clockwise | CR: Counterclockwise

	23mm	25mm
● 10	DFSEF-B-L2310A2-CW(CR)	DFSEF-B-L2510A2-CW(CR)
○ 15	DFSEF-B-L2315A4-CW(CR)	DFSEF-B-L2515A4-CW(CR)
● 20	DFSEF-B-L2320A6-CW(CR)	DFSEF-B-L2520A6-CW(CR)
● 25	DFSEF-B-L2325A4-CW(CR) DFSEF-B-L2325A6-CW(CR)	DFSEF-B-L2525A4-CW(CR) DFSEF-B-L2525A6-CW(CR)
● 30	DFSEF-B-L2330A4-CW(CR) DFSEF-B-L2330A6-CW(CR)	DFSEF-B-L2530A4-CW(CR) DFSEF-B-L2530A6-CW(CR)
● 40	-	DFSEF-B-L2540A5-CW(CR)
● 50	-	DFSEF-B-L2550A5-CW(CR)

**Denflex**  
www.denflex.kr

**Denflex Co., Ltd.**

#708, Byucksan Digital Valley I, 24, Digital-ro 27-gil,  
Guro-gu, Seoul, Republic of Korea

**Tel.** +82-2-866-0617 **Fax.** +82-2-866-0618 **Email.** info@denflex.kr

