



**YOUR  
IMAGINATION  
IS THE  
LIMIT.**

# MiYO - MAKE IT YOUR OWN!

With the paste-like MiYO liquid ceramic, monolithic restorations can be designed with aesthetics and depth effects that are only known from layered restorations - with considerably less time effort. The minimal layering thickness of only 0.1 - 0.2 mm prevents chipping.

Another special feature of MiYO restorations is their unique light dynamics. Due to the scattering of light, there is no shadowing as known from stains. Therefore, MiYO restorations look natural in the mouth and not only on the model.

Versatile applications make MiYO a real game changer for your restorations: Either MiYO pure - colour and structure - or the smart hybrid technique - the combination of MiYO and InSync ZR layering ceramic.

**Discover MiYO for your lab today!**





Only 0.1 - 0.2 mm layering thickness, therefore no risk of chipping.

Natural aesthetics and depth effects with unique light dynamics - no shadowing.



Detailed design of natural effects for teeth and gingiva. MiYO Colors can be applied wet-on-wet without blurring.



Colour, shape and effects are visible directly when applying MiYO and do not change after firing.



Raise the brightness value by 1-2 levels with the Value Enhancers Lumin or Lumin plus, lower it with Smoke or Slate without losing any characteristics.




Smart hybrid technique - MiYO Color masses in combination with veneering ceramics.

Pre-fab, homogeneous pastes for reproducible results and easy handling.



Achieve dynamic 3D effects with MiYO Structrue without cut-back. Design individual surface texture, tooth shape and occlusal contact points simply with a brush.





# MiYO & MiYO pink I Kit configurations

Choose either the MiYO „Standard“ Kit or for even more creativity the MiYO Kit Color+.

The MiYO pink Kit offers everything you need for gingiva restorations.

All kits include Glaze Paste and Stain & Glaze Liquid.



# Our SERVICE for you

SERVICE is not only written in capital letters for us, but also lived in this way. Because we want our customers to become real MiYO fans. So we offer the right support depending on your needs:

- Free online trainings where the practical application of MiYO and MiYO pink is shown live.
- Training videos, information on furnace programming and much more in the customer portal at [www.miyoworld.eu](http://www.miyoworld.eu)
- Technical support for quick help via mail or telephone.



## Accessories

High quality brushes and tools with exchangeable tips.



## Learn more ...



[www.miyoworld.eu](http://www.miyoworld.eu)



[shop.jensendental.de](http://shop.jensendental.de)



Jensen Dental GmbH  
Gustav-Werner-Str. 1  
72555 Metzingen, Germany  
[info@jensendental.de](mailto:info@jensendental.de)  
+49 (0) 7123 92260  
[www.jensendental.de](http://www.jensendental.de)