



**NEW**

**YOUR  
IMAGINATION  
IS THE  
LIMIT.**

**Ormu**  
**ON**

# 850 MiYO

## **New esthetics for lithium disilicate restorations.**

Especially for crystallization.

The 850 MiYO ceramic system was specifically developed for the characterization of lithium disilicate restorations before the crystallization firing. By eliminating the need for an additional characterization firing, 850 MiYO enables the creation of highly esthetic and naturally looking restorations faster than ever before.

Additionally, its enhanced firing stability allows the use with high-fusing layering ceramics in hybrid techniques. The high color stability of 850 MiYO provides reliability, ensuring that colors do not fade during firing.

Another advantage is the unique MiYO light dynamics, which also distinguishes 850 MiYO restorations. The diffuse light scattering prevents shadows as seen with traditional stains.

Reproducible results without surprises and with natural esthetics are the result.

**Discover 850 MiYO for your lab today!**



## Product features

### Easy to use


- 850 MiYO is easy to apply, even wet-on-wet.
- 850 MiYO doesn't drift, puddle or slump.
- 850 MiYO shines after crystallization firing without additional glaze firing.
- 850 MiYO is optimized for existing crystallization cycles.

### Natural esthetics

- MiYO 850 offers the translucency and vitality of natural dentition.
- MiYO 850 crowns have a finished surface sheen like natural tooth enamel.
- 850 MiYO creates natural color and depth effects.
- 850 MiYO Value Enhancer to control the brightness value.

### Consistent & Predictable

- 850 MiYO is color-stable - no graying after firing.
- With 850 MiYO, the result after firing corresponds to the design before firing.
- MiYO 850 perfectly matches all Vita® Classic shades.



# System components

850 MiYO Translucent Body Shades   fluor	Shade A	Shade B	Shade C	Shade D
850 MiYO Translucent Value Enhancer   fluor	Smoke	Lumin		
850 MiYO Translucent Effect Material   fluor	Straw	Clementine	Storm	Sage
850 MiYO Color   fluor	Snow	Fissure (powder)		
MiYO Structure   fluor	Enamel			
850 Glaze Paste fluor				

## Accessories

High quality brushes and tools with exchangeable tips.



## Learn more ...



[www.miyoworld.eu](http://www.miyoworld.eu)



[shop.jensendental.de](http://shop.jensendental.de)



Jensen Dental GmbH  
Gustav-Werner-Str. 1  
72555 Metzingen, Germany  
[info@jensendental.de](mailto:info@jensendental.de)  
+49 (0)7123 92260  
[www.jensendental.de](http://www.jensendental.de)

