


RAMO MEDICAL

ENDODONTIC CATALOGUE



Make Endodontic Simple Safe Efficient



CONTACT US


Our team's vision is to be a team capable of continuous innovation and continuous progress. We hope that through our efforts, we can bring better services and products to endodontists and make our own contributions to the development of the industry.



Alisa Huo

[Brand Founder]

WhatsApp: 008618855951157
Email: Alisa@ramomedical.com



Cassie Wang

[Marketing Manager]

WhatsApp: 008615955571773
Email: market02@ramomedical.com



Wendy Zhang

[Global Sales Manager]

WhatsApp: 008617802162869
Email: Wendy@ramomedical.com



Beatriz Luo (Spanish)

[Area Sales Manager]


WhatsApp: 008615907388101
Email: Beatriz@ramomedical.com



Jack Liu (English)

[Area Sales Manager]


WhatsApp: 008613013033162
Email: Jack@ramomedical.com



Lizzy Li (English)

[Area Sales Manager]


WhatsApp: 008617856002256
Email: lizzy@ramomedical.com



Mandy Zhang (English)

[Area Sales Manager]

WhatsApp: 008618205601493
Email: Sales001@ramomedical.com



Rebecca Lin (French)

[Area Sales Manager]

WhatsApp: 008618105659671
Email: sales003@ramomedical.com




Icey Jiao (English)

[Area Sales Manager]

WhatsApp: 008619856039223
Email: Sales05@ramomedical.com

FEATURES

- ◎ Flexible.
- ◎ Upgraded slim S-shaped cross section.
- ◎ The first nano coated rainbow file in the world. Maximum increases surface rigidity while improving cyclic fatigue resistance.
- ◎ Advanced polishing technology increases cyclic fatigue resistance by 30%.
- ◎ Greater torsional resistance and flexural strength.

Size	2004	2006	2504	2506	3004	3006	3504	3506
Product								
Length	21/25/31(mm)							
Packing	SX+3PCS 2004/Blister SX+3PCS 2006/Blister		SX+3PCS 2504/Blister SX+3PCS 2506/Blister		SX+3PCS 3004/Blister SX+3PCS 3006/Blister		SX+3PCS 3504/Blister SX+3PCS 3506/Blister	

PRODUCT DETAILS



- Slim S cross section



- Speed&Torque
400RPM/2.0N




IMAGE&PACKING



FEATURES

- ◎ Two safe sides, multi cross section symmetrical design to make whole file dint even, making the cutting more stable and efficient.
- ◎ Upgraded slim S-shaped cross section with two safe sides to release more space for debris with higher efficiency.
- ◎ Increases patient safety with enhanced flexibility and resistance to cyclic fatigue.
- ◎ Respect the original anatomy of the tooth.
- ◎ Unique design.

Size	2004	2006	2504	2506	3004	3006	3504	3506
Product								
Length	21/25/31(mm)							
Packing	SX+3PCS 2004/Blister SX+3PCS 2006/Blister		SX+3PCS 2504/Blister SX+3PCS 2506/Blister		SX+3PCS 3004/Blister SX+3PCS 3006/Blister		SX+3PCS 3504/Blister SX+3PCS 3506/Blister	

PRODUCT DETAILS



- Slim S cross section
with two safe sides



- Speed&Torque
500RPM/2.5N










IMAGE&PACKING



FEATURES

- ◎ Rotary path file.
- ◎ Suitable for curved canals.
- ◎ More flexible than K-files.
- ◎ Faster instrument time than manual K-files.
- ◎ 7 sizes available.

Size	1302	1602	1902	1003	1203	1503	1903
Product							
Length	21/25/31(mm)						
Packing	4PCS/Blister				6PCS/Blister		

PRODUCT DETAILS



- Square cross section



- Speed&Torque
350RPM/1.5N








IMAGE&PACKING



FEATURES

- ◎ Multi-Taper system.
- ◎ SX,S1,S2 are flexible;F1, F2, F3 are extremely flexible, both have memory control.
- ◎ Ensures that shaping file has good cutting efficiency and finishing file can save more dentin.

Size	SX	S1	S2	F1	F2	F3
Product						
Length	21/25/31(mm)					
Packing	6PCS/Blister					

PRODUCT DETAILS



- Triangle cross section



- Speed&Torque
350RPM/2.0N









IMAGE&PACKING



FEATURES

- © Extremely flexible.
- © Suitable for all root canal treatment, including very curved canal.
- © Traditional six files system, easy to operate.

Size	1708	1303	2004	2504	3004	3504
Product						
Length	21/25/31(mm)					
Packing	6PCS/Blister					

PRODUCT DETAILS




- Triangle cross section



- Speed&Torque
350RPM/2.0N

**IMAGE&PACKING****FEATURES**

- © Control Memory.
- © Reciprocating motion, CW30, CCW150.
- © Consistent taper, more stable.
- © The cutting angles and the improved raw material provide high efficiency and cutting performance.
- © Non-cutting tips.

Size	2006	2506	4006
Product			
Length	21/25/31(mm)		
Packing	3PCS/Blister		

PRODUCT DETAILS

- Triangle cross section







- Speed&Torque
CCW150/CW30

**IMAGE&PACKING**

FEATURES

- ◎ Control Memory.
- ◎ The perfect combination of resistance to cyclic fatigue and cutting efficiency.





Size	1708	1805	2506	3504
Product				
Length	21/25/31(mm)			
Packing	4PCS/Blister			

PRODUCT DETAILS


• Triangle cross section		• Speed&Torque 350RPM/2.0N	
--------------------------	---	-------------------------------	---

IMAGE&PACKING**FEATURES**

- ◎ Super flexible.
- ◎ Short length design is more suitable for use on children.
- ◎ Improved cutting efficiency and resistance of cyclic fatigue.

Size	1708	2004	2504	3004
Product				
Length	10mm	16mm		
Packing	4PCS/Blister			


PRODUCT DETAILS

• Triangle cross section		• Speed&Torque 350RPM/2.0N	
--------------------------	---	-------------------------------	---

IMAGE&PACKING

FEATURES

- ⊗ High shaping to remove gutta-percha and paste fillers fast.
- ⊗ 5 sizes to meet all retreatment needs.

Size	1005	1506	2007	2508	3009
Product					
Length	16mm	18mm	22mm	18mm	16mm
Packing	3PCS/Blister			5PCS/Blister	

PRODUCT DETAILS

- Rhombus cross section



- Speed&Torque
350RPM/2.0N

**IMAGE&PACKING****FEATURES**

- ⊗ Cleaning file.
- ⊗ Up to 98% bacterial reduction in the main canal;
- ⊗ 2X more debris removal.
- ⊗ Superior resistance to cyclic fatigue.
- ⊗ Be compatible with all rotary or reciprocating endo systems.
- ⊗ Gently abrade canal walls and dislodge debris in difficult to reach areas.
- ⊗ Non-cutting and non-active tip.

Length	Taper&Size	Packing
25(mm)	1%/#25	4PCS/Pack



- Triangle cross section









- Speed&Torque
800RPM/1.0N

**IMAGE&PACKING**

FEATURES


- © 02 taper design.
Large size files can also safely enter curved root canals, reducing the risk of instrument separation.
- © Protect dentin.
Enlarge the apical hole enough to remove bacteria, and retain the dentin in the middle and upper section of the root canal to avoid excessive. cutting.
- © Effectively solve the problem of insufficient root tip preparation.

Size	3002	3502	4002	4502	5002	6002
Product						
Length	21/25/31(mm)					
Packing	6PCS/Blister					


PRODUCT DETAILS



- Square cross section



- Speed&Torque
300RPM/1.5N



IMAGE&PACKING



Ultrasonic NiTi File

Ultrasonic NiTi File use the principle of ultrasonic vibration to clean root canals and remove calcification in the root canals, and are used in conjunction with ultrasonic root canal treatment instruments.



Size	32.5mm
#15	S18 0101 00015 32.5
#20	S18 0101 00020 32.5
#25	S18 0101 00025 32.5
#30	S18 0101 00030 32.5
#35	S18 0101 00035 32.5
#40	S18 0101 00040 32.5




PRODUCT DETAILS



STEP-BY-STEP INSTRUCTIONS FOR USE OF RCS GP FILES


- (1) Establish the working length and apical stop of the root canal.
- (2) Insert the gutta-percha point and the RCS GP FILE into the canal, and press the RCS GP FILE against the gutta-percha.
- (3) Descend the RCS GP FILE into the canal, rotating clockwise at 8'000 min-1. Do not penetrate it to a distance of less than 2 mm from the apex.
- (4) Back the files out of the canal, without opposing strong resistance.



HAND FILE

RCS SPREADER

RCS Spreader File is tapered and is used to move the gutta-percha laterally to make room for additional gutta-percha cones to be placed in the root canal space.



- For lateral filling for root canal filling materials.
- Hand instrument with a smooth and tapered working end.
- Provide a better sense of tactile and less likely to cause a fracture in the root of the tooth.


RCS C FILE

- Gold Wire Tech.
- Extremely sharp.
- Special for negotiating through a calcified root canal.



RCS SHAPE D

- Stainless steel, rigid.
- Shaped D cross section.
- Use for enlarging root canal and negotiating calcified and narrowed root canal.



RCS H FILE

- High quality stainless steel.
- Comfortable grip, clear product texture.
- High-quality fine steel, high temperature resistance.
- high pressure resistance and high corrosion resistance.
- Not easy to loose, adopt integrated design.



OPTIONAL MODEL

