



**Vonflex S Putty**  
**Proven and Reliable Precision Impression Material**  
 • Kit : Base 280ml x 1bottle, Catalyst 280ml x 1bottle,  
 2measuring scoops



**Vonflex S Heavy**  
**Excellent Fixation Characteristic Material**  
 • Type : Injection(50ml)  
 • Package : - Injection(50ml x 2cartridges), Mixing tips  
 - Injection(50ml x 12cartridges), Mixing tips



**Vonflex S Supermix**  
**Excellent Fixation Characteristic Material**  
 • Type : Heavy Type\_Injection(380ml)  
 • Package : - Injection(380ml x 1cartridge)  
 - Injection(380ml x 2cartridges)  
 - Injection(380ml x 4cartridges)



**Vonflex S Medium**  
**Proven One Step Material Based on A-silicone**  
 • Type : Injection(50ml)  
 • Package : Injection(50ml x 2cartridges), Mixing tips



**Vonflex S Light**  
**Excellent Hydrophilic Impression Material**  
 • Type : Injection(50ml)  
 • Package : - Injection(50ml x 2cartridges), Mixing tips  
 - Injection(50ml x 12cartridges), Mixing tips



**Vonflex S Light XLV**  
**Excellent Hydrophilic Impression Material**  
 • Type : Injection(50ml)  
 • Package : - Injection(50ml x 2cartridges), Mixing tips  
 - Injection(50ml x 12cartridges), Mixing tips



**Vonflex S Bite**  
**Thixotropic Bite Registration Material**  
 • Type : Injection(50ml)  
 • Package : Injection(50ml x 2cartridges), Mixing tips

DISINFECTION SERIES



**Dual Clean**  
 Highly Active  
 Enzymatic Cleaner

• 1L / 2L



**Dual Clean Plus**  
 Highly Active Clean  
 With Disinfectant

• 1L / 2L



**Dual Vac**  
 Highly Active Cleaner,  
 Disinfectant & Deodorant  
 for All Suction Systems

• 1L / 2L



**V-varnish Premium**  
**Clear Cavity varnish with TCP & Xylitol**  
 • Standard pack : 50 single-unit-dose application sticks  
 • Bulk pack : 200 single-unit-dose application sticks  
 • Flavor : Bubblegum, Strawberry, Melon



**Eco-S**  
**Light-Cured Pit & Fissure Sealant**  
 • Economic kit : 1.2ml x 2syringes, DenFil Etchant-37 1.5g x  
 2syringes, Disposable tips  
 • Normal kit : 1.2ml x 4syringes, Disposable tips



**Eazigel**  
**Topical Fluoride Gel**  
 • Excellent flowability for minute area  
 with 1.23% APF fluoride ion  
 • Containing xylitol for healthy teeth  
 • Soft bottle for easy pouring  
 • Hardly overflow out of tray  
 • 480ml(Strawberry and peach flavors)



Eazitray

OTHER PRODUCTS



**V-Ject**  
**Superior Quality Disposable Dental Needles**  
 • 100pcs / Box  
 • Gauge : 27G, 30G



Visit our YouTube channel  
 to find out more videos.

VERICOM CO., LTD.

**Head Office / Factory**  
 48, Toegyegongdan 1-gil, Chuncheon-Si, Gangwon-Do, 24427, Korea

**Sales office / R&D Center**  
 15, Jeonpa-ro 62beon-gil, Manan-gu, Anyang-Si, Gyeonggi-Do, 14086, Korea  
 T. +82-31-441-2881 F. +82-31-441-2883 H. www.vericom.co.kr  
 E-mail. vericom@vericom.co.kr / overseas@vericomdental.com

Dental Products  
 Catalog

Restoratives / Adhesives & Etchants / Preventives /  
 Endodontics / Impressions / Disinfections /  
 Other Products





### DenFil N

#### Light-Cured Restorative Nano Hybrid Composite Resin

- Shade : A1, A2, A3, A3.5, B1, B2, C1, A2O, A3O, BL
- Type : Syringe Kit, Refill (4g x 1syringe), Capsule (0.25g x 20ea)
- Package : **5 Syringe kit (4g x 5syringes)**
  - A pack (A1, A2, A3, A3.5, B2)
  - B pack (A2, A3, A3.5, B1, B2)
  - C pack (A1, A2, A3, A3.5, C1)
- 8 Syringe kit (4g x 8syringes)**
  - A pack (A1, A2, A3, A3.5, B2, A2O, A3O, BL)
  - B pack (A2, A3, A3.5, B1, B2, C1, A2O, BL)
  - C pack (A1, A2, A3, A3.5, B1, C1, A2O, A3O)



### DenFil

#### Light-cured Restorative Hybrid Composite Resin

- Shade : A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D3, A2O, A3O, I, P, E
- Type : Syringe Kit, Refill (4g x 1syringe), Capsule (0.25g x 20ea)
- Package : **5 Syringe kit (4g x 5syringes)**
  - A pack (A1, A2, A3, A3.5, B2)
  - B pack (A2, A3, A3.5, B1, B2)
  - C pack (A1, A2, A3, A3.5, C2)



### DenFil Flow

#### Light-cured Radiopaque Flowable Composite Resin

- Shade : A1, A2, A3, A3.5, B2, B3, C2, A2O, A3O, A2L, A3L, I
- Type : Syringe Kit, Refill(2g x 1syringe)
- Package : **4 syringe kit (2g X 4syringes)**
  - A pack (A2, A3, A3.5, B2)
  - B pack (A2)
  - C pack (A3)
  - D pack (A3.5)



### DenFil Bulk Flow

#### Light-cured Flowable Bulk Fill Restorative Material

- Shade : U, A1, A2, A3, B1
- Type : Syringe(2g x 1syringe)
- Package : 2g x 1syringe, 2g x 2syringes



### Dual Core Automix type

#### Dual-Cured Core Build-up Material

- Mini-mix type : 9g x 2syringes, Mixing tips, Intraoral tips
- Color : Natural, Blue



### U-Cem Premium

#### Self-Adhesive Universal Resin Cement

- Double dispenser type : 9g x 1 or 2syringes, Mixing pad, Spatula
- Automix type : 9g x 1 or 2syringes, Mixing tips, Intraoral tips
- Color : Universal, Clear, Opaque



### Care Cem Temp

#### Dual-cured temporary resin cement

- High bond strength(Strong retention)
- Easy to remove(Easy removal of temporary prosthesis)
- Flow and mixed easily(The viscosity that is suitable for using)
- Automix Type : 5ml x 1, 2 or 3ea(Mixing tips, Mixing pad)
- Handmix Type : 5ml x 1ea, Mixing pad



### BC-plus

#### Single Component Bonding Agent

- 5ml x 1bottle



### U-Bond

#### One Component Self-etching Light-cured Adhesive

- Kit type : 5ml x 1bottle, Brush holder, Brushes, Mixing well
- Refill type : 5ml x 1bottle



Refill 5ml



### DenFil Etchant-37

#### Semi-gel Type Etching agent with Phosphoric Acid

- 37% Phosphoric acid
- Eco-Kit: 3ml x 5syringes, Disposable tips
- Eco-Kit: 3ml x 12syringes, Disposable tips
- Bulk Kit: 5ml x 5syringes, Disposable tips



### V-Etch Viscous solution

#### Novel and Optimum Etchant

- 35% Phosphoric acid
- Eco-Kit: 1.2ml x 4syringes, Disposable tips
- Bulk Kit: 30ml x 1syringe, Disposable tips



V-Etch 30ml



### Well-Root ST

#### Premixed Bioceramic Sealer

- 2g x 1syringe
- Disposable tips, Extension holder



### Well-Root PT

#### Premixed Bioceramic Putty

- 0.25g x 10capsules



### Well-Paste

#### Water-Based Calcium Hydroxide Paste with Barium Sulfate

- 2g x 1syringe
- Disposable tips, Disposable tip cap, Extension holder



### Well-Pex

#### Oil-Based Calcium Hydroxide Paste with Iodoform

- 2g x 1syringe
- Disposable tips, Disposable tip cap, Extension holder



### Well-Prep

#### Creamy EDTA Root Canal Cleaner

- 6g x 3syringes, Disposable tips



### Spacer

#### Light-Cured Temporary Filling Material

- Color : Blue, Yellow
- 3g x 3syringes



### Care C&B

#### Dual-cured temporary crown and bridge material

- Shade : A1, A2, A3, A3.5, B1, BL
- Low shrinkage and Low setting and Low Exothermic Heat
- Fast setting and High Flexural Strength
- More Natural Looking
- Automix Type : 50ml x 1ea, Mixing tips  
20ml x 1ea, Mixing tips



### U-Bond Ortho

#### Light-cured Orthodontic Adhesive

- Color : Sky, Clear
- Kit : U-Bond Ortho 4g x 2syringes, U-Bond Primer 6ml x 1bottle, Instrument 1ea, Brushes
- Refill : U-Bond Ortho 4g x 1syringe
- Refill (Primer) : U-Bond Primer 6ml x 1bottle, 3ml x 1bottle



U-Bond Primer

# MAZIC C-200

## Uniformity Curing Unit

Possible to adjust power & time by optimized curing schedule per each resin.



Visit our YouTube channel to find out more videos.

## Dental Laboratory Products Catalog

Infinite possibilities of CAD/CAM systems



# MAZIC L-200

## MONO 4K LCD Printer

Up to x3 faster and 4K LCD for precise and detailed reproducibility. Setting values are provided for all kinds of Vericom 3D resin materials.



### Heat/Cool

Uniform shade and physical properties are provided through constant internal temperature control



### Light uniformity


Minimum deviation of light inside the equipment




### 4 Curing Setting Option



### Compact Design



Light uniformity



Heat / Cool

- Uniform shade and physical properties are provided through constant internal temperature control
- Fast curing with a high-power light source
- Time by optimized curing schedule per each resin
- 4 Curing Setting Option
- Compact Design

- Light source : 405nm LED
- LED Power Output : 190W
- Product Size : L220 x W240 x H390mm

## VERICOM CO., LTD.

Head Office / Factory  
48, Toegyegongdan 1-gil, Chuncheon-si, Gangwon-do, Korea

Sales Office / R&D Center  
15, Jeonpa-ro 62beon-gil, Manan-gu, Anyang-si, Gyeonggi-do, Korea

Tel +82-31-441-2881 Fax +82-31-441-2883 Website [www.vericom.co.kr](http://www.vericom.co.kr)  
Email [vericom@vericom.co.kr](mailto:vericom@vericom.co.kr) / [overseas@vericomdental.com](mailto:overseas@vericomdental.com)

Dental Laboratory\_ENG\_240527\_ver 08

VERICOM



### Heater



### Uniform light irradiation



### Wide Build Plate

L132 x W78



### Fascinating & Modern Design

- Fast output and high thermal stability with Heater function
- High quality output with uniform light irradiation
- Convenient Replacement
- Wide build plate (132 x 78mm)
- Convenient use as a magnetic build plate & electromagnet resin tank
- Compact Design

- Light source : 405nm LED
- XY Resolution : 35µm
- Printing speed : 40mm/hour
- Product Size : L220 x W240 x H390mm



Heater



Uniform light irradiation  
High-quality print

